


# Bending of a leaf spring with Cobra SMARTsense



Physics

Mechanics

Forces, work, power &amp; energy



Difficulty level

easy



Group size

2



Preparation time

10 minutes



Execution time

10 minutes


This content can also be found online at:




<https://www.curriculab.de/c/604f1a2e3b95820003807eff>

PHYWE

## Teacher information



## Application



Experiment set-up

In this experiment, a leaf spring is to be clamped horizontally and a force applied at the outer end. The applied force then results in a bending moment  $M_b$ , which stresses the spring. The bending moment is greatest at the loading point and decreases towards the support of the leaf spring until it finally becomes zero in the bearing itself.

The bending moment results from the product of the acting force  $F$  and the lever arm  $l$ :

$$M_b = F \cdot l \sim [\text{Nm}]$$

## Teacher information (1/2)

### Prior knowledge



Students should have a basic understanding of forces. Ideally, the students have already carried out the experiment on Hooke's law and know the concept of the spring constant as well as the correlations between the deflection of a spring under a specific force effect.

### Learning objective




The students are to investigate the behaviour or deformation of a leaf spring under load (by the acting force), present the measurement results in the form of a diagram and, in particular, determine the spring constant  $D$ .

## Teacher information (2/2)

PHYWE

### Scientific principle



Measured values and fit function

### About the evaluation

Generate fit line:

Using the straight line tool of the app, a balance line can be laid through the first 5 measuring points. It is important that the first measurement actually corresponds to 1 cm deflection, the 2nd measurement to 2 cm and so on. The equation for the straight line can then be read directly in the output.

## Notes on procedure

PHYWE

- The force sensor must be tared - this is done automatically when the sensor is switched on and off.

Alternatively: Select the sensor in the app. Then select "Set to zero" and save.

- Attention must be paid to the correct height of the force sensor: The height of the free end of the leaf spring decreases with increasing deflection. If the height of the force sensor is not adjusted accordingly, the measurement results will be falsified. Usually the height of the force sensor should be adjusted according to the following deflections: 5 cm, 8 cm, 10 cm, 12 cm, 13 cm, 14 cm, 15 cm.

## Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE



## Student Information

### Motivation

PHYWE



<https://www.youtube.com/watch?v=egwjc5axZjo>

Have you ever wondered how it is possible to roll up the cord on a vacuum cleaner at the push of a button? The mechanism is called a scroll spring drive, where 2 spiral springs roll up and unroll against each other, moving 2 reels that wind up or unwind the power cable.

A spiral spring is a leaf spring that is rolled up in a spiral and thus very strongly bent. In this experiment, you will take a closer look at such a leaf spring.

Although it is not as strongly bent as a spiral spring, it still has similar properties.

## Equipment

| Position | Material   | Item No. | Quantity |
|----------|--|----------|----------|
| 1        | <a href="#">Cobra SMARTsense - Force and Acceleration, <math>\pm 50\text{N}</math> / <math>\pm 16\text{g}</math> (Bluetooth + USB)</a> | 12943-00 | 1        |
| 2        | <a href="#">Support base, variable</a>   | 02001-00 | 1        |
| 3        | <a href="#">Support rod, stainless steel, <math>l = 250\text{ mm}</math>, <math>d = 10\text{ mm}</math></a>                            | 02031-00 | 1        |
| 4        | <a href="#">Support rod, stainless steel, <math>l = 600\text{ mm}</math>, <math>d = 10\text{ mm}</math></a>                            | 02037-00 | 2        |
| 5        | <a href="#">Boss head</a>  | 02043-00 | 2        |
| 6        | <a href="#">Leaf spring</a>  | 02228-00 | 1        |
| 7        | <a href="#">Glass tube holder with tape measure clamp</a>  | 05961-00 | 1        |
| 8        | <a href="#">Measuring tape, <math>l = 2\text{ m}</math></a>  | 09936-00 | 1        |
| 9        | <a href="#">Fishing line, <math>l. 20\text{m}</math></a>   | 02089-00 | 1        |
| 10       | <a href="#">Support rod, stainless steel, <math>l = 250\text{ mm}</math>, <math>d = 10\text{ mm}</math></a>                            | 02031-00 | 1        |
| 11       | <a href="#">measureAPP - the free measurement software for all devices and operating systems</a>                                       | 14581-61 | 1        |

## Task

PHYWE



- Deflect a leaf spring from its resting position and measure the restoring force with a force gauge.
- Display the measurement results in a diagram.
- In doing so, contrast the restoring force and the deflection.

## Set-up (1/5)

PHYWE

For measurement with the **Cobra SMARTsense sensors** the **PHYWE measureAPP** is required. The app can be downloaded free of charge from the relevant app store (see below for QR codes). Before starting the app, please check that on your device (smartphone, tablet, desktop PC) **Bluetooth** is **activated**.



iOS



Android



Windows

## Set-up (2/5)

PHYWE

Connect the separate tripod rods to form two long tripod rods, with a length of 600 mm each.

Attach the two tripod base halves to the two ends of the long tripod rods.

By raising the locking levers, the respective tripod rod is fixated in the tripod base.



Tripod rods with thread



Assembly of the bases



Fixing the stand rod

## Set-up (3/5)

PHYWE

Now insert the short tripod rod vertically into a tripod base and screw it tight.

Attach a double socket to the center of one of the horizontal tripod rods.


Then attach the tape measure with the glass tube holder to the second horizontal rod.



Assembly of the base



Measuring tape in glass tube holder




Fasten tape measure



## Set-up (4/5)

PHYWE




Experiment set-up


- Complete the experimental set-up as shown in the accompanying diagram.
- To do this, clamp the force sensor to the vertical tripod rod.
- Refer to the illustrations on the next page for securing the leaf spring in the double socket and for securing the fishing line to the leaf spring.

## Set-up (5/5)

PHYWE



Assembly of leaf spring in double socket




Fastening fishing line to leaf spring

- Do make sure that the leaf spring is correctly attached to the double socket (top illustration). The leaf spring is clamped vertically upwards.
- Using an approx. 5 cm-10 cm long loop, connect the leaf spring to the hook of the force sensor according to the bottom illustration.
- Then move the double socket with the leaf spring to the right until the loop is taut and the leaf spring is minimally bent by the loop in the direction of the force sensor.
- Adjust the tape measure so that the "0 cm end" is just at the left side of the double socket with the leaf spring!

## Procedure (1/2)

PHYWE




Example of measurement of the restoring force

- Turn on your Cobra SMARTsense Force Sensor. Open the measureAPP and select the force sensor as the sensor.
- When switched on, the force sensor is tared, i.e. at the start it shows a weight force of 0 N. Alternatively, the sensor must be tared manually in the start position. In this case, select "Set to zero".
- The measurement is performed as a point-by-point measurement. For this purpose, set "Measurement to keystroke".
- Now move the double socket with the leaf spring 1 cm to the right and start the measurement.

## Procedure (2/2)

PHYWE




Regression line

- The first measuring point thus corresponds to 1 cm of deflection. Before each new measurement, the spring is pushed further by 1 cm.
- Proceed accordingly until the total deflection is 15 cm. Then finish the measurement and save it. Make sure that the hook of the force sensor is approximately at the same height as the free end of the leaf spring. You may have to adjust the height of the force sensor accordingly during the measurement.
- Evaluate your data with the help of the questions in the protocol. As an observation, use your measurement on the tablet.

PHYWE

# Report



## Task 1

PHYWE

[www.giphy.com](http://www.giphy.com)

When a force is applied to the metal plate, it deforms. When the force is removed, the deformation decreases. This is characteristic of...

☐ an elastic and inelastic deformation☐ an inelastic deformation☐ an elastic deformation

## Task 2

PHYWE

Drag the words to the right places

For  deflections the course of the measured values deviates  
 from the linear course. This happens approximately from a deflection of  
.

Not needed:  (adjective),  (deflection).

large

7 cm

small

14 cm

increasingly


✓ Check

## Task 3

PHYWE

The spring constant  $D$  of a leaf spring is defined as the quotient of the force pulling perpendicularly to the end of the leaf spring  $F$  and the corresponding deflection  $\Delta s$

$$D := F / \Delta s$$



What has the spring constant  $D$  to do with the straight line gradient that was determined?

The spring constant  $D$  has nothing to do with the straight line gradient!

The spring constant  $D$  equals the gradient of the straight line!

The spring constant  $D$  equals the doubled straight line gradient!

| Slide  | Score / Total |
|--|---------------|
| Slide 19: Type of deformation  | 0/3           |
| Slide 20: Deviation of the measured values                             | 0/5           |
| Slide 21: Relationship of the straight line slope to the spring con... | 0/2           |

Total  0 / 10

 Solutions

 Repeat