

# Damping with Cobra SMARTsense



In this experiment, students are asked to observe the amplitudes of vibration over long periods of time and to know and measure their decrease with time.

Physics

Mechanics

Vibrations &amp; waves



Difficulty level

easy



Group size

2



Preparation time

10 minutes



Execution time

10 minutes


This content can also be found online at:



<http://localhost:1337/c/6123a85ce3e7380003d48737>

PHYWE

## Teacher information



## Application

PHYWE



### Damping with SMARTsense

So far, the students have conducted experiments on idealized oscillations in which influences due to friction, damping and air resistance are neglected. In fact, however, every oscillation experiences damping, which causes the amplitude of the oscillation to decrease over time.

Here the free damped oscillation can be described as a decreasing e-function as follows:

$$x(t) = x_0 \cdot e^{-\delta t} \cdot \cos(\omega t + \varphi_0)$$

With initial deflection  $x_0$  the angular frequency  $\omega$  and the phase shift  $\varphi_0$  at the beginning of the oscillation.

## Teacher information (1/2)

PHYWE

## Prior knowledge



Students should already be familiar with the helical spring pendulum and the acceleration due to gravity.  $g$  and its averaged value  $9,81 \text{ m/s}^2$  as this value plays a major role in the oscillation behaviour of each pendulum.

## Scientific Principle



The amplitude of a free, damped oscillation decreases with time. This phenomenon is sometimes due to the damping of the system. This results, for example, from the air resistance of the coil spring pendulum or from the damping in the material of the spring.

Note: The damped oscillation can be described with the help of an e-function. In this case, the e-function should have been treated mathematically in advance. It does not necessarily have to be introduced in the context of this experiment.

## Teacher information (2/2)

PHYWE

## Learning objective



The students should observe the amplitudes of vibration over longer times and learn about and measure their decrease with time. This leads to the concept of damping, which is to be introduced here primarily qualitatively and not as the commonly used logarithmic damping decrement.

## Tasks



1. The students let a thread pendulum swing freely in air and measure the deflection after different times. They also compare the deflection with the original deflection, then attach a disc to the thread pendulum and repeat the measurement.
2. The students immerse the mass of the pendulum in water and investigate the deflection after different times.

## Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons to be applied to this experiment.

PHYWE

## Student Information



## Motivation

PHYWE



### Damping with SMARTsense

#### Do vibrations occur in everyday life?

Whether when swinging on the playground, playing a musical instrument or generating clocks for electronic devices. Vibrations can be intentional (e.g. oscillating quartz in watches) or unintentional (e.g. vibrations when driving a car).



In this experiment you will learn to observe the amplitudes of oscillations over longer periods of time and to know and measure their decrease with time.

## Task

PHYWE



1. Let a thread pendulum oscillate freely in air and measure the deflection after different times. Compare the deflection with the original deflection. Then attach a disc to the thread pendulum and repeat the measurement.
2. Immerse the mass of the pendulum in water and examine the deflection after different times.



## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Cobra SMARTsense - Force and Acceleration, <math>\pm 50\text{N}</math> / <math>\pm 16\text{g}</math> (Bluetooth + USB)</a>	12943-00	1
2	<a href="#">Support base, variable</a>	02001-00	1
3	<a href="#">Support rod, l = 600 mm, d = 10 mm, split in 2 rods with screw threads</a>	02035-00	1
4	<a href="#">Boss head</a>	02043-00	1
5	<a href="#">Holding pin</a>	03949-00	1
6	<a href="#">Helical spring, 3 N/m</a>	02220-00	1
7	<a href="#">Weight holder, 10 g</a>	02204-00	1
8	<a href="#">Slotted weight, black, 10 g</a>	02205-01	4
9	<a href="#">Glass tube holder with tape measure clamp</a>	05961-00	1
10	<a href="#">Digital stopwatch, 24 h, 1/100 s and 1 s</a>	24025-00	1
11	<a href="#">Measuring tape, l = 2 m</a>	09936-00	1
12	<a href="#">Beaker, 250 ml, plastic (PP)</a>	36013-01	1
13	<a href="#">measureAPP - the free measurement software for all devices and operating systems</a>	14581-61	1

## Additional material

PHYWE

Position	Equipment	Quantity
1	Drawing cardboard	approx. DIN A4

## Set-up (1/3)

PHYWE

For measurement with the **Cobra SMARTsense sensors** the **PHYWE measureAPP** is required. The app can be downloaded free of charge from the relevant app store (see below for QR codes). Before starting the app, please check that on your device (smartphone, tablet, desktop PC) **Bluetooth** is **activated**.



iOS



Android




Windows

## Set-up (2/3)

**PHYWE**  
excellence in science

- Rotate the two-part tripod rod together (left figure).
- Assemble the tripod base and the tripod rod to a tripod (right picture).



## Set-up (3/3)

**PHYWE**

- Attach the force sensor in the double socket (left figure).
- Hang the coil spring on it (right figure).





## Procedure (1/8)

PHYWE



Switch on





Select sensor in measureAPP

- Turn on the force sensor by pressing the power button for several seconds.
- After successful switching on you will see a flashing LED (left picture).
- Start the measureAPP. Tap on the tab "Sensor" and select the force sensor (right figure).

## Procedure (2/8)

PHYWE  
excellence in science

- Tap the "Settings" tab and select "Set to zero" (left figure). In the following window, tap on the force sensor.
- Exit the window by clicking on save (right image).



## Procedure (3/8)

PHYWE



Execution - weight plate

- Hang the weight plate with a total mass of 150g on the eyelet of the coil spring.
- The spring should be completely at rest and not vibrate. Therefore, steady the system with your hand.

**Note** To attach the slotted weights to the weight plate, slide them over the top of the weight plate (see illustration).

## Implementation (4/8)

PHYWE



- Set the mass in oscillation. Make sure that the spring is only deflected vertically and does not move sideways.
- Make sure of this by observing the vibration for a few seconds.
- If the oscillation is not smooth, stop the spring and try again.


## Procedure (5/8)



Execution - Start measurement

- Start the measurement (figure).
- Stop the measurement after three minutes (180 seconds).
- Use the auto-zoom function (figure).
- Find the mean value of the curve and write it down (figure).

## Procedure (6/8)



Execution - Maximum amplitude

- Read off the maximum deflection of the force at the start of the measurement and note it down (Figure). To do this, zoom with your fingers into the upper left corner of the trace.
- Use the auto-zoom function again. Read off the maximum displacement at the end of the measurement and note it down. To do this, zoom with your fingers into the upper right corner of the trace.
- Save the measurement.


## Procedure (7/8)

PHYWE



- Attach a circular cardboard disc to the center of the weight plate and repeat the entire measurement.
- Place the cup filled with water underneath the spring pendulum and immerse the weight plate so that the weight is in the middle underneath the water surface.
- Deflect the spring pendulum downwards from the equilibrium position until the weight reaches the bottom of the vessel. Then release the mass and simultaneously start the measurement for a duration of 10 seconds.

## Procedure (8/8)

PHYWE  
excellence in science


Feedthrough - Tripod base

- To disassemble the tripod base, press the buttons in the middle and pull both halves apart.



PHYWE

# Report



## Task 1

PHYWE

Give for the case of **free** oscillation displays the determined values:

mean force	<input type="text"/>
Force at max	<input type="text"/>
Deflection	<input type="text"/>
(start)	<input type="text"/>
Force at max	<input type="text"/>

### Determination of mean value and deflection

The mean value of the force corresponds to the  of the spring pendulum. Via the relation  can therefore be used to determine the mass. The mean value must be divided by  $g$  .

The difference between the forces at maximum deflection and the mean value corresponds to the spring force . If one divides this difference by the spring constant  $D$  one obtains the geometric

.

## Task 2

PHYWE

Give for the oscillation case with **Cardboard disc** displays the determined values:

mean force

Force at max

Deflection

(start)

Force at max

Calculate the geometric deflections of the spring pendulum at the beginning and end of the measurement. Indicate by what percentage the amplitude has decreased and compare with the previous result. What stands out?

In task 1, the amplitude has decreased by . In task 2, it has decreased by about . This value is  than before. The reason for this is that the resistance in task 1 is  than the resistance caused by the cardboard disc.

☒ Check

## Task 3

PHYWE

Give for the vibration case **underwater** displays the determined values:

mean force

Force at max

Deflection

(start)

Force at max

Calculate the geometric deflections in the same way as before. What stands out? What does the limit value 0 mean for the amplitude of the oscillation?

The amplitude dropped significantly  than the drops in air. The water provides a  resistance which slows the pendulum down quickly (only 10 seconds measurement). The amplitude of the oscillation has dropped , i.e. the oscillation has come to a .

☒ Check

Slide	Score / Total
Slide 23: Mean value and deflection	0/5
Slide 24: Geometric deflections of the spring pendulum	0/4
Slide 25: Geometric deflections of the spring pendulum	0/4

Total  0 / 13

[Solutions](#)[Repeat](#)[Export text](#)