


Instantaneous and average speed with timer 2-1



Physics

Mechanics

Dynamics & Motion



Difficulty level

medium



Group size

2



Preparation time

10 minutes



Execution time

20 minutes


This content can also be found online at:



<http://localhost:1337/c/5fdbb8485098f00003f1ede2>


PHYWE

Teacher information



Application

PHYWE



Experiment set-up

By using light barriers, both the instantaneous and the average speed can be determined.

This technology is used, among other things, in traffic monitoring. Here, for example, the speed of road users can be determined using the difference in travel time of transmitted radar waves or light pulses. For correct scaling, there are often standardized markings on the road.

Other teacher information (1/2)

PHYWE

Prior knowledge



Students should know how a photoelectric sensor works.

Scientific principle



A body experiences a constant acceleration in an inclined path due to the component of gravity acting on it parallel to the path. Therefore, the laws of motion apply to a uniformly accelerated motion.

Other teacher information (2/2)

PHYWE

Learning objective



In this experiment, students will quantitatively investigate the differences between uniform and non-uniform motions. Specifically, the average velocity $v = \Delta s / \Delta t$ should be differentiated from the momentary speed $v = \dot{s}$.

Tasks



1. Determine the average speed: measure the time taken by the car to cover a given distance using two light barriers at the beginning and end of the respective distance.
2. Determine the instantaneous velocity: Measure the time it takes for the aperture on the car to pass the light barrier after such a distance.

Safety instructions


PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student Information



Motivation

PHYWE



Measuring device for radar control

Since everyday movements are particularly non-uniform, it is important to distinguish between the so-called instantaneous and average speeds. In road traffic, the instantaneous speeds are generally determined with radar traps as average speeds for very short periods of time. The longer this time span is, the more the average speed can deviate from the intermediate instantaneous speeds.

In this lesson, you will learn the difference between these two quantities by using two photoelectric sensors in different measurement modes.

Tasks

PHYWE



1. Measure the time taken by the cart to travel a certain distance using two light barriers at the beginning and end of the respective distance. Calculate the average speed from the time measured between the interruption of one and the other light barrier and the length of the track.
2. Measure the time required for the diaphragm on the cart to pass the light barrier after such a distance. Calculate approximately the instantaneous velocity with this shading time of the light barrier and the aperture width.


Equipment

Position	Material	Item No.	Quantity
1	Cart for measurements and experiments	11060-00	1
2	Shutter plate for cart	11060-10	1
3	Holding pin	03949-00	1
4	Slotted weight, black, 50 g	02206-01	3
5	PHYWE Timer 2-1	13607-99	1
6	Light barrier, compact	11207-20	2
7	Adapter plate for Light barrier compact	11207-22	2
8	Connecting cord, 32 A, 1000 mm, red	07363-01	2
9	Connecting cord, 32 A, 1000 mm, yellow	07363-02	2
10	Connecting cord, 32 A, 1000 mm, blue	07363-04	2
11	Track, I 900 mm	11606-00	1

Set-up (1/3)

PHYWE

To tilt the carriage, screw the adjustable foot of the carriage all the way down and place it on two stacked 50 g slotted weights. Then attach the shade screen to the cart using the holding pin and weigh it down with a 50 g slotted weight.



Set-up (2/3)


PHYWE

Screw the spacer bolts and adapter plates to the forked light barriers in such a way that they can be easily positioned next to the carriageway and the aperture on the cart fits through the respective light barrier without hitting it. Connect both light barriers to the timing device. Set the slide switch above the field labelled "Start" on the timing device to the right position.



Set-up (3/3)

PHYWE




Direction of travel of the cart

Position the cart so that the end of the cart is flush with the end of the track. Position the first light barrier so that the screen on the cart interrupts it as soon as the cart is released.

Position the second light barrier at a distance of 20 cm from the first one.

Procedure (1/2)

PHYWE




Direction of travel of the cart

- Set the rotary switch on the timing device to the third position from the left. The device then displays the time that elapses between the interruption of the first and the second light barrier. In this experiment, this is the time that the cart takes to pass through the distance Δs between the two light barriers.
- Press the "Reset"-button on the timing device before each measurement.
- Let go of the car without nudging it and note the time t required for the distance by the car $\Delta s = 20\text{ cm}$ in Table 1 in the report. Repeat the measurement for the distances $\Delta s = 30, 50, 70\text{ cm}$.

Procedure (2/2)

PHYWE



Removing the first light barrier

- Now remove the first light barrier from the track far enough away so that it is not interrupted by the aperture of the car.
- Set the rotary switch to the second position from the left. The device then displays the shading time. This is the time period in which a light barrier is interrupted by an aperture.
- Repeat the measurement for all positions of the second light barrier from the first part of the experiment. Start each measurement with the "Reset" button and let the car roll down the track. Enter the measured times t also in Table 1 in the report.

PHYWE

Report




Table 1

PHYWE

Enter the calculated values for the travel times Δt and the shading times t for the respective distances Δs into the table. Then calculate the average speed $v_d = \Delta s / \Delta t$ and the instantaneous speed $v_m = b / t$ with the aperture width $b = 5 \text{ cm}$.


Distance Δs [cm]	Δt [s]	v_d [cm/s]	t [s]	v_m [cm/s]
20				
30				
50				
70				

Task 1

PHYWE


Now take a sheet of paper and draw a diagram on it. In this diagram you place the two velocities v_d and v_m (y -axis) as a function of the distance travelled Δs (x -axis).

Draw both curves in a diagram.



Task 2

PHYWE



Cart

How do the speeds compare v_m to each other?

- ☐ The speeds are all the same.
- ☐ The speeds increase as the length of the route increases.
- ☐ The speeds drop with increasing distance.

☒ Check

Task 3

PHYWE

Drag the words to the right places.

The v_m is at the end of the line s is always
 than the v_d over the same
distance. Thus, the speed over the course of the
distance.
not needed: (adjective),
(noun).

increases

maximum speed

greater

instantaneous speed

average speed

smaller

☒ Check

Task 4

PHYWE



Cart

Can we speak of a uniform movement here?

- ☐ Yes, there is uniform motion because the car is moving at a constant speed.
- ☐ No, there is no uniform movement, because the car is continuously accelerating.

[✓ Check](#)

Slide

Score/Total

Slide 18: Comparison of $\backslash(v_m\backslash)$	0/1
Slide 19: Shapes of the velocities	0/6
Slide 20: Nature of the movement	0/1

Total  0/8[Solutions](#)[Repeat](#)[Export text](#)