


# Determining the magnification of a concave lens



Physics

Light &amp; Optics

Optical devices &amp; lenses



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/62dbce3da52f910003dffb99>


PHYWE



## Teacher information

## Application

PHYWE



Experimental setup

Convex lenses, also called converging lenses, can produce a magnified image. They are an important element of ray optics and are therefore widely used in optical instruments and in camera lenses.

## Other teacher information (1/4)

PHYWE

## Principle



Incident light that is parallel to the optical axis is focused by the convex lens at the focal point. This can produce a magnified real image.

## Learning objective



Students should observe the optical effect of a convex lens and calculate the magnification of the lens  $b/g = B/G$ .

## Other teacher information (2/4)

PHYWE

## Task



- Students should investigate the relationship between the width of the object  $g$ , the image width  $b$ , the object size  $G$  and the image size  $B$  in the case of images through convex lenses.
- For this purpose, for different object widths  $g$  the other variables are measured and noted in the table in the report.

## Other teacher information (3/4)



Obtaining the equation for the magnification by theoretical means generally presents few difficulties. It is therefore advisable to use the experiment as a confirmation experiment. In this way, the students have sufficient goal orientation and know from the outset that the important thing is to find the two quotients  $B/G$  and  $b/g$  by calculation and to compare them after the four occurring variables have been determined experimentally.

## Other teacher information (4/4)

PHYWE

### Notes on set-up and procedure

- With the available scale for the stand bench, it is possible to measure all four sizes. However, it is easier to measure if the scale is applied to the front stand rod and only the distances  $g$  and  $b$  are measured with it. Item size  $G$  and image size  $B$  can be measured more conveniently with an additional ruler.
- The choice of convex lens  $f = +100 \text{ mm}$  ensures that the entire length of the optical bench can be used for image generation without producing images that are too large or too small. In addition, with this lens you also get magnified images that are still very sharp, so that quite accurate measurements become possible. When using the lens  $f = 5 \text{ mm}$ , this is not the case, as lens aberrations have a very strong effect.

## Safety instructions


PHYWE



- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

## Student information



## Motivation

PHYWE




Spectacle lenses as an example of a convex lens

Convex lenses, also called converging lenses, can produce a magnified image. They are a central component of optics and are therefore often found in everyday devices, such as telescopes, camera lenses or even eyeglasses.

## Tasks

PHYWE



Experimental setup

- Investigate the relationship between the width of the object  $g$ , the image width  $b$ , the object size  $G$  and the image size  $B$  in the case of images through convex lenses.
- For this purpose, you should measure the other sizes for different object widths  $g$  and note them in the table in the report.

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Optical profile-bench for student experiments, l = 600 mm</a>	08376-00	1
2	<a href="#">Light box, halogen 12V/20 W</a>	09801-00	1
3	<a href="#">Bottom with stem for light box</a>	09802-20	1
4	<a href="#">Lens on slide mount, f=+100mm</a>	09820-02	1
5	<a href="#">Slide mount for optical bench</a>	09822-00	1
6	<a href="#">Screen, white, 150x150 mm</a>	09826-00	1
7	<a href="#">Object -L-, glass bead</a>	11609-00	1
8	<a href="#">PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A</a>	13506-93	1

## Set-up (1/3)

PHYWE

- Assemble the optical bench from the two tripod rods and the variable tripod foot and place the scale on the front tripod rod.
- Place the base with stem under the light box.



## Set-up (2/3)

PHYWE


- Clamp the light box in the left part of the tripod base so that the lens side faces away from the optical bench.
- Slide an opaque shade in front of the lens and the Perl-L into the shaft at the other end of the luminaire.





## Set-up (3/3)

PHYWE




Experimental setup

- Place the lens and, with the help of the slide, the screen on the optical bench.

## Procedure (1/3)

PHYWE






- Connect the lamp to the power supply unit (12 V~) and switch it on.
- Place the lens at a distance of 150 mm from the Perl-L and move the screen until the image of the Perl-L is sharp.

## Procedure (2/3)


PHYWE

- Measure the distance of the lens from the screen, the image width  $b$  and record for an object width  $g = 150 \text{ mm}$  in Table 1 in the report.
- Measure the item size  $G$  and the image size  $B$ .  $G$  and  $B$  should be the distances between the centres of the top and bottom beads or their images.



## Procedure (3/3)


PHYWE



Displacement of the lens

- Create another enlarged and two reduced sharp images of the Perl-L by moving the lens to the right.
- Measure the object width for each of these pictures  $g$ , the image width  $b$ , the object size  $G$  and the image size  $B$ . Enter the measured values in Table 1.
- Switch off the power supply unit.

PHYWE



Report

Table 1

PHYWE

Enter your measurements in the table.  
 $g$  in mm  $b$  in mm  $G$  in mm  $B$  in mm  $g/b$   $G/B$


## Task 1

PHYWE

Compare the quotients in the individual rows of Table 1. What do you find?

- ☐ The quotients  $b/g$  and  $B/G$  are (almost) the same in the individual lines.
- ☐ The quotients  $b/g$  are greater in each row than the quotients  $B/G$ .
- ☐ The quotients  $b/g$  are smaller in each row than the quotients  $B/G$ .

☒ Check

## Task 2

PHYWE

What is the result of your considerations in mathematical form?

- ☐  $g/b > G/B$
- ☐  $g/b = G/B$
- ☐  $g/b < G/B$

☒ Check

## Task 3

PHYWE

With which devices do you know that the magnification is important?

☐ Camera☐ Windows☐ TV screen☐ Writing projectors☐ Slide projectors☒ Check

Slide

Score / Total

Slide 20: Comparison of the quotients

0/1

Slide 21: mathematical relation

0/1

Slide 22: Applications

0/3

Total

 0/5 Solutions Repeat Export text