

Charge distribution in a Faraday cup



Physics

Electricity & Magnetism

Electrostatics & electric field



Difficulty level

easy



Group size

-



Preparation time

10 minutes



Execution time

10 minutes


This content can also be found online at:



<http://localhost:1337/c/6426c878ab58420002f62a76>

PHYWE

Teacher information



Application

PHYWE



Faraday cup

A Faraday cup is - according to its name - a cup-shaped object with the help of which one can determine electric charge (especially ion beams).

The inside of the metal cup is basically field-free and electrical charges can only be taken out via the outer wall of the cup. However, the electric charge can also be supplied to the cup via the inner wall.

This therefore means that no charge can migrate from the inner wall of the Faraday cup to an influence plate, but it is possible from the outer wall of the cup.

Other teacher information (1/2)

PHYWE

Prior knowledge



The students should already have studied electric charge and its effects in detail. The previous experiments, in which the electric charge is examined using a simple electroscope, provide basic knowledge for this. In addition, the experiment "Conductors as charge storage" provides a good basis for successfully carrying out this experiment.

Principle



The inside of the Faraday cup is basically field-free. Thus, no electrical charges can migrate from the inner wall of the cup to other objects. The outer wall of the cup, on the other hand, is able to transfer charges to other bodies.

Other teacher information (2/2)

PHYWE

Learning objective



The students realise that in a charged Faraday cup, the charges are on the outside of the cup and can only be removed from there. However, charges can be applied to both the outside and the inside of the cup.

Tasks



In this experiment, the students are to investigate the properties of the Faraday cup in detail. To do this, they have to determine the charge distribution on the Faraday cup.

Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.


Notes on set-up and procedure:

Experimentally critical is the proof that there are no charges on the inside of the Faraday cup. When using the influence plate, this can only be proven if the influence plate is safely inserted into the beaker without contact. Charge transfer can already occur when approaching the upper rim. Therefore, the plate must be handled with care.

Since the Faraday cup, in addition to its function as a charge storage device, also plays a major role in a modified form as a Faraday cage, the shielding of electrostatic fields can be investigated in the additional task. The protective effect of a car body during a thunderstorm can be pointed out.

PHYWE

Student information



Motivation

PHYWE



Faraday cup

A Faraday cup is - according to its name - a cup-shaped object with the help of which one can determine electric charge (especially ion beams).

Due to the special geometry of the Faraday cup, the inside of the metal cup is basically field-free, which has special effects on the charge distribution in relation to the cup.

This experiment shows you exactly how the charges are distributed on the Faraday cup and what properties result from this.

Tasks

PHYWE



In this experiment, you will take a close look at the so-called Faraday cup and its properties.

First set up an electroscope and make sure that the Faraday cup is insulated as well as possible.

Then investigate the distribution of the electrical charge on a metal cup.

Equipment

Position	Material	Item No.	Quantity
1	Electroscope w. metal pointer	13027-01	1
2	Faraday pail, d. 40mm, h. 75mm	13027-03	1
3	Polypropylene rod, l=175mm, d=10 mm	13027-09	1
4	Acrylic resin rod, l=175 mm, d=8 mm	13027-08	1
5	Electrostatic influence plate, 30 x 60 mm	13027-12	1
6	Film, transparent, DIN A4, 100 sheets	08186-10	1
7	Rubber stopper,d=49/41mm, 1 hole	39263-01	1

Additional Equipment

PHYWE

Position	Equipment	Quantity
1	Dry, rough paper	DIN A4
1	Adhesive tape	

Set-up

PHYWE




Assemble the electroscope. The pointer should hang vertically (one side is slightly longer and therefore minimally heavier) without bumping, the axis is in the notch.

Then attach the influence plate to the polypropylene rod and the Faraday cup to the acrylic rod. Then insert the acrylic rod into the rubber stopper as shown in the illustration opposite.



Procedure (1/6)

PHYWE



Charging the Influence Plate

Experiment 1:

- Rub the influence plate over the transparent film lying on the table.

Procedure (2/6)

PHYWE
excellence in science

Experiment 1:

- Then touch the outside of the Faraday cup with the Influence Plate.
- Repeat this process several times so that the Faraday cup is well charged.
- Now unload the influence plate and the electroscope by touching them with your hand.



Procedure (3/6)

PHYWE



Alternately touching the Faraday cup and the electroscope

Experiment 1:

- Now alternately touch the outside of the Faraday cup and the electroscope with the Influence Plate
- Keep a close eye on the pointer.

Procedure (4/6)

PHYWE



Examine the inside of the cup

Experiment 2:

- Check whether charges can also be transferred to the electroscope from the inside of the Faraday cup.
- As in the first part of the experiment, charge the Faraday beaker well and discharge the electroscope.
- Touch the bottom or the inner wall of the beaker (but not the upper rim) with the influence plate and then the electroscope.
- Finally, touch the outside of the cup again and then touch the electroscope to check.

Procedure (5/6)

PHYWE




Discharge Faraday Cup



Electroscope discharged

Experiment 3:

- Now check whether charges can be transferred to the inside of the Faraday cup.
- First discharge the Faraday beaker and the electroscope.
- Then charge the Influence Plate by rubbing it on the transparent film.



Charging the Influence Plate on Transparent Film

Procedure (6/6)

PHYWE



Examine the inside of the cup

Experiment 3:


- Touch the inside of the cup with the charged Influence Plate.
- Be careful not to touch the edge or the outside.
- Then check with the influence plate and the electroscope whether charges can be picked up from the inside of the cup.
- Then also check whether charges can be taken from the outside.



Examine the outside of the cup

PHYWE

Report



Task 1

PHYWE



Alternately touching the Faraday cup and the electrostatic demonstrator

How does the position of the pointer change with the individual touches of the electrostatic demonstrator in the first partial experiment?

With each further movement, the pointer deflected further. Strongly at first, then less and less.

With each movement, the pointer moved briefly and then returned to its original position.

Task 2

PHYWE



Alternately touching the Faraday cup and the electroscope

What does the changing pointer deflection in the first partial experiment tell us about the size of the charges on the Faraday cup and the electroscope?

The the increase in the pointer deflection, the more electrical charge was removed from the cup.

The the increase in the pointer deflection, the less residual charge was in the cup.

Task 3

PHYWE



Examine the inside of the cup

When did a pointer deflection occur in the second part of the experiment with the electroscope?

Both after touching the inside and the outside of the cup, a deflection occurred on the electroscope.

While no pointer deflection occurs after previous contact with the inside, the pointer deflects after contact with the outside wall.

Task 4

PHYWE



Examine the inside of the cup

Where on the Faraday cup were the transferred charges located during the second partial experiment?

The charge in the Faraday cup was on the of the cup, so it cause the electroscope to move.

The of the cup was free of charge, so it cause the electroscope to move.

☒ Check

Task 5

PHYWE



Examine the inside of the cup

Where could charges be taken from in the third partial experiment?



Examine the outside of the cup

Task 6

PHYWE

What is the difference between the possibilities of applying and removing charges on the Faraday cup?

The Faraday cup is made of metal and is therefore electrically .
Therefore, it can be charged at point during contact, especially on the . Due to the force effect of charges of the same name, the charges are distributed at the distance from each other, i.e. especially on the . Accordingly, the charges can only be reduced there.

inside

outside

conductive

any

repelling

largest possible

 Check

Slide	Score / Total
Slide 19: Observation: Experiment 1	0/1
Slide 20: Conclusion: Experiment 1	0/2
Slide 21: Explanation: Experiment 2	0/1
Slide 22: Explanation: Experiment 2	0/4
Slide 23: Observation: Experiment 3	0/1
Slide 24: Conclusion	0/6

Total  ★ 0 / 15 Solutions Repeat