

Discharging at points



Physics

Electricity & Magnetism

Electrostatics & electric field



Difficulty level

easy



Group size

-



Preparation time

10 minutes



Execution time


10 minutes

This content can also be found online at:



<http://localhost:1337/c/64280d855e30a7000275e9de>

PHYWE



Teacher information

Application

PHYWE



Lightning in a thunderstorm

For the exchange of electrical charge between two bodies to take place, they do not necessarily have to come into contact. It is even sufficient if they approach each other at a certain distance (in this case a few centimetres).

This transition of the electrical charge then usually becomes visible in the form of a spark.

A classic example of the visible exchange of electrical charges over long distances is lightning in a thunderstorm. Here, the transition of the charge from the cloud to the ground is very clearly visible.

Other teacher information (1/2)

PHYWE

Prior knowledge



Students should already have thoroughly studied and understood electric charge and its effects. In addition, they should have understood that electric charge can be transferred between two bodies by touching them.

Principle



If the electric charge of a body is sufficiently large, it can also be transferred to other bodies that are close to it without touching it!

Other teacher information (2/2)

PHYWE

Learning objective



The students should realise that charges can be transported through the air with the help of sparks and that these sparks mainly occur at the tips and corners of objects.

Tasks



In this experiment, students are asked to discharge an electrically charged electroscope without touching it.

They are to show that charge can also be transported without contact between two bodies and that the shape of the bodies can play a role in this.

Safety instructions

PHYWE




The general instructions for safe experimentation in science lessons apply to this experiment.

Notes on set-up and procedure:

This experiment is not mainly about the spark itself, but about the transfer of charges between two bodies at points. Point discharge on bodies can only be shown in demonstration experiments, as higher voltages are needed. The comparison of distances in the first sub-experiment is limited by the fact that the electroscope would have to carry the same amount of charge each time. However, this is difficult to realise. For the series of events in the third part of the experiment, the spark is already easily visible in moderately darkened rooms. An explanation of the mechanism of conduction through air can be omitted if the students' prior knowledge is not sufficient for this.

PHYWE



Student information

Motivation

PHYWE



Lightning in a thunderstorm

Electrically charged objects can transfer this charge when in direct contact with conductors.

But how does the exchange of charge work in the event of a thunderstorm between the cloud and the ground? In this case, the discharge is very clearly visible as a flash of lightning, when the excess charges are transferred from the cloud to the ground. However, the principle also works in the same way on a smaller scale, for example with the Tesla transformer or with plasma spheres.

You will investigate the phenomenon of contactless charge transfer more closely in the next experiment.

Tasks

PHYWE



In this experiment you will again deal with the discharge of different bodies.

To do this, first set up an electroscope.

Then examine the electrically charged electroscope without touching it!

Equipment

Position	Material	Item No.	Quantity
1	Electroscope w. metal pointer	13027-01	1
2	Faraday pail, d. 40mm, h. 75mm	13027-03	1
3	Polypropylene rod, l=175mm, d=10 mm	13027-09	1
4	Acrylic resin rod, l=175 mm, d=8 mm	13027-08	1

Additional Equipment

PHYWE

Position	Equipment	Quantity
1	Dry, rough paper	DIN A4

Set-up (1/2)

PHYWE



Structure of the electroscope

First assemble the electroscope.

The pointer should hang vertically (one side is slightly longer and thus minimally heavier) without bumping, the axis is in the notch.

Set-up (2/2)

PHYWE



Attach Faraday cup to acrylic rod

Then mount the Faraday cup at the end of the acrylic rod.

Procedure (1/5)

PHYWE

- Rub the polypropylene rod vigorously with paper to charge it.
- Then charge the electroscope with the polypropylene rod rubbed against the paper. Repeat the rubbing and charging a few times to fully charge the electroscope.



Procedure (2/5)

PHYWE



Bring your finger close to one corner of the electroscope (without touching it!)

Experiment 1, Part 1:

- Slowly move a finger towards a corner on the upper arm of the electroscope.
- Note the distance between the finger and the corner of the electroscope when a small crack is heard.
- Also observe the behaviour of the pointer.

Procedure (3/5)

PHYWE



Bring your finger close to the upper surface of the electroscope (without touching it!).

Experiment 1, part 2:

- Repeat the experiment, but this time move your finger towards a surface on the electroscope.
- Again, note the distance when a crack is heard.

Procedure (4/5)

PHYWE



Approach Faraday cup to electroscope

Experiment 2:

- Recharge the electroscope as much as possible.
- Move the Faraday cup on the acrylic rod to one corner of the electroscope until you hear a crack, but do not touch the electroscope.
- Now devise an experiment to find out whether charges were transferred from the electroscope to the Faraday cup.
- Carry out this experiment.

Procedure (5/5)

PHYWE



Charge mug




Approach cup to electrohead

Experiment 3:

- Charge the Faraday cup strongly on the acrylic rod with the polypropylene rod rubbed on the paper.
- Then move it slowly to a corner of the electroscope, which should be discharged beforehand.
- Listen carefully to the crackling and watch the hand.

PHYWE



Report

Task 1

PHYWE



Bring your finger close to one corner of the electrostatic demonstrator (without touching it!)

What were your observations of part 1 of the 1st experiment?

The pointer did not move at any time during the experiment.

The pointer was deflected even further after cracking.

The pointer has clearly moved back towards the starting position after the crackling.

Task 2

PHYWE



Bring your finger close to the upper surface of the electroscope (without touching it!).

What were your observations of part 2 of the 1st experiment?

The pointer was deflected even further after cracking.

The pointer has clearly moved back towards the starting position after the crackling.

The pointer did not move at any time during the experiment.

Task 3

PHYWE



Experiment 1, Part 1



Experiment 1, Part 2

In which part of experiment 1 was the distance between the finger and the electroscope greater than the crackling sound?

In the 1st part, the distance was greater (finger approaching the corner of the electroscope).

In the 2nd part, the distance was greater (fingers approaching the surface of the electroscope).

The distance was the same in both sub-tests.

Task 4

PHYWE



Approach Faraday cup to electroscope

What were your observations of the second experiment?

The pointer did not move at any time during the experiment.

The pointer was deflected even further after cracking.

The pointer has clearly moved back towards the starting position after the crackling.

Task 5

PHYWE



Approach Faraday cup to electroscope

What were your observations of the third experiment?

The pointer did not move at any time during the experiment.

The hand briefly left its initial position, but then returned to it.

A repeated crackling sound is heard on approach, and the pointer deflection increases with each crack.

Task 6

PHYWE



Experiment 1, Part 1



Experiment 1, Part 2

What influence does the shape of the surface have on the distance at which the observed phenomenon occurred?

The the shape, the easier it is to transfer charge.

The the shape, the closer you have to get to the body for the charge to be transferred.



Check

Task 7

PHYWE



Faraday cup attached to electroscope

How can the spark be made more visible or the distance be increased?

☐ Approach oppositely charged object☐ Increase charge quantity (through Faraday cup)☐ Darken room☐ Approaching a conductive object at the tip

Check

Task 8

PHYWE

In the individual experiments, the bodies were always negatively charged. Think about whether you would observe or obtain the same phenomena and results with positive charges?

No, the observations and results would only be partially the same.

No, the observations and results would always be opposite.

Yes, the observations and results would be the same.

Slide	Score/Total
Slide 19: Observation: Experiment 1, Part 1	0/1
Slide 20: Observation: Experiment 1, Part 2	0/1
Slide 21: Comparison Experiment 1, Part 1 & 2	0/1
Slide 22: Observation: Experiment 2	0/1
Slide 23: Observation: Experiment 3	0/1
Slide 24: Conclusion 2	0/2
Slide 25: Sparks clarify	0/4
Slide 26: Positive charges	0/1

Total  0/12

 Solutions

 Repeat