


Shadow formation by a point light source



Physics

Light & Optics

Dispersion of light



Difficulty level

easy



Group size

-



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/642814775e30a7000275e9ee>


PHYWE



Teacher information

Application

PHYWE



Experimental set-up:

Adhesive luminaire with point-shaped light source, divergent light beam and shadow body

Light spreads out in a straight line. If a ray of light hits an opaque object, a shadow is created.

The straight-line spread of light, from a point light source, creates a sharp shadow that reproduces the contour of the shadow-casting object on a scale.

Other teacher information (1/2)

PHYWE

Prior knowledge



Students need prior theoretical knowledge about the straight-line, ray-shaped propagation of light.

Principle



It is to be shown that a shadow is created behind an opaque body, the area of which depends on the distance of the body from the light source.

Other teacher information (2/2)

PHYWE

Learning objective



Students should gain knowledge about the principles of shadowing.

Furthermore, they should recognise that the sharp formation of shadows is the result of the straight-line propagation of light from a point-shaped light source.

Tasks



The students should collect observations and insights into the effect of sharp shadows.

Safety instructions

PHYWE



- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Motivation

PHYWE



Shadow play of a group of people

When the sun shines, our bodies and objects in our surroundings cast a shadow.

But how is such a shadow created, and how is the sharp outline formed?


This is the question this experiment aims to clarify.

Equipment

PHYWE

Set-up and Procedure

PHYWE




Adhesive luminaire with point-shaped light source,
divergent light beam and shadow body

- Place the adhesive luminaire on the adhesive panel, push up the trim of the rear wall of the luminaire and create a divergent light beam with the point-shaped light source
- Bring shadow body earth into the light beam so that it still partially passes the shadow body on both sides; observe shadow
- Place the shadow body towards and away from the light source while keeping an eye on the shadow

PHYWE

Report



Task 1

PHYWE

Fill in the blanks

Behind the body, an area remains unlit; a is created. The shadow area is the larger, the the distance between the body and the light source. The shadow is limited.

☒ Check

Task 2

PHYWE

Answer the questions below to check that you have understood the experiment correctly.

Tip: Although our sun is very far away, our shadow is often almost as big as we are.

Behind an opaque body, an area remains unlit!

☐ True☐ False☒ Check

The smaller the distance between the body and the light source, the smaller the shadow area!

☐ True☐ False☒ Check

Slide	Score / Total
Slide 11: Shadow area	0/4
Slide 12: Multiple tasks	0/2

Total   0/6



Solutions



Repeat