

Aberrations with a concave mirror (catacaustics)



Physics

Light & Optics

Reflection & refraction of light



Difficulty level

easy



Group size

-



Preparation time

10 minutes




Execution time

10 minutes

This content can also be found online at:

<http://localhost:1337/c/642827335e30a7000275ea16>


PHYWE



Teacher information

Application

PHYWE



Experimental set-up:

Concave mirror

Light propagates in a straight line. If a light beam hits a reflecting object, the light beam also spreads out in a straight line from there.

A concave mirror can only sharply reproduce the incoming rays to a certain extent. If the rays are too far away from the optical axis, they can no longer be perceived.

The experiment is intended to show that it is only possible to work properly with such a mirror at a certain distance.

Other teacher information (1/2)

Prior knowledge



The students need theoretical knowledge about the straight-line, ray-shaped propagation of light and that objects reflect light rays.

Furthermore, students should know the basics about a concave mirror and its reflection behaviour.

Principle



It is to be shown that aberrations occur at the concave mirror if the light rays do not hit the mirror close to the axis.

Other teacher information (1/2)

PHYWE

Prior knowledge



The students need theoretical knowledge about the straight-line, ray-shaped propagation of light and that objects reflect light rays.

Furthermore, students should know the basics about a concave mirror and its reflection behaviour.

Principle



It is to be shown that aberrations occur at the concave mirror if the light rays do not hit the mirror close to the axis.

Other teacher information (2/2)

PHYWE

Learning objective



The students should gain knowledge about the principles of light reflection. Furthermore, they should realise that a concave mirror does not function as desired from every position.

Tasks



The students should observe how a concave mirror reflects the different incident rays. And understand that images on a concave mirror appear blurred or distorted when it is used disproportionately.

Additional teacher information

PHYWE

Note



If the concept of catacaustics is to be dealt with, then one can sketch the envelope curve for the reflected rays. This envelope curve can be obtained even more elegantly if the adhesive luminaire is used last without a diaphragm and the wide parallel light beam is shifted up and down parallel to the optical axis.

Safety instructions


PHYWE



- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Equipment

Position	Material	Item No.	Quantity
1	PHYWE Demo Physics board with stand	02150-00	1
2	Halogen lamp for experiments, 12V/50W, with magnetic base	08270-20	1
3	Diaphragm w. holder, magnet held	08270-10	2
4	Concave/convex mirror, magnet held	08270-12	1
5	PHYWE Multitap transformer DC: 2/4/6/8/10/12 V, 5 A / AC: 2/4/6/8/10/12/14 V, 5 A	13533-93	1
6	G-clamp	02014-00	2

Equipment

PHYWE

Position	Material	Item No.	Quantity
1	PHYWE Demo Physics board with stand	02150-00	1
2	Halogen lamp for experiments, 12V/50W, with magnetic base	08270-20	1
3	Diaphragm w. holder, magnet held	08270-10	2
4	Concave/convex mirror, magnet held	08270-12	1
5	PHYWE Multitap transformer DC: 2/4/6/8/10/12 V, 5 A / AC: 2/4/6/8/10/12/14 V, 5 A	13533-93	1
6	G-clamp	02014-00	2

Set-up and Procedure (1/3)

PHYWE




Fig.1:

Concave mirror with 3-slit aperture

- Draw optical axis
- Draw an arc of a circle with radius $r = 200 \text{ mm}$ on the board using a template or compass.
- Set mirror on circular arc
- Position the luminaire with the 3-slit aperture so that the central light beam runs along the optical axis; readjust the mirror if necessary; mark the focal point.

Set-up and Procedure (2/3)

PHYWE




Fig.2:

Concave mirror with 5-slit aperture

- Replace 3-slit aperture with 5-slit aperture
- Use apertures with a holder to alternately block the two parallel beams further away from the axis or the parallel beams closer to the axis.
- The illustration shows the course of the light rays further away from the axis

Set-up and Procedure (3/3)

PHYWE





Fig.3:

Concave mirror with 5-slit aperture

- Shift the luminaire downwards and upwards so that the parallel beams are partly very far away from the axis.
- Observe reflected rays and trace some as far as possible if necessary
- Remove mirror and luminaire and complete beam paths

PHYWE

Report




Task 1

PHYWE



Write down your observations on the different experimental set-ups.
Complete the ray traces of the third set-up by drawing.



Task 2

PHYWE



Fill in the blanks!

Only parallel rays [] the axis are reflected by the concave mirror in such a way that they intersect at a point on the [] axis, the [], after reflection at the concave mirror. Parallel rays that are far away from the axis run after [] at the concave mirror in such a way that they intersect the optical axis [] to the apex of the mirror the further they are from the optical axis.

near

optical

closer

reflection

focal point

☒ Check

Task 3

PHYWE



Fill in the blanks!

Images on a concave mirror become [] or distorted if the object is relatively [] in relation to the [] of the mirror. For then [] rays emanating from object points are partly [] rays.

blurred

radius of curvature

axis-distant

axis-parallel

large

☒ Check