


Reflection of light by a convex mirror



Physics

Light & Optics

Reflection & refraction of light



Difficulty level

easy



Group size

-



Preparation time

10 minutes



Execution time


10 minutes

This content can also be found online at:

<http://localhost:1337/c/64283c585e30a7000275ea4e>


PHYWE

Teacher information



Application

PHYWE



Experimental set-up:

Domed mirror

Light propagates in a straight line. If a light beam hits a reflecting object, the light beam also spreads out in a straight line from there.

With a concave mirror, the surface of the mirror is curved. The radius is similar to that of a concave mirror.

Therefore, a convex mirror reflects the rays exactly as if plane mirrors were aligned tangentially to the circular arc.

Other teacher information (1/2)

Prior knowledge



The students need theoretical knowledge about the straight-line, ray-shaped propagation of light and that objects reflect light rays.

Principle



The aim is to investigate how parallel and divergent light rays are reflected at the camber mirror.

Other teacher information (1/2)

PHYWE

Prior knowledge



The students need theoretical knowledge about the straight-line, ray-shaped propagation of light and that objects reflect light rays.

Principle



The aim is to investigate how parallel and divergent light rays are reflected at the camber mirror.

Other teacher information (2/2)

PHYWE

Learning objective



The students should gain knowledge about the principles of light reflection. This experiment is about the reflection behaviour of a convex mirror.

Tasks



Ask the students to observe how a convex mirror reflects the different rays that fall on it.

Additional teacher information

PHYWE

Note



The distance between luminaire and mirror should not be too large. The course of the reflected light rays can be observed particularly well when the light falls at an angle. To make the circle template, the pattern can be copied, glued onto thin cardboard and cut to size.

Safety instructions


PHYWE



- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Equipment

Position	Material	Item No.	Quantity
1	PHYWE Demo Physics board with stand	02150-00	1
2	Halogen lamp for experiments, 12V/50W, with magnetic base	08270-20	1
3	Opt. block, planoconcave, magn. held	08270-03	1
4	Concave/convex mirror, magnet held	08270-12	1
5	PHYWE Multitap transformer DC: 2/4/6/8/10/12 V, 5 A / AC: 2/4/6/8/10/12/14 V, 5 A	13533-93	1
6	G-clamp	02014-00	2

Equipment

PHYWE

Position	Material	Item No.	Quantity
1	PHYWE Demo Physics board with stand	02150-00	1
2	Halogen lamp for experiments, 12V/50W, with magnetic base	08270-20	1
3	Opt. block, planoconcave, magn. held	08270-03	1
4	Concave/convex mirror, magnet held	08270-12	1
5	PHYWE Multitap transformer DC: 2/4/6/8/10/12 V, 5 A / AC: 2/4/6/8/10/12/14 V, 5 A	13533-93	1
6	G-clamp	02014-00	2

Set-up and Procedure (1/2)

PHYWE




Fig.1:

Convex mirror with 5-slit aperture

- Draw optical axis
 - Draw an arc of a circle with radius $r = 200 \text{ mm}$ on the board using a template or compass.
 - Set mirror on circular arc
- Mark centre of curvature ($r = 200 \text{ mm}$)
- Position the luminaire with the 5-slit aperture so that the light beams run in the direction of the optical axis.

Set-up and Procedure (2/2)

PHYWE





Fig.2:

Convex mirror with 5-slit aperture


- Change the position of the luminaire so that the light beams fall on the convex mirror from different directions; observe the course of the beams.
- Move the adhesive luminaire to its original position and place the plane-concave model body in front of the luminaire ; cover the slit openings above the model body.
- Allow divergent beams from different positions and also from different distances of the luminaire to strike the convex mirror; observe the course of the beams.

PHYWE

Report




Task 1



Write down your observations of the different experimental set-ups.

Task 1

PHYWE

Write down your observations of the different experimental set-ups.