curricuLAB<sup>®</sup> PHYWE

# **Refraction of light at the air-water boundary**



Physics

ght & Optic

Difficulty level

easy

Group size

22

Preparation time



10 minutes

10 minutes

This content can also be found online at:



http://localhost:1337/c/64286c5c5e30a7000275eaea





# **Teacher information**

## **Application**

### **PHYWE**



Experimental set-up: Refraction of light in water Light rays not only reflect, they are refracted in the transition between different media.

With each new transition between two media ( air / glass ) or ( air / water ) the light beam is slightly deflected.

These distractions are technically exploited, e.g. to correct a visual impairment with the help of glasses.









#### **PHYWE**

#### Note

It is recommended to colour the water in the cuvette, e.g. with Paten Blue-V (48376-04).

Then the course of the light beam in the water is even more clearly visible.

If this experiment is also to demonstrate the reversibility of the light path during refraction, then it is advisable to additionally use the light box so that its light beam runs in the opposite direction to the light beam of the adhesive luminaire, and to operate both luminaires alternately.

## Additional teacher information (2/2)

#### **PHYWE**

Note

If you want to work out the different refraction at the air-glass and air-water transitions, then the angle of refraction is determined for the same angle of incidence (e.g. 45°) by marking the refracted beam.

Transition air-water:  $lpha=45\degree; eta=32\degree$ 

Transition air-glass:  $lpha=45\degree; eta=28\degree$ 

Therefore, it is realised that light is refracted more strongly at the air-glass transition than at the air-water transition.



## Safety instructions

#### **PHYWE**



• The general instructions for safe experimentation in science lessons apply to this experiment.

# **PHYWE**

# Student information



## **Motivation**

#### **PHYWE**



Straw in water glass

Some of you may have noticed that a straw in a glass of water seems to deform below the water surface.

Of course, the straw doesn't really deform, it just looks that way because the water refracts the light rays differently than the air.

The present experiment aims to explain how this "kink in the optic" occurs.



www.phywe.de

## Equipment

Position	Material	Item No.	Quantity
1	PHYWE Demo Physics board with stand	02150-00	1
2	Halogen lamp for experiments, 12V/50W, with magnetic base	08270-20	1
3	Cuvette, magnet held, 230x75 mm	08270-08	1
4	PHYWE Multitap transformer DC: 2/4/6/8/10/12 V, 5 A / AC: 2/4/6/8/10/12/14 V, 5 A	13533-93	1
5	G-clamp	02014-00	2



www.phywe.de

**PHYWE** 

#### Set-up and Procedure



1-slit orifice with water-filled cuvette

- Draw a horizontal line to mark the interface approximately in the middle of the lower half of the adhesive board; erect incidence plumb bob
- Place the cuvette filled with water in such a way that the water surface is level with the marking of the interface.
- Position the luminaire with the 1-slit aperture so that the beam runs in the direction of the incident plumb bob and does not experience any refraction ( $\alpha\alpha = \beta = 0^{\circ}$ ); readjust the cuvette if necessary.
- Vary the angle of incidence while observing the course of the light beam

**PHYWE** 



## Report



Task 1	PHYWE
	Write down your observations on the refractive behaviour of the light beam.





Task 3			PHYWE
Fill in the blanks!	When the light beam passe is refracted and perpendicular lie in then the light beam is away air the sam Check	es obliquely from the water into by the perpendicular. Incident ray, ref plane. If the angle of incid refracted. ne not	, it fracted ray dence is 0°,

Slide			Score / Total
Slide 14: Refraction			0/5
Slide 15: Angle of incidence			0/4
		Total	0/9
Sol	ons 🛛 🤁 Repea	at 🖹 Export text	