



# Refraction by a prism





This content can also be found online at:



http://localhost:1337/c/642877cd5e30a7000275eb14






# **PHYWE**



# **Teacher information**

# **Application** PHYWE



Refraction of light at prism

After the students have learned how a light beam behaves when passing through two parallel transitions, this experiment is intended to draw attention to the behaviour at different angles.

The refractive index always describes the refractive behaviour perpendicular to the transition surface.

Due to the different spatial directions of the prism walls, the light beam is also refracted differently.





## Other teacher information (1/2)

#### **PHYWE**

### Prior knowledge



**Principle** 



Students need prior theoretical knowledge about the straight-line, ray-shaped propagation of light. They should have learned about light refraction and refractive indices.

The aim is to show how a beam of light passes through a prism.

# Other teacher information (2/2)

**PHYWE** 

# Learning objective



**Tasks** 



refracted twice. if the final medium is also the initial medium, the beam runs parallel.

The students should observe that when light crosses two media transitions, it is also

The students should observe the experiment and realise that passing through mutually inclined planes also leads to a tilting of the light beam.





## **Additional teacher information**

#### **PHYWE**

Note



In this experiment the refractive angle is 45°. If you want to show that the total deflection is also dependent on the refractive angle, then the experiment can be repeated with the model body trapezoid (08270-05), which contains an angle of 60°.

The fact that after the refraction of the light beam at the prism, colour appearances also occur due to dispersion should not be emphasised at this point.

## **Safety instructions**





 $\circ\,$  The general instructions for safe experimentation in science lessons apply to this experiment.





# **PHYWE**



# **Student information**

## Motivation PHYWE



The picture opposite is a symbolic image.

This experiment deals with the refraction of a light beam passing through a prism.

Unlike the parallel glass pane, the light is refracted here depending on the angle of the prism side.






## **Equipment**

| Position | Material  | Item No. | Quantity |
|----------|---|----------|----------|
| 1        | PHYWE Demo Physics board with stand   | 02150-00 | 1        |
| 2        | Halogen lamp for experiments, 12V/50W, with magnetic base                         | 08270-20 | 1        |
| 3        | Opt. block,triangular,magnet held   | 08270-06 | 1        |
| 4        | PHYWE Multitap transformer DC: 2/4/6/8/10/12 V, 5 A / AC: 2/4/6/8/10/12/14 V, 5 A | 13533-93 | 1        |
| 5        | G-clamp   | 02014-00 | 2        |





## Set-up and Procedure PHYWE



1 slit diaphragm directed towards prism

- Place the prism (model body right triangle) on the adhesive board approximately as shown.
- Aim the adhesive luminaire with 1-slit diaphragm at the prism; mark the point of incidence of the light beam (P)
- Move the luminaire to another position, whereby the light beam should always hit P, or rotate the prism around the point P.
- Observing the course of the light beam at different angles of incidence


# **PHYWE**



# Report







# The prism used has an angle of .... 45° 180° 60°





