


Equivalence of inertial mass and gravitational mass with the demonstration track and the timer 4-4



Physics

Mechanics

Energy conservation & impulse



Difficulty level

medium



Group size

2



Preparation time

20 minutes



Execution time

10 minutes

This content can also be found online at:



<http://localhost:1337/c/600497c427aa1c00038a12df>


PHYWE



General information

Application

PHYWE



Experiment set-up

The mass of a body can be determined in different ways; by weighing the heavy mass or by measuring the acceleration at a known constant accelerating force, resulting in the inertial mass. It is assumed that heavy mass and inertial mass are equivalent.

In this experiment, a car with different heavy masses is accelerated on the roller track, its inertial mass is determined and the relationship between heavy and inertial mass is derived.

Other information (1/2)

PHYWE

Prior knowledge



Students should have a good understanding of the classical equations of motion and Newton's axioms.

Scientific principle



Up til now, the term mass has always been used in connection with the weight of a body and its movement.

In principle, these are two different properties, the heavy mass and the inertial mass.

In general, they are not distinguished because of their equivalence, but it should be noted that they are basically independent properties.

Other information (2/2)

PHYWE

Learning objective



It should be possible to distinguish between the two physical properties of weight and mass of a body.

Tasks



1. Determination of the total heavy mass and the driving mass by weighing.
2. Determination of the acceleration of the trolley and from this the inertial mass of the system.
3. Comparison of heavy and inertial mass.

Safety instructions

PHYWE

The general instructions for safe experimentation in science lessons apply to this experiment.

Equipment

Position	Material	Item No.	Quantity
1	Demonstration track, aluminium, 1.5 m	11305-00	1
2	Cart, low friction sapphire bearings	11306-00	1
3	Shutter plate for low friction cart, width: 100 mm	11308-00	1
4	Needle with plug	11202-06	1
5	Tube with plug	11202-05	1
6	Plasticine, 10 sticks	03935-03	1
7	End holder for demonstration track	11305-12	1
8	Starter system for demonstration track	11309-00	1
9	Magnet w.plug f.starter system	11202-14	1
10	Pulley for demonstration track	11305-10	1
11	Holder for pulley	11305-11	1
12	Weight for low friction cart, 400 g	11306-10	1
13	Slotted weight, black, 50 g	02206-01	3
14	Slotted weight, black, 10 g	02205-01	4
15	Weight holder, silver bronze, 1 g	02407-00	1
16	Slotted weight, blank, 1 g	03916-00	20
17	Silk thread, l = 200 m	02412-00	1
18	Light barrier, compact	11207-20	4
19	Holder for light barrier	11307-00	4
20	PHYWE Timer 4-4	13604-99	1
21	Connecting cord, 32 A, 1000 mm, red	07363-01	4
22	Connecting cord, 32 A, 1000 mm, yellow	07363-02	5
23	Connecting cord, 32 A, 1000 mm, blue	07363-04	5
24	Portable Balance, OHAUS CR2200	48914-00	1


PHYWE



Set-up and Procedure

Set-up (1/6)

PHYWE



Launching device without shock


1. In order to compensate for minor friction effects, the track must be set at a slight angle using the adjusting screws on the feet so that the measuring carriage just does not start to roll to the right.

2. A launching device shall be installed at the left end of the runway.

Note; to start the trolley with initial pulse, the starting device must be mounted so that the punch moves away from the measuring trolley when triggered.

Set-up (2/6)

PHYWE




End bracket with plasticine

3. A tube filled with plasticine is attached to the end bracket at the right end of the track to slow the car down without hard impact.

4. The deflection roller is attached to the right end of the track with the bracket for deflection roller and the incremental wheel is inserted.

Set-up (3/6)

PHYWE




Fastening the thread to the trolley

5. The measuring trolley is equipped with the holding magnet with plug as well as the cover for measuring trolley ($b = 100 \text{ mm}$).

6. The beginning of the thread is inserted from above into the vertical hole of the carriage end cap and fixed by inserting the needle with plug from the front.

Set-up (4/6)

PHYWE




Positioning of the weight plate

7. The thread is placed over the incremental wheel of the deflection pulley and knotted at the end to the weight plate so that it hangs freely directly below the wheel, as shown in the illustration. The weight plate together with the 5-20 slotted weights (1g each) lying on it serve as a constant accelerating force. The thread must run parallel to the track.

8. The mass of the car can be varied by means of the black painted weights.

Set-up (5/6)

PHYWE



Release of the light barriers after aperture passage

9. The four forked photoelectric sensors are mounted on the roadway with the photoelectric sensor holders and distributed evenly over the measurement section.


When rolling the trolley, make sure that all light barriers can be passed through by the rear part of the screen before the weight plate touches the ground.

10. The forked light barriers are connected from left to right in sequence to the sockets in fields "1" to "4" of the timing device.

The yellow sockets of the light barriers are connected to the yellow sockets of the measuring device, the red sockets to the red sockets and the blue sockets of the light barriers to the white sockets of the time measuring device.

Set-up (6/6)

PHYWE




Connecting the light barriers and the starting device

11. The starting device must be connected to the two "Start" connection sockets of the timing device.

Make sure that the polarity is correct.



The red socket of the starting device is connected to the yellow socket of the timing device.

12. The two slide switches on the timing device are set to the right-hand position "falling edge" () to select the trigger edge.

Procedure (1/2)

PHYWE

1. The measuring carriage is released by the starter and experiences constant acceleration until the weight plate touches the ground. It then continues to roll at a constant speed.

2. First the times $t_1 \dots t_4$ are measured in mode 2 (), which are the times required from the start until the respective light barriers are reached. Subsequently, a measurement is made in mode 1 () to determine the corresponding speeds.

When performing this measurement, the shading times are $\Delta t_1 \dots \Delta t_4$ of the four forked light barriers; from these in turn, the average speed during the corresponding passage is later calculated via the aperture length (100 mm).

Procedure (2/2)

PHYWE

3. The measuring times are recorded up to five times for different trolley masses with constant mass of the weight plate. Before each measurement, press the "Reset" key to reset the displays. Before starting a new series of measurements, the total heavy mass is measured with a balance. $m_s = m_W + m_A$ (car mass and accelerating mass of the driving body).

4. The heavy mass of a body is proportional to the weight force F_g acting on it. The drive body (m_A , weight plate with weights) can be determined with a balance.

Evaluation (1/8)

PHYWE

Measurement example for light barrier 1

F_g in N	m_s in kg	t_{1m} in s	Δt_{1m} in s	v_1 in m/s	a_1 in m/s ²
0,0981	0,41	1,235	0,293	0,341	0,276
0,0981	0,55	1,436	0,344	0,291	0,202
0,0981	0,81	1,781	0,43	0,233	0,131

Measurement example for light barrier 2

F_g in N	m_s in kg	t_{2m} in s	Δt_{2m} in s	v_2 in m/s	a_2 in m/s ²
0,0981	0,41	1,776	0,22	0,454	0,256
0,0981	0,55	2,067	0,258	0,387	0,187
0,0981	0,81	2,573	0,323	0,309	0,12

Evaluation (2/8)

PHYWE

Measurement example for light barrier 3

F_g in N	m_s in kg	t_{3m} in s	Δt_{3m} in s	v_3 in m/s	a_3 in m/s ²
0,0981	0,41	2,202	0,185	0,541	0,245
0,0981	0,55	2,566	0,217	0,46	0,179
0,0981	0,81	3,2	0,273	0,367	0,115

Measurement example for light barrier 4

F_g in N	m_s in kg	t_{3m} in s	Δt_{3m} in s	v_3 in m/s	a_3 in m/s ²
0,0981	0,41	2,555	0,163	0,613	0,24
0,0981	0,55	2,98	0,192	0,52	0,175
0,0981	0,81	3,722	0,243	0,412	0,111

Evaluation (3/8)

PHYWE

Mean values of the measurement examples

F_g in N	m_s in kg	a_m in m/s ²	$1/a$ in s ² /m
0,0981	0,41	0,254	3,94
0,0981	0,55	0,186	5,38
0,0981	0,81	0,119	8,4

Evaluation (4/8)

PHYWE


1. From the five measurements each of $t_1 \dots t_4$ and $\Delta t_1 \dots \Delta t_4$, the mean values $t_{1m} \dots t_{4m}$ and $\Delta t_{1m} \dots \Delta t_{4m}$ can be determined.
2. Calculate the velocities from the shading times $v_{im}(t_{im}) = b / \Delta t_{im}$ and the accelerations $a_i = v_{im}(t_{im}) / t_{im}$ with the aperture length $b = 0.1$ m.
3. The four accelerations belonging to a series of measurements a_i are averaged. From the resulting values a_m the reciprocal values of the acceleration $1/a_m$ is determined (see measurement example 1).

Evaluation (5/8)

PHYWE

4. Finally, one carries the reciprocal accelerations $1/a_m$ against the total heavy mass m_s from the trolley and drive mass in a coordinate system and determines the gradient of the straight line through these points.

For comparison, the values of the measurement example are shown in the figure.



Evaluation (6/8)

PHYWE

5. Every body has the property of gravity. Due to the earth's gravity, it experiences an accelerating force, the weight force F_g which is proportional to its heavy mass m_s (Newton's law of gravity).

In this experiment, the heavy mass consists of $m_s = m_W + m_A$ from the mass of the car and the mass of the driving body.

6. If no force acts on a body, it remains in its state of motion (Newton's 1st axiom). When an external force acts on the body, it opposes the change of its state of motion with a resistance, the property of inertia.

The greater its inertial mass m_t the more force must be applied to accelerate the body to a certain value (Newton's 2nd axiom).

Evaluation (7/8)

PHYWE

7. From the values of the $(1/a, m_s)$ coordinate system shown in Fig. 8, a straight line through the origin is obtained. At constant accelerating force, the heavy mass of the system is inversely proportional to the acceleration of the carriage:

$$m : s \propto \frac{1}{a} \text{ for } F_g = \text{const.}$$

With larger heavy mass M_s the acceleration decreases a .

8. According to Newton's second axiom $F_g = m_t \cdot a$, the mass m_t of a body is inversely proportional to its acceleration a .

If both heavy and inertial mass are inversely proportional to acceleration, then heavy and inertial mass must be proportional to each other:

$$m_s = k \cdot m_t.$$

Evaluation (8/8)

PHYWE

9. The determination of the proportionality factor k is done via a gradient triangle on the graph from slide 20.

The gradient gives a value of $10.1 \text{ s}^2/(\text{kg}\cdot\text{m})$, the reciprocal of which is $0.0990 \text{ kg}\cdot\text{m}/\text{s}^2$ of the weight force $F_g = 0.0981 \text{ kg}\cdot\text{m}/\text{s}^2$ (cf. measurement example 1). Heavy and inertial mass must therefore be identical:

$$m_s = m_t.$$

Notes

PHYWE

1. To reduce the distance between the weight adjuster and the incremental wheel, the thread length can be shortened by turning the needle with plug several times on the carriage, thus winding up the thread.

2. The from Δt_i calculated speeds v_i are, strictly speaking, not instantaneous velocities, since acceleration continues to act on the carriage as it passes through the light barrier.

The velocities thus arise from a secant slope, but not from a tangent slope, of the graph of $s(t)$. Using $\Delta s = 0.1 \text{ m}$, a systematic error of about 2 % must be expected.