


# Two-sided lever



P1253100

Physics

Mechanics

Forces, work, power &amp; energy



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time


20 minutes

This content can also be found online at:

<http://localhost:1337/c/661cdb5c9cc71c0002761f33>

PHYWE

## General information



## Application

PHYWE




Fig. 1: Two-sided lever using the example of pliers

Levers are force-transforming devices. They are often used to generate larger forces with small forces

With a **two-sided lever** both forces act on different sides as seen from the axis of rotation.

## Other information (1/2)

PHYWE

### Prior knowledge



No prior knowledge required.

### Principle



The aim is to demonstrate that a two-sided lever is in equilibrium if the product of the applied force and the lever arm is the same on both sides.

## Other information (2/2)

PHYWE

### Learning objective



The experiment is designed to help students understand when a two-sided lever is in equilibrium.

The relationship between force, length of the lever arms and torque is also shown.

### Tasks



The students are supposed to carry out observations and measurements to determine the conditions for equilibrium with a two-sided lever.

## Safety instructions


PHYWE

The general safety instructions for experimentation in science lessons apply.

## Equipment

Position	Equipment	Item no.	Quantity
1	<a href="#">PHYWE Demo Physics board with stand</a>	02150-00	1
2	<a href="#">Rod on fixing magnet</a>	02151-02	1
3	<a href="#">Torsion dynamometer</a>	03069-03	1
4	<a href="#">Scale for demo board</a>	02153-00	1
5	<a href="#">Weight holder, 10 g</a>	02204-01	1
6	<a href="#">Slotted weight, silver-bronze, 10 g</a>	02205-03	2
7	<a href="#">Slotted weight, silver-bronze, 10 g</a>	02205-03	2
8	<a href="#">Slotted weight, silver-bronze, 50 g</a>	02206-03	1
9	<a href="#">Slotted weight, silver-bronze, 50 g</a>	02206-03	1
10	<a href="#">Lever</a>	03960-00	1
11	<a href="#">Pointer for demo lever</a>	03963-00	1
12	<a href="#">Marker, black</a>	46402-01	1
13	<a href="#">Screw clamp</a>	02014-00	2

PHYWE



## Setup and procedure

### Setup

PHYWE

- Place the rod on the fixing magnet on the upper part of the demo board and place the lever in the centre on the rod.
- Draw a line vertically downwards from the axis with the marker.
- Attach the pointer for the demo levers (its tip lies exactly on the line drawn when the lever is balanced).

## Procedure (1/2)

PHYWE




Fig. 2: Beam balance with weight and dynamometer


- Put on the dynamometer and measure the weight force - in the following called  $F_1$  - for the weight holder with all slotted weights; note  $F_1$  in Table 1, upper part.
- Hook the pull cord of the dynamometer on the right, at mark no. 10, and the weight holder with the slotted weights at no. 10, on the left.
- Move the dynamometer until the lever is horizontal and the pull cord is perpendicular to it (Fig. 2)
- Read  $F_2$  on the dynamometer and note it in Table 1.

## Procedure (2/2)

- Shorten the lever arm  $I_1$  step by step; measure and record the force required to balance the lever in each case  $F_2$  (see specifications in Table 1, upper part).
- Remove two 50g slotted weights from the weight holder, measure and record the weight force  $F_1$ .
- Hook the weight holder at mark no. 9, left, and leave it there for the following procedure.
- Shorten lever arm  $I_2$  step by step; to do this, hook the pull cord of the dynamometer at marks no. 10, 9, ..., 6 (see Table 1, lower part), measure and note down  $F_2$  each time.

PHYWE

# Evaluation



## Evaluation (1/4)

PHYWE

### Observation:

Marke Nr.	l1 / cm	F1 / N	(F1*l1) / (N*cm)	Marke-Nr.	l2 / cm	F2 / N	(F2*l2) / (N*cm)
10	20	1.95	39.0	10	20	1.96	39.2
8	16	1.95	31.2	10	20	1.55	31.0.
6	12	1.95	23.4	10	20	1.19	23.8
4	8	1.95	15.6	10	20	0.79	15.8
2	4	1.98	7.8	10	20	0.39	7.8
9	18	0.98	17.6	10	20	0.88	17.6
9	18	0.98	17.6	9	18	0.99	17.8
9	18	0.98	17.6	8	16	1.10	17.6
9	18	0.98	17.6	7	14	1.27	17.8
9	18	0.98	17.6	6	12	1.47	17.6

Table 1: Example measurements

## Evaluation (2/4)

PHYWE

The distance between the marks is 2 cm. The resulting  $I_1$  and  $I_2$  are to be entered in Table 1.

After calculating the products  $F \cdot I$  for the left and right lever side  $F_1 \cdot I_1 = F_2 \cdot I_2$  can be recognised.

A two-sided lever is in equilibrium if the products of the forces acting to the left and right of the centre of rotation and their lever arms are equal. Here the lever arms are the distances between the points of application of the forces and the centre of rotation of the lever.

## Evaluation (3/4)

PHYWE

The forces  $F_1$  were selected in order to utilise the measuring range of the dynamometer as much as possible. The pointer for demo levers promotes the association with an application of the lever, the beam balance.

If the suggestion of this association is not desired at this point, leave the pointer out and draw a horizontal line behind the lower edge of the lever.

If the term torque can be used, the law found is generalised into the torque theorem:

At the two-sided lever, the sum of all torques is zero in the case of equilibrium.

The following applies with two forces:

$$\vec{M} = \vec{F}_1 \vec{I}_1 + \vec{F}_2 \vec{I}_2 = 0$$

Clockwise and anti-clockwise torques have different signs:

$$\vec{F}_1 \vec{I}_1 = -\vec{F}_2 \vec{I}_2$$

## Evaluation (4/4)

PHYWE

In the event that  $\vec{F}$  and  $\vec{I}$  form an angle of  $90^\circ$ , is

$$|\vec{M}| = M = |\vec{F} * \vec{I}| = F \cdot I$$

In the general case, however

$M = F \cdot I \cdot \sin \alpha$  is also known as the effective length of the force arm or effective lever length (see Fig. 2).

The unit of  $M$  is the Newton metre (Nm), i.e. the same as the unit for mechanical work, which is often a source of confusion for students.




Fig. 3: Calculation of the lever arm

## Task 1

PHYWE

What is the effective length of the lever arm or effective lever length?

☐  $I \cdot \sin \alpha$

☐  $I \cdot \cos \alpha$

☐  $F \cdot I \cdot \sin \alpha$

☐  $F \cdot \sin \alpha$

☒ Check

Slide

Score / Total

Slide 16: Effective lever length

0/1

Total score



Show solutions



Repeat