


# Pressure distribution in flowing liquids



P1297500

Physics

Mechanics

Mechanics of liquids &amp; gases



Difficulty level

medium



Group size

-



Preparation time

-



Execution time

-

This content can also be found online at:

<http://localhost:1337/c/66e3f607b9fd8b000221b2b5>

PHYWE



## General information

### Application (1/2)

PHYWE




Fig. 1: Experimental setup

Static pressure is the pressure in a flowing liquid or in a flowing gas that is measured perpendicular to the direction of flow. This static pressure depends on the flow velocity. The relationship is described by Bernoulli's law: The greater the flow velocity of a liquid or gas, the lower the static pressure.

## Application (2/2)

PHYWE

This principle has already been used in everyday life, e.g. the atomiser. In such an atomiser, air flows through a nozzle. The flow velocity is relatively high in the area of the nozzle and the static pressure is therefore low according to Bernoulli's law. Due to this low pressure, the liquid in a vessel is literally sucked in. It emerges from an opening and is atomised by the fast-flowing air.

## Other information (1/2)

PHYWE

### Prior knowledge



Students should have prior knowledge of static pressure and "Bernoulli's law".

### Principle



The aim is to investigate the behaviour of pressure in a flowing liquid.

## Other information (2/2)

PHYWE

### Learning objective



The students are to be shown and taught about the pressure distribution in flowing liquids in this experiment. The students are supposed to understand "Bernoulli's law".

### Tasks



Students are to observe how pressure behaves in a flowing liquid. The students should be familiar with static pressure and understand Bernoulli's law.

## Safety instructions

PHYWE



The general safety instructions for experimentation in science lessons apply.

## Theory

PHYWE

Pressure distribution in flowing liquids refers to the static pressure that can be measured perpendicular to the direction of flow in a flowing liquid or in a flowing gas. This static pressure depends on the flow velocity. This in turn increases as the pipe cross-section is reduced. The relationship between flow velocity and pressure is described by Bernoulli's law. It states that the greater the flow velocity of a liquid or gas, the lower the static pressure.


With a constantly flowing liquid, the reduced diameter can lead to an increased flow velocity and also to a decreased static pressure.

## Equipment

Position	Equipment	Item no.	Quantity
1	<a href="#">PHYWE Demo Physics board with stand</a>	02150-00	1
2	<a href="#">Hook on fixing magnet</a>	02151-03	1
3	<a href="#">Clamping holder, 0-13 mm, fixing magnet</a>	02151-07	2
4	<a href="#">Scale for demonstration board</a>	02153-00	1
5	<a href="#">Pointers f. Demonst.Board, 4 pcs</a>	02154-01	1
6	<a href="#">Flow pipe w.uniform cross-section</a>	02765-00	1
7	<a href="#">Flow pipe w.varying cross-section</a>	02766-00	1
8	<a href="#">Storage tray, 413 x 120 x 100 mm</a>	47325-01	1
9	<a href="#">Rubber hose, various diameters (8 mm)</a>	39283-00	1
10	<a href="#">Beaker, 100 ml, plastic (PP)</a>	36011-01	1
11	<a href="#">Microspoon, steel</a>	33393-00	1
12	<a href="#">Pinchcock, width 20 mm</a>	43631-20	1
13	<a href="#">Pipette with rubber bulb, long</a>	64821-00	1
14	<a href="#">Patent blue-V, 25 ml</a>	48376-05	1
15	<a href="#">Screw clamp</a>	02014-01	2

PHYWE

## Setup and procedure



### Setup and procedure(1/3)

PHYWE




Fig. 2

#### 1st experiment

- Attach the flow pipe with a constant diameter horizontally to the demo board using the clamps and secure it against tilting using the hook (Fig. 2)
- Place the bowl under the flow pipe (to catch overflowing water while the flow is being regulated to a desired height of the water column and the air bubbles are being removed).
- Push a tube of suitable length onto both ends of the flow pipe.

## Setup and procedure (2/3)

PHYWE

- Connect one hose to the tap, fit the other with the hose clamp and lead it to the drainage basin.
- With the hose clamp loosened, carefully turn on the tap.
- Now use the hose clamp to regulate the flow so that the first riser does not overflow and the last one has a water column of several centimetres.
- Colour some water in the beaker and add a few drops of the solution to each riser using the pipette.
- Mark the height of the water columns in the risers with arrows.
- Formulate a statement about the height of the water columns.

## Setup and procedure (3/3)

PHYWE




Fig. 2

### 2nd experiment

- Now replace the pipe with the constant diameter from the first setup with a pipe with a varying diameter.
- Proceed as in experiment 1.



## Result

PHYWE

### 1st result

The height of the water columns decreases from the beginning of the flow pipe.

### 2nd result

The height of the water columns - seen from the beginning of the flow tube - is much lower in the last two risers than in the first two. It is very low in the centre riser, which is attached to the narrowing of the flow tube.

## Evaluation (1/2)

PHYWE

### 1st evaluation

After removing the flow pipe from the demo board, the scale is placed at the tips of the pointers. This shows that the height of the water columns in the riser tubes decreases linearly. The height of the water columns is a measure of the pressure ( $p = \rho * g * h$ ), which prevails perpendicular to the direction of flow and is called static pressure.

So the following applies: In a flowing liquid, the static pressure decreases linearly in the direction of flow if the flow cross-section is the same everywhere.

The cause of the linear decrease in static pressure are frictional forces that occur between the liquid and the pipe wall as well as within the liquid.

## Evaluation (2/2)

PHYWE

### 2nd evaluation

Within the narrowing of the flow pipe, the static pressure is much lower than before and after it. In order for the same amount of liquid to flow through all pipe cross-sections per unit of time, the liquid particles must be accelerated. This happens at the expense of the static pressure.

## Notes

PHYWE

If the friction could be eliminated, then all water columns would have the same height in experiment 1 and all but the column above the narrowing would have the same height in experiment 2.

The Bernoulli equation applies to the steady flow of incompressible, frictionless liquids:

Total pressure = static pressure + dynamic pressure or

$$p_{ges} = p + \frac{1}{2}\rho v^2 \text{ with}$$

$v$  = Flow velocity and

$\rho$  = Density of the flowing liquid .

The dynamic pressure acts in the direction of the flow. It is called dynamic pressure in analogy to static pressure..