






Conductivity of aqueous solutions of electrolytes with Cobra SMARTsenseors of electrolytes



In this experiment, the students should find out why an undissolved (or not melted) electrolyte and distilled water are not or almost not conductive.

Physics		Electricity & Magnetism		Electric current & its effects			
	Difficulty level		Group size		Preparation time		Execution time
	medium		2		10 minutes		10 minutes

This content can also be found online at:



<https://www.curriculab.de/c/6808abc33e1f6d00029c34cf>


PHYWE

Teacher information



Application

PHYWE



Experimental setup

Salts, acids, and bases are electrolytes. In their pure form, they conduct electricity poorly or not at all, as they contain no — or only very few — freely mobile ions.

When dissolved in water, electrolytes dissociate into positive and negative ions.

Other teacher information (1/2)

PHYWE

Prior knowledge



No prior knowledge is required for this experiment.

Principle



If a voltage is applied to two electrodes immersed in the aqueous solution of an electrolyte, the ions move in the direction of the electrode with the opposite electrical polarity. Aqueous solutions of electrolytes are therefore electrically conductive.

Other teacher information (2/2)

PHYWE

Learning objective



In this experiment, the students should find out why an undissolved (or not melted) electrolyte and distilled water are not or almost not conductive.

Tasks



Investigate whether water in which substances are dissolved conducts electricity.

Safety instructions


PHYWE



- Put on safety goggles!
- Wear gloves!
- Please refer to the relevant safety data sheets for the H and P phrases.
- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Motivation

PHYWE

As soon as a thunderstorm approaches, you should leave the pool. It is said that it is also not safe to swim in another body of water when there is thunder and lightning in the air.

But why is it so dangerous to go swimming during a thunderstorm? This experiment investigates the conductivity of substances dissolved in water and thus allows conclusions to be drawn about the answer to this question.



Thunderstorm over the sea

Tasks

PHYWE



1. Set up the experiment and measure the conductivity of distilled water.
2. Repeat the experiment and determine the conductivity of several more substances including:
 - salt
 - salt diluted in water
 - drinking water
 - a diluted acid
 - a diluted base

Material

Position	Material	Item No.	Quantity
1	Cobra SMARTsense Current - Sensor for measuring electrical current \pm 1 A (Bluetooth + USB)	12902-01	1
2	Cobra SMARTsense Voltage - Sensor for measuring electrical voltage \pm 30 V (Bluetooth + USB)	12901-01	1
3	Angled connector module, SB	05601-02	2
4	Straight connector module, SB	05601-01	1
5	Interrupted connector module with sockets, SB	05601-04	2
6	Junction module, SB	05601-10	2
7	Angled connector module with socket, SB	05601-12	2
8	On-off switch module, SB	05602-01	1
9	Trough, grooved, w/o lid	34568-01	1
10	Copper electrode, 76 mm x 40 mm	45212-00	2
11	Alligator clips, bare, 10 pcs	07274-03	1
12	Connecting cord, 32 A, 250 mm, red	07360-01	2
13	Connecting cord, 32 A, 250 mm, blue	07360-04	2
14	Connecting cord, 32 A, 500 mm, red	07361-01	2
15	Connecting cord, 32 A, 500 mm, blue	07361-04	2
16	PHYWE Power supply, 230 V,DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
17	Sulphuric acid, 10%, tech.gr., 1000 ml	31828-70	1
18	Sodium hydroxide sol.,10%, 1000ml	31630-70	1
19	Water, distilled 5 l	31246-81	1
20	Emery paper, medium	01605-00	1
21	Spoon, with spatula end, 180 mm, plastic	38833-00	1
22	measureAPP - the free measurement software for all devices and operating systems	14581-61	1

Setup (1/3)

PHYWE

To measure with the **Cobra SMARTsense sensors**, the **PHYWE measureAPP** is required. The app can be downloaded free of charge from the respective app store (QR codes below). Please check that **Bluetooth is enabled** on your device (smartphone, tablet, desktop PC) before starting the app.



iOS



Android



Windows

Setup (2/3)

PHYWE

- Set up the experiment as shown in Figs. 1 and 2. The switch is open. Clean the grooved trough and the copper electrodes carefully before inserting the electrodes into the outer grooves of the trough.
- Fill the grooved trough about half full with distilled water.




Fig. 1





Fig. 2

Setup (3/3)

PHYWE




Screenshot of the app without power supply switched on

- Switch on both Cobra SMARTsense sensors by pressing and holding the power button on each device for about three seconds.
- Then open the measureAPP and connect it to both sensors. Adjust the display so that the measured values are shown as numbers. To do this, tap the "0.0" icon at the top of the app. A visual example is shown on the left-hand side of the screen.

Procedure (1/4)

PHYWE

- Set the power supply unit to 0 V and switch it on.
- Close the switch, increase the voltage at the power supply unit until the measureAPP displays 2 V. Measure the current and note the measured value in Table 1 in the log.
- Open the switch, empty and dry the groove tray.




Experimental setup

Procedure (2/4)

PHYWE

- Place the electrodes back in the grooved trough and fill the trough with a layer of saline salt about 2 cm high.
- Close the switch and measure the current at $U = 2 \text{ V}$. Note the measured value.
- Pour distilled water slowly onto the salt in the grooved trough and monitor the displayed current values in the measureAPP while doing so.



Procedure (3/4)

PHYWE

- Stir the saline solution with the spoon and measure the strength of the current at the end.
- Open the switch and note the measured value for I .
- Empty the groove tray, wash it thoroughly - just like the electrodes - and dry both. Put the electrodes back into the groove tray.
- Fill the grooved trough about half full with drinking water.
- Close the switch and measure the current again at $U = 2 \text{ V}$. Note the measured value.
- Empty the grooved trough with the switch open.
- Close the switch, carefully pour diluted acid into the groove trough, measure the current and make a note of the measured value.


Procedure (4/4)

PHYWE

- Open the switch, dispose of the aqueous solution properly, rinse and dry the groove tray and the electrodes with water.
- Proceed in the same way with the diluted base.
- Set the power supply unit to 0 V and switch it off.
- Dispose of the aqueous solution properly, rinse and dry the groove tray and the electrodes with water and finally wash your hands.

PHYWE

Report



Observation

PHYWE

Test part no.	Substances	Current I [mA]
1	Distilled water	<input type="text"/>
2	Salt	<input type="text"/>
3	Aqueous solution of a salt	<input type="text"/>
4	Drinking water	<input type="text"/>
5	Aqueous solution of an acid	<input type="text"/>
6	Aqueous solution of a base	<input type="text"/>

Task (1/4)

PHYWE

Summarise the results of the individual parts of the experiment in words.

Task (2/4)

PHYWE

Drag the words into the correct boxes!

In liquids, a only takes place if freely movable (migratory) are present due to dissociation. When a voltage is applied and thus an is present, the ions move in a directed manner. Electrical energy is converted into . A special feature of conduction processes in liquids that is important for applications is that the ions not only transport charges, but also transport substances.

thermal energy

electric field

ions

conduction process

☒ Check

Task (3/4)

PHYWE

Why, for example, does table salt not conduct electricity in its solid form? Why does distilled water conduct electricity only poorly — if at all — while ordinary tap water does conduct electricity, even if not very well?

Moving charged particles such as ions are required to conduct electricity.

does not have these, which is why it does not conduct electricity.

Although consists of electrically charged ions, these are not freely mobile. In contrast, contains small quantities of dissolved salts, which now contribute freely movable ions and thus lead to electrical conductivity.

Distilled water

drinking water

cooking salt

☒ Check

Task (4/4)

PHYWE

In electrical engineering, the ground is often used as a conductor. How can this be explained?

Slide


Score/Total

Slide 20: Functionality

0/4

Slide 21: Reason for conductivity

0/3

Total score  0/7 Show solutions Repeat Export text