






Galvanisation with Cobra SMARTsense



In a model experiment, the students are asked to apply a metallic copper coating to part of a sheet of iron.

Physics		Electricity & Magnetism		Electric current & its effects	
					
Difficulty level		Group size		Preparation time	
medium		2		10 minutes	
					
					Execution time
					10 minutes

This content can also be found online at:



<https://www.curriculab.de/c/6811e5a83e87bd0002fc904f>


PHYWE

Teacher information



Application

PHYWE



Experimental setup

For aesthetic purposes and to protect against corrosion, everyday objects are often coated with a thin layer of nickel, chrome, silver, or gold. This is achieved through an electrochemical process known as electroplating.

Other teacher information (1/2)

PHYWE

Prior knowledge



For this experiment, the students should be familiar with the fact that aqueous solutions conduct electricity.

Principle



If a voltage is applied to two electrodes immersed in the aqueous solution of an electrolyte, the ions move in the direction of the electrode with the opposite electrical polarity. As the charge transport is realised via mass transport, this can be used to coat substances with other substances.

Other teacher information (2/2)

PHYWE

Learning objective



In a model experiment, the students are asked to apply a metallic copper coating to part of a sheet of iron.

Tasks



Demonstrate in a model experiment how a sheet of iron can be coated with a layer of copper.

Safety instructions


PHYWE



- Wear safety goggles and protective gloves.
- Please refer to the relevant safety data sheets for the H and P phrases.
- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student information



Motivation

PHYWE

For aesthetic reasons and to protect against corrosion, everyday objects are often coated with a layer of nickel, chrome, silver, or gold. This is done electrochemically, and the process is called electroplating.



Galvanised objects

Tasks

PHYWE



1. Set up the experiment using the photos
2. Clean the electrodes
3. Start the galvanizing process using the instructions
4. Record your observations

Materials

Position	Material	Item No.	Quantity
1	Cobra SMARTsense Current - Sensor for measuring electrical current ± 1 A (Bluetooth + USB)	12902-02	1
2	Cobra SMARTsense Voltage - Sensor for measuring electrical voltage ± 30 V (Bluetooth + USB)	12901-02	1
3	Straight connector module, SB	05601-01	1
4	Angled connector module, SB	05601-02	4
5	Interrupted connector module with sockets, SB	05601-04	2
6	Junction module, SB	05601-10	2
7	On-off switch module, SB	05602-01	1
8	Trough, grooved, w/o lid	34568-01	1
9	Copper electrode, 76 mm x 40 mm	45212-00	2
10	Iron electrode, 76 x 40 mm	45216-00	2
11	Alligator clips, bare, 10 pcs	07274-03	1
12	Connecting cord, 32 A, 250 mm, red	07360-01	2
13	Connecting cord, 32 A, 250 mm, blue	07360-04	2
14	Connecting cord, 32 A, 500 mm, red	07361-01	1
15	Connecting cord, 32 A, 500 mm, blue	07361-04	1
16	PHYWE Power supply, 230 V,DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
17	Sulphuric acid, 10%, tech.gr., 1000 ml	31828-70	1
18	Water, distilled 5 l	31246-00	1
19	Emery paper, medium	01605-00	1
20	Copper-II sulphate,cryst. 250 g	30126-25	1
21	Denaturated alcohol (spirit for burning), 1000 ml	31150-70	1
22	Spoon, with spatula end, 180 mm, plastic	38833-00	1
23	measureAPP - the free measurement software for all devices and operating systems	14581-61	1

Setup (1/4)

PHYWE

To measure with the **Cobra SMARTsense sensors**, the **PHYWE measureAPP** is required. The app can be downloaded free of charge from the respective app store (QR codes below). Please check that **Bluetooth is enabled** on your device (smartphone, tablet, desktop PC) before starting the app.



measure App for Android operating system




measureAPP for iOS operating systems



measureAPP for tablets / PCs with Windows 10

Setup (2/4)

PHYWE




Experimental setup

- Set up the experiment as shown, initially leaving the switch open. Insert the electrodes into the grooved trough and connect them using the crocodile clips and short connecting leads so that the iron electrode is connected to the negative terminal, i.e. the cathode. Add some diluted sulphuric acid to the solution and stir gently.
- Clean the grooved trough, then carefully sand the electrodes with sandpaper. Wipe them with alcohol (e.g. spirit), and avoid touching them with your fingers afterwards to ensure the surfaces remain grease-free.

Setup (3/4)

PHYWE


- Fill the grooved trough two-thirds full with distilled water and then add about two spoons of copper sulphate to the water while stirring until the solution is saturated.
- Set the power supply unit to 0 V and switch it on.



Experimental setup

Setup (4/4)

PHYWE




Example screenshot of the measureAPP

- Start both Cobra SMARTsense sensors by pressing and holding the on/off button on each device for about three seconds.
- Then launch the measureAPP and connect to both sensors. Configure the display so that the measured values appear as numbers. To do this, tap on "0.0" at the top of the app. A visual example is shown on the left-hand side.

Procedure (1/2)

PHYWE

- Close the switch and increase the voltage at the power supply unit until a current of approx. 150 mA is reached.




Procedure (2/2)

PHYWE

- Observe the processes at the electrodes and record your observations in the lab log.
- After approximately 3 minutes, open the switch, set the power supply unit to 0 V, and switch it off.
- Rinse the iron electrode with water, examine the part that was immersed, and note your observations in the lab log.
- Dry the copper electrode, dispose of the aqueous solution properly, clean the grooved trough, and wash your hands thoroughly with soap.

PHYWE



Report

Observation (1/2)

PHYWE

Note down your observations on the processes during the current flow:

a) at the anode:

b) at the cathode:

Observation (2/2)

PHYWE

Note down your observations on the condition of the iron electrode after the experiment:

Task (1/2)

PHYWE

Drag the words into the correct boxes!

If contacts are in a conductive liquid, electricity flows. Liquids, e.g. acids, alkalis or water with dissolved salts, conduct electricity. These substances are also called . If a is applied, the positive cations move to the and the negative anions to the . takes place at the anode. Cations that get close enough to the cathode pick up from the cathode. They are therefore and accumulate as metal on the cathode.

anode

Oxidation

reduced

electrolytes

electrons

cathode

voltage

☒ Check

Task (2/2)

PHYWE

The process in which an electric current through a salt solution is used to coat the surfaces of conductive materials with a layer of metal is called electroplating. Give examples of electroplated objects.

Slide

Score/Total

Slide 19: Principle

0/7

Total score   0/7 Show solutions Repeat Export text