


Characteristic curve of a Zener diode



The students should use the experiment to see how a Z-diode differs from an ordinary rectifier diode.

Physics	Electricity & Magnetis	sm Electronics	3
Difficulty level	R Group size	Preparation time	Execution time
medium	2	10 minutes	10 minutes

This content can also be found online at:



http://localhost:1337/c/63170a3d9ebaee00039a3aee






PHYWE



Teacher information

Application PHYWE



Experimental setup

Z-diodes are silicon diodes with high doping of the p- and n-region. They behave like ordinary Si diodes in the forward direction. When a voltage is applied in the reverse direction, a high electric field strength is generated in the very narrow boundary layer. When a voltage dependent on the selected doping, the breakdown voltage, is exceeded, charge carrier pairs are released under the influence of the electric field, resulting in a strong increase in current intensity. This reduces the diode resistance. If the external voltage is increased further, the diode current rises sharply and produces an increasing voltage drop at the series resistor, while the voltage at the diode rises only very slightly.





Other teacher information (1/2)

PHYWE

Prior knowledge



Principle



Students should be familiar with how an ordinary diode works.

A Z-diode (also Zener diode) is a diode that is designed to be operated permanently in the reverse direction in the breakdown voltage range. The level of this breakdown voltage UBR is the main characteristic of a Z-diode and is specified in the data sheet. This is achieved by a heavily doped p+ and a heavily doped n+ layer. The strong recombination of both layers leads to a very small junction thickness and thus to high field strengths in the junction region.

Other teacher information (2/2)

PHYWE

Learning objective



The students should use the experiment to see how a Z-diode differs from an ordinary rectifier diode.

Tasks



Investigate the relationship between current and voltage on a Z-diode in the forward and reverse directions.





Safety instructions

PHYWE



• The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE



Student information





Motivation PHYWE

A diode is an electronic component that allows current to pass in one direction and blocks the flow of current in the other direction.

A Z-diode (also Zener diode) is a diode that is designed to be operated permanently in the reverse direction in the breakdown voltage range. The level of this breakdown voltage UBR is the main characteristic of a Z-diode and is specified in the data sheet. This is achieved by a heavily doped p+ and a heavily doped n+ layer. The strong recombination of both layers leads to a very small junction thickness and thus to high field strengths in the junction region.



Electronic components





Equipment

Position	Material	Item No.	Quantity
1	Straight connector module, SB	05601-01	3
2	Angled connector module, SB	05601-02	3
3	T-shaped connector module, SB	05601-03	2
4	Interrupted connector module with sockets, SB	05601-04	2
5	Straight connector module with socket, SB	05601-11	1
6	Angled connector module with socket, SB	05601-12	1
7	Resistor module 50 Ohm, SB	05612-50	1
8	Resistor module 100 Ohm, SB	05613-10	2
9	Z-diode module ZF4.7, SB	05652-00	1
10	Connecting cord, 32 A, 250 mm, red	07360-01	1
11	Connecting cord, 32 A, 250 mm, blue	07360-04	1
12	Connecting cord, 32 A, 500 mm, red	07361-01	2
13	Connecting cord, 32 A, 500 mm, blue	07361-04	2
14	PHYWE Power supply, 230 V, DC: 012 V, 2 A / AC: 6 V, 12 V, 5 A	13506-93	1
15	PHYWE Analog multimeter, 600V AC/DC, 10A AC/DC, 2 MΩ, overload protection	07021-11	2
16	Flashlight, w/o battery,medium	08164-00	1
17	Battery Type C 1.5 V - Pack of 2 pieces	07400-00	2





Set-up PHYWE


1st experiment

 Set up the circuit as shown in Figs. 1 and 2. Pole the Z-diode in forward direction. Select the measuring range 1 V- and 30 mA-. Pay attention to the correct polarity and the correct connection of the measuring instruments.

2nd experiment

• Change the circuit as in fig. 3 and 4.

Hint: Pay attention to the changed polarity of the Z-diode! Switch the measuring range for the voltage measurement to 10 V-.



Procedure PHYWE

1st experiment

- \circ Switch on the power supply unit. Increase the voltage U_N on the power supply unit in steps of 1 V from 0 V to 10 V. Measure the forward voltage in each case U_D and the forward current I_D . Enter the measured values in Table 1 in the report.
- Set the power supply unit to 0 V and then switch it off.

2nd experiment

- \circ Switch on the power supply unit. Increase the voltage again U_N in steps of 1 V each from 0 V to 10 V. Measure the blocking voltage in each case U_{Sp} and the reverse current I_{Sp} . Due to the changed polarity of the Z-diode, enter the measured values as negative in Table 2.
- Switch off the power supply unit.





PHYWE




Report


Result PHYWE

U_N [V]	U_D [V]	I_D [mA]	U_N [V]	U_D [V]	$I_D \left[mA ight]$	U_N [V]	U_D [V]	I_D [mA]	U_N [V]	U_D [V]	I_D [mA]
-10			-4			2			8		
-9			-3			3			9		
-8			-2			4			10		
-7			-1			5					
-6			0			6					
-5			1			7					
















Task (4/4) Read from the characteristic curve how large the change in voltage in the reverse direction is when the current increases from 20 mA to 30 mA. How could this property of the Z-diode be used?





			Score/Total
Slide 13: How it works			0/5
Slide 15: Series resisto	r		0/5
		Total score	0/10

