


Ohm's law



Students should recognise the proportional relationship between voltage and amperage and develop an understanding of the concept of resistance.



This content can also be found online at:



http://localhost:1337/c/6474e16c21530f000293d973



Tel.: 0551 604 - 0

Fax: 0551 604 - 107




PHYWE



General information

Application PHYWE



Experimental setup

This experiment is about the experimental derivation of Ohm's law. The law was named after Georg Simon Ohm. In 1825, he was looking for a mathematical relationship to calculate the "effect of flowing electricity" (current strength) as a function of the dimension and material of the electrical conductor.. After many experiments, he was able to establish the connection between electrical voltage U, electrical current I and electrical resistance R, by means of $U = R \cdot I$.

This experiment is about the experimental derivation of Ohm's law. The law was named after Georg Simon Ohm. In 1825 he was looking for a mathematical connection to calculate the "effect of flowing electricity" (current



Other information (1/2)

PHYWE

Prior knowledge



The students should be familiar with the functioning of the simple electric circuit and know the terms voltage and amperage.

Principle



By means of simple experiments, it can be demonstrated that the electric voltage is proportional to the electric current.

Other information (2/2)

PHYWE

Learning objective



Students should recognise the proportional relationship between voltage and amperage and develop an understanding of the concept of resistance.

Tasks



First, the current intensity is measured for two concrete resistors, at different voltages. In the second part, the resistor is replaced by an incandescent lamp.





Safety instructions

PHYWE

The general instructions for safe experimentation in science lessons apply to this experiment.

Theory

Ohm's law states that the electrical voltage U which decreases across an object is proportional to the electrical current intensity I. The proportionality factor is the material-specific resistance R.

The following therefore applies: $U = R \cdot I$.

The assumption of a constant resistance is only correct in certain cases, especially for metals of constant temperature. In this experiment, however, it is sufficient to consider the resistance as ideal and thus cosntant. In order to avoid irritations when using the word resistance before defining the physical term resistance, the components R_1 and R_2 should be first discussed. The evaluation should then clarify the function of the resistance (inhibiting or impeding the current).

The set of magnetic electrical symbols for the demo board is used for demonstrative labelling of the circuit.





Equipment

Position	Material	Item No.	Quantity
1	PHYWE Demo Physics board with stand	02150-00	1
2	Connector, angled, module DB	09401-02	4
3	Connector, T-shaped, module DB	09401-03	2
4	Connector interrupted, module DB	09401-04	3
5	Switch on/off, module DB	09402-01	1
6	Socket for incandescent lamp E10 ,module DB	09404-00	1
7	Resistor 50 Ohm, module DB	09412-50	1
8	Resistor 100 Ohm,module DB	09413-10	1
9	Electr.symbols f.demo-board,12pcs	02154-03	1
10	Connecting cord, 32 A, 1000 mm, red	07363-01	3
11	Connecting cord, 32 A, 1000 mm, blue	07363-04	3
12	PHYWE Power supply, universal, analog display DC: 18 V, 5 A / AC: 15 V, 5 A	13503-93	1
13	PHYWE Demo Multimeter ADM 3: current, voltage, resistance, temperature	13840-00	2
14	Filament lamps 12V/0.1A, E10, 10 pieces	07505-03	1
15	G-clamp	02014-01	2






PHYWE



Set-up and Procedure

Set-up PHYWE



 \circ Set up the experiment as shown in the illustration on the left and first adjust the resistance to $R_1=50\Omega$. Set the measuring ranges to 10V and 300mA.





Procedure PHYWE

- \circ Close the switch, switch on the power supply unit and starting from 0V increase the voltage in 2V steps; measure the respective amperage and note down the measured values.
- \circ Voltage to 0V and instead of R_1 , put the resistance $R_2=100\Omega$ into the circuit.
- Repeat as before and note it down.
- \circ Switch the power supply unit to 0V and install the incandescent lamp in the circuit instead of the technical resistor.
- \circ As with the 1st attempt, starting from 0V increase the voltage in 2V steps, measure the current reached in each case and note it down; observe the bulb during the series of measurements.

Evaluation (1/5)

PHYWE

Based on the test series for the first part of the test, the following table of measured values can be drawn up as an example:

U	I		U/I		
$\overline{\mathbf{V}}$	$\overline{\mathbf{A}}$		$\overline{V/A}$		
	(R ₁)	(R ₂)	(R ₁)	(R ₂)	
0	0	0	-	-	
2	0,039	0,020	51	100	
4	0,079	0,039	51	102	
6	0,120	0,060	50	100	
8	0,161	0,081	50	98	
10	0,201	0,100	50	100	





Evaluation (2/5)

PHYWE

If you now replace the resistor with the glow lamp, you will get the following measured values:

U	I	U/I	
$ \overline{\mathbf{v}} $	$\overline{\mathbf{A}}$	$\overline{\text{V/A}}$	
0	0	-	
2	0,037	54	
4	0,053	75	
6	0,069	87	
8	0,082	98	
10	0,093	108	

The brightness of the light bulb increases with increasing voltage or current.

\circ	Tı	ru	ie





Evaluation (3/5)


PHYWE

What can be seen from the graph of the measured values from the 1st part of the experiment?

There is a linear relationship between voltage and amperage $U \sim I$.

There is an exponential relationship between voltage and current.

This relationship is called Ohm's law. After the introduction of the term electrical resistance, it is obvious to use the quotient U/I to define resistance R.



Graphical representation of the measured values for the resistors R_1 and R_2





Evaluation (4/5)

PHYWE


What does the graph on the right show?

Ohm's law does apply to the incandescent lamp. The resistance of the incandescent lamp with an increase in current, which causes an increase in temperature. Obviously, the resistance of the incandescent lamp is temperature-dependent.

not

grows





Evaluation (5/5)

PHYWE

So what is Ohm's law?

$$U = R \cdot I$$

$$R = U \cdot I$$

$$R = I/U$$

$$U = R/I$$

What other insight was gained during the experiment?

Ohm's law is universally valid.

The resistance is not temperature-dependent.

Ohm's law only applies to components whose resistance is constant. R=constant is the validity condition for Ohm's law.



Slide	Score/Total
Slide 12: Proportionality of light bulb brightness and voltage	0/1
Slide 13: Proportionality between voltage and current	0/1
Slide 14: Measurement with incandescent lamp	0/2
Slide 15: Multiple tasks	0/2
Total score	0/6

