

Diffraction at a slit



If monochromatic light hits a narrow slit, intensity minima and maxima appear behind it on a screen, from which positions the slit width can be determined at a known wavelength.

Physics

Light & Optics

Diffraction & interference



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time

20 minutes

This content can also be found online at:



<http://localhost:1337/c/647c850ed59fde0002e27e64>

PHYWE

General information



Application

PHYWE



Experimental setup

When monochromatic light hits a narrow slit, an interference pattern with intensity maxima and minima appears behind it on a screen.

The slit width can be determined from their positions at a known wavelength.

Other teacher information (1/2)

PHYWE

Prior knowledge



To understand this experiment, students should already be familiar with the wave behaviour of light. For illustration purposes, it can be helpful to show interference of water waves beforehand.

Principle



A laser beam shines through an aperture with single slits and creates an interference pattern on a screen behind it.

The interference pattern can be used to read off the intensity minima and maxima and, if the wavelength is known, the original slit width can also be determined.

Other teacher information (2/2)

PHYWE

Learning objective



If we look at the interference pattern of a slit, we see that the maxima and minima are equally spaced.

In addition, the central intensity maximum is twice as wide as the other maxima. If one compares the interference patterns of the different slits with each other, it can be seen that the distances between the maxima and minima increase with decreasing slit width.

Tasks



- Observing the interference patterns on the screen.
- Marking the positions of the maxima and minima on the screen.
- Determination of the slit width of a slit.

Safety instructions


PHYWE



It is essential to ensure that you do not look directly into the laser beam.

The general instructions for safe experimentation in science lessons apply to this experiment.

Theory (1/3)



If a laser beam falls on a slit of width b , the beam in the slit area can be thought of as divided into two equal sub-beams, as shown on the left.

If after diffraction by the angle α the path difference $\Delta l/2$ between an edge and a central ray is an integer multiple of $\lambda/2$, then they interfere destructively. This path difference also exists between each ray of one half of the slit and a corresponding ray of the other half of the slit.


Theory (2/3)

There is always darkness in the interference pattern if the relation for Δl applies:

$$\Delta l = k * \lambda = b * \sin \alpha; k = \pm 1, \pm 2, \pm 3, \dots \quad (1)$$

If r is the distance between the slit and the sufficiently distant collecting screen S and if x_k is also the distance of the k th minimum from the center, then:

$$\sin \alpha_k = \frac{x_k}{\sqrt{x_k^2 + r^2}} \approx \frac{x_k}{r} \text{ for } x_k \ll r \quad (2)$$



Theory (3/3)

From (1) and (2) it follows for x_k :


$$x_k = k \frac{\lambda * r}{b} \quad (3)$$

Brightness maxima result for the angles

$$\sin \alpha_m = \frac{2m+1}{2} * \frac{\lambda}{b}; m = \pm 1, \pm 2, \pm 3, \dots \quad (4)$$

Accordingly, the following applies to x_m

$$x_m = m \frac{\lambda * r}{b} + \frac{\lambda * r}{2b} \quad (5)$$



Equipment

Position	Material	Item No.	Quantity
1	Optical profile-bench, l = 1000 mm	08370-00	1
2	Slide mount for optical bench	09822-00	3
3	Plate mount for three objects	09830-00	1
4	Diaphragm, 3 single slits	08522-00	1
5	Screen, metal, 300 x 300 mm	08062-00	1
6	Barrel base expert	02004-00	1
7	Measuring tape, l = 2 m	09936-00	1
8	Diodelaser, red, 1 mW, 635 nm	08761-99	1


Additional Equipment

PHYWE

Position	Equipment	Quantity
1	Sellotape	1
2	white sheet of paper	1

PHYWE

Set-up and Procedure



Set-up

PHYWE

The experimental setup is as shown in Fig. 1.

The tally marks of the tabs for holding the components have the following positions on the optical bench.

- Rider with diode laser at 2cm
- Rider with panel holder and inserted panel with gaps at 11cm

The barrel base with umbrella is located at a distance $r \geq 3\text{m}$ to the slit diaphragm.




Fig. 1

Procedure

A sheet of typewriter paper is attached to the screen with its surface normal pointing in the direction of the optical axis using adhesive tape. The aperture with slits is moved in the plate holder so that one slit at a time is completely irradiated by the laser light.

Use a water-soluble felt-tip pen to mark the positions of the maxima and minima of several diffraction orders. The distance r between the slit diaphragm and the screen must be determined with the tape measure. Use a ruler to determine with an accuracy of $0,5\text{mm}$ the distances $2x$ of the various intensity maxima and minima.




Procedure

PHYWE


A sheet of typewriter paper is attached to the screen with its surface normal pointing in the direction of the optical axis using adhesive tape. The aperture with slits is moved in the plate holder so that one slit at a time is completely irradiated by the laser light.

Use a water-soluble felt-tip pen to mark the positions of the maxima and minima of several diffraction orders. The distance r between the slit diaphragm and the screen must be determined with the tape measure. Use a ruler to determine with an accuracy of $0,5\text{mm}$ the distances $2x$ of the various intensity maxima and minima.



PHYWE

Evaluation



Evaluation (1/4)

$\pm k$	$2x_k/mm$
1	20,5
2	40,0
3	60,0
4	80,0
5	99,0

$\pm m$	$2x_m/mm$
1	30,5
2	49,5
3	70,0
4	89,5
5	111,0

If we look at the interference pattern of a slit, we see that the maxima and minima are equally spaced.

In addition, the central intensity maximum is twice as wide as the other maxima. Comparing the interference patterns of the different slits with each other, it can be seen that as the slit width decreases, the distances between the maxima and minima increase (Fig.2).




Fig. 2

Evaluation (2/4)

For a gap of the width $b = 0,2mm$ The distance values given in the table were used as examples. $2x_k$ resp. $2x_m$ and the corresponding value for the slit width was calculated with (3) or (5). From the data sheet of the diode laser, the value for the wavelength was taken as $\lambda = 635nm$ taken over.

The mean value for the gap width is obtained from the individual values:

$b = (0,203 \pm 0,002)mm$; $\Delta b/b \approx 1\%$ An inaccuracy of r at $\pm 5mm$ can be neglected when considering errors.

Calculated results from the sample data:

$\pm k$	b/mm
1	0,199
2	0,204
3	0,204
4	0,204
5	0,206

$\pm k$	b/mm
1	0,200
2	0,206
3	0,204
4	0,205
5	0,202

(Table: $r = 3205$ mm, gap $b = 0.2$ mm)

Evaluation (3/4)

PHYWE

What does this stand for b ?

Gap width

Wavelength

Gap centre distance

If we look at the interference pattern of a slit, we see that the maxima and minima are unequally spaced.

☐ True

☐ False

✓ Check

Evaluation (4/4)

PHYWE

What must be taken into account so that the formula gives the correct result?

The wavelength λ must not be smaller than 700nm.

All given values must be converted into the base units.

