


# Diffraction at an edge



If a plane wave meets a straight edge of an obstacle, the edge acts as the starting point of secondary waves according to Huygens' principle. These overlap with the wave part passing the edge undisturbed. Interference of both bundles results in a sequence of light and dark interference fringes in the passband.

Physics

Light &amp; Optics

Diffraction &amp; interference



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time

20 minutes

This content can also be found online at:



<http://localhost:1337/c/6157071cbcbea20003ba0ef5>


PHYWE



## General information

## Application

PHYWE



Test setup

If a plane wave meets a straight edge of an obstacle, the edge acts as the starting point of secondary waves according to Huygens' principle. These overlap with the part of the wave passing the edge undisturbed.

Interference of both beams results in a sequence of light and dark interference fringes in the passband.

The phenomenon of diffraction at the edge applies not only to light waves (transverse waves) but also to longitudinal waves (sound waves). Thus, when sound hits an edge, loudness minima and maxima can be observed behind this edge depending on the location.

## Other information (1/2)

PHYWE

### Previous



It should already be known, that diffraction and interference occur at openings as well as at obstacles. For this purpose it is useful to perform the experiment "Diffraction at a bar - Babinet's theorem" beforehand.

### Principle



If a wave meets an obstacle, the edge of this obstacle is the starting point of elementary waves according to Huygens' principle.

This circular wave emanating from the edge then interferes with the wave front flowing past the obstacle.

## Other information (2/2)

PHYWE

### Learning



Since edges can always be regarded as starting points of elementary waves, interference also takes place where a wave partly meets an obstacle and partly flows past it. Thus, only one edge is needed for interference.

### Tasks



- Observing the interference pattern
- Understanding why interference occurs at an edge
- Answering the questions in the protocol

## Safety instructions

PHYWE



It is imperative that you do not look directly into the laser beam.


The general instructions for safe experimentation in science lessons apply to this experiment.

## Theory

PHYWE

If a light bundle with plane wave fronts meets an edge, a part of the bundle is cut off while the other part can pass it unhindered. According to Huygens' principle, elementary waves emanate from the edge, which also propagate into the shadow area with decreasing intensity.

In the passband, these scattered waves overlap with the unobstructed plane waves. These two waves can then either cancel or amplify each other on the screen, depending on the path difference.



schematic propagation of the waves

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Optical profile-bench, l = 1000 mm</a>	08370-00	1
2	<a href="#">Slide mount for optical bench</a>	09822-00	2
3	<a href="#">Plate mount for three objects</a>	09830-00	1
4	<a href="#">Diaphragm with single slit, bar and edge</a>	08521-00	1
5	<a href="#">Lens, mounted, f +50 mm</a>	08020-01	1
6	<a href="#">Mount with scale on slide mount</a>	09823-00	1
7	<a href="#">Screen, metal, 300 x 300 mm</a>	08062-00	1
8	<a href="#">Barrel base expert</a>	02004-00	1
9	<a href="#">Diodelaser, red, 1 mW, 635 nm</a>	08761-99	1


## Equipment

PHYWE

Position	Equipment	Quantity
1	Adhesive tape	1
2	white sheet of paper	1


PHYWE

## Structure and implementation



## Set-up

PHYWE



Experimental setup

The experimental setup is as shown in the figure, with the tally marks of the tabs to hold the components on the optical bench to the following positions:

- Slide mount with diode laser at 2cm
- Mount with scale and inserted converging lens  $f = +50mm$  at 13cm
- Slide mount with plate holder and diaphragm at 25cm

The barrel base with screen is located at a distance of between 1.5 and 2 metres from the bar screen.

## Procedure

PHYWE




Set-up with screen

A white sheet of paper is attached to the screen with adhesive tape, and the point of impact of the unobstructed laser beam is marked with a felt-tip pen. A vertical line is drawn through the point of impact.

Afterwards, the mount with the inserted converging lens is mounted to expand the laser beam.

Finally, insert the diaphragm with bar into the plate holder and position it so that about 2/3 of the beam is covered by the diaphragm and the remaining part can hit the screen past the edge.

PHYWE




# Evaluation

## Evaluation (1/3)

PHYWE

The nearly parallel light beam of a laser  $L$  is directed through a converging lens ( $f = +50\text{ mm}$ ) is made divergent and then hits a straight edge of an obstacle.

On a screen  $S$  at a distance  $r$  on the edge  $K$  an enlarged image of bright and dark interference fringes can be observed in the free passband.




The very demanding theoretical treatment of diffraction at an edge must be omitted here.

However, it would show that the maximum is not at the geometric shadow boundary, but already in the passband.



## Evaluation (2/3)

PHYWE



Where is the brightest intensity maximum?

Perpendicular behind the edge.

Shifted a little into the shadows.

Shifted a little into the passband.

## Evaluation (3/3)

PHYWE

At which of these objects does interference occur?

☐ Polarizing Filter.

☐ Obstacles whose diameter is small enough.

☐ Openings whose diameter is small enough.

☐ Edges that cut off some of the light.

✓ Check


Which waves interfere in this experiment?

The waves cut off by the obstacle with the continuous wave front.

The continuous wave front with the elementary waves generated at the edge.

Slide	Score / Total
Slide 14: Intensity maximum	0/3
Slide 15: Multiple tasks	0/4

Total score  0/7

 Show solutions

 Repeat