


# The series motor (DEMO) with Cobra SMARTsense



Physics

Electricity &amp; Magnetism

Electromagnetism &amp; Induction

Physics

Electricity &amp; Magnetism

Electric generator, motor, transformer



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time

20 minutes

This content can also be found online at:



<http://localhost:1337/c/647641eba4c6900002803dd9>


PHYWE

# Teacher information



## Application

PHYWE



Experimental setup

Electric motors are installed in many machines. Be it the electric car or the electric toothbrush. In addition to a permanent magnet, an electric motor can also be operated with an electromagnet. If armature coils and field coils are connected in series, then it is a main shunt motor.

The properties of this motor are investigated by observing the direction of rotation and measuring the current. In this experiment, the principle of the main shunt motor is clarified.

## Other teacher information (1/2)

PHYWE

## Prior knowledge



No prior knowledge is required.

## Principle



The attraction and repulsion of magnetic fields causes the motor to rotate. The external magnetic field is generated by the coils connected in series. The T-armature also forms a magnetic field, which is reversed at the right time with the help of a commutator.

## Other teacher information (2/2)

PHYWE

## Learning objective



The students should understand how a main shunt motor works.

## Tasks



Investigate how a main shunt motor works with direct current.

PHYWE



## Student information

### Motivation

PHYWE

Electric motors are installed in many machines. Be it the electric car or the electric toothbrush. In addition to a permanent magnet, an electric motor can also be operated with an electromagnet. If armature coils and field coils are connected in series, then it is a main shunt motor.

The properties of this motor are investigated by observing the direction of rotation and measuring the current. In this experiment, the principle of the main shunt motor is clarified.



An electric car

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">PHYWE Power supply, universal, analog display DC: 18 V, 5 A / AC: 15 V, 5 A</a>	13503-93	1
2	<a href="#">Cobra SMARTsense High Current - Sensor for measuring electric current <math>\pm 10</math> A (Bluetooth + USB)</a>	12925-00	1
3	<a href="#">Bench clamp</a>	02012-00	1
4	<a href="#">Plate holder, opening width 2 - 35 mm</a>	06509-00	1
5	<a href="#">Iron core, U-shaped, laminated</a>	06501-00	1
6	<a href="#">Coil, 300 turns</a>	06513-01	2
7	<a href="#">Motor set</a>	06550-00	1
8	<a href="#">Rotor coil, Double-T armature</a>	06554-00	1
9	<a href="#">Cord pulley</a>	06558-01	1
10	<a href="#">Connecting cord, 32 A, 750 mm, red</a>	07362-01	3
11	<a href="#">Connecting cord, 32 A, 750 mm, blue</a>	07362-04	2
12	<a href="#">measureAPP - the free measurement software for all devices and operating systems</a>	14581-61	1

## Set-up (1/3)

PHYWE

For measurement with the **Cobra SMARTsense sensors** the **PHYWE measureAPP** is required. The app can be downloaded free of charge from the relevant app store (see below for QR codes). Before starting the app, please check that on your device (smartphone, tablet, desktop PC) **Bluetooth** is **activated**.



iOS



Android



Windows

## Set-up (2/3)

- Assemble the motor attachment according to Fig. 1.
- Push the axle [1] of the double T-anchor into the bearing hole [3] of the motor attachment and screw it tight with the cord washer [2].
- Place the abrasive brushes [4] of the motor attachment against the interrupted copper ring [7], tighten the knurled screws [5] slightly upwards so that the spring of the lever arms is tensioned. This presses the brushes firmly onto the copper ring. The electrical contact between the armature coils and the connection sockets [6] is established.




Fig. 1

## Setup (3/3)

PHYWE

- Complete the set-up according to Fig. 2 and Fig. 3.
- Clamp the iron core with holder in the bench clamp.
- Place the coils and motor attachment on the iron core.
- Set the DC voltage at the power supply unit to 0 V-.
- Connect field coils and armature coil in series and connect the motor to the power supply unit via the meter.




Fig. 2

## Procedure (1/2)

PHYWE



Cobra SMARTsense

- Switch on the SMARTsense sensor and make sure that the end device can connect to Bluetooth devices.
- Open the PHYWE measureApp and select the sensor "High Current".
- Select the sampling rate of your choice. The higher the sampling rate, the more accurate the measurement.

## Procedure (2/2)

PHYWE





Fig. 3

- Set the voltage to approx. 6 V-, you may have to "start" the motor by turning it.
- Change the voltage. Observe the speed and the recorded measurement signal.
- Set the voltage to 0 V-. Reconnect the operating voltage at the power supply unit. Increase the voltage and observe the direction of rotation.
- Set the voltage to 0 V-. Reverse the voltage at the terminals of the armature coil. Increase the voltage and observe the direction of rotation.
- Load the motor by pressing on the pulley with your finger. Observe the speed and the recorded measurement signal.

PHYWE

## Report





**Task (1/6)****PHYWE**

How does the speed of the motor and the amperage change as the voltage increases?

The speed changes little, the amperage increases.

The speed increases, the current changes little.

The speed and the current increase.

The speed and the amperage change little.

**Task (2/6)****PHYWE**

By reversing the polarity of the operating voltage...

... the direction of rotation remains constant.

... the engine stops.

... the direction of rotation changes.

**Task (3/6)****PHYWE**

If the direction of current only changes in the armature coil,...

... the engine stops turning.

... the direction of rotation remains constant.

... the direction of rotation changes.

**Task (4/6)****PHYWE**

With increased load...

... the speed of the motor decreases and the current increases.

... the speed of the motor increases and the current decreases.

... the speed of the motor decreases and the amperage decreases.

... the speed of the motor increases and the amperage increases.

## Task (5/6)

PHYWE

Drag the words into the correct boxes!

If an  is used to operate an , it must generate a sufficiently large  in the vicinity of the armature. Therefore, a U-shaped iron core with two field coils is used, between whose poles the armature runs. The  and field coils are connected in  in a main shunt motor (Fig. 3).

series

electric motor

electromagnet

armature coils

magnetic field

☒ Check

## Task (6/6)

PHYWE

Drag the words into the correct boxes!

When the polarity of the  is reversed, both the field of the armature coil and that of the  are reversed, so that the  remains the same. If, on the other hand, only the  of the armature coil changes, then only this magnetic field changes its  and thus also the direction of rotation.

current direction

direction

field coils

sense of rotation

operating voltage

☒ Check

Slide	Score / Total
Slide 14: Engine speed	0/1
Slide 15: Reversing the polarity of the operating voltage	0/1
Slide 16: Current direction of the armature coil	0/1
Slide 17: Behaviour under load	0/1
Slide 18: Main shunt motor	0/5
Slide 19: Reverse polarity behaviour	0/5

Total score  0 / 14

 Show solutions

 Repeat