


Generation of induced voltage with a permanent magnet (DEMO)



Physics

Electricity & Magnetism

Electromagnetism & Induction

Physics

Electricity & Magnetism

Electric generator, motor, transformer



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time

20 minutes

This content can also be found online at:



<http://localhost:1337/c/6492f42c3759e10002069160>


PHYWE



Teacher information

Application

PHYWE



Experimental setup

Electromagnetic induction (also Faraday induction, after Michael Faraday, induction for short) refers to the creation of an electric field when the magnetic flux changes.

In many cases, the electric field can be detected directly by measuring an electric voltage.

Other teacher information (1/2)

PHYWE

Prior knowledge



No prior knowledge is required.

Principle



The change in magnetic flux acting on an electrical conductor induces an electrical voltage and thus an electrical current flow in that conductor.

Other teacher information (2/2)

PHYWE

Learning objective




Students should understand the principle behind electromagnetic induction.

Tasks



Investigate how induction voltages can be generated with the help of a permanent magnet and which conditions influence their level.

PHYWE




Student information

Motivation

PHYWE

Electromagnetic induction (also Faraday induction, after Michael Faraday, induction for short) refers to the creation of an electric field when the magnetic flux changes.

In many cases, the electric field can be detected directly by measuring an electric voltage.



One coil

Equipment

Position	Material	Item No.	Quantity
1	Cobra SMARTsense Voltage - Sensor for measuring electrical voltage $\pm 30\text{ V}$ (Bluetooth + USB)	12901-01	1
2	Coil, 300 turns	06513-01	1
3	Coil, 1200 turns	06515-01	1
4	magnet, $l = 72\text{mm}$, rodshaped, colored poles	07823-00	1
5	Connecting cord, 32 A, 750 mm, red	07362-01	1
6	Connecting cord, 32 A, 750 mm, blue	07362-04	1
7	measureAPP - the free measurement software for all devices and operating systems	14581-61	1

Set-up (1/2)

PHYWE

For measurement with the **Cobra SMARTsense sensors** the **PHYWE measureAPP** is required. The app can be downloaded free of charge from the relevant app store (see below for QR codes). Before starting the app, please check that on your device (smartphone, tablet, desktop PC) **Bluetooth** is **activated**.



iOS



Android



Windows

Set-up (2/2)

- Set up the experiment according to Fig. 1.
- Carry out the following experimental steps one after the other and observe the recorded voltage signal in each case.
- Enter your observation in the prepared table in the report.




Fig. 1

Procedure (1/2)

PHYWE



Cobra SMARTsense


- Switch on the SMARTsense sensor and make sure that the end device can connect to Bluetooth devices.
- Open the PHYWE measureApp and select the sensor "Voltage".
- Select the sampling rate of your choice. The higher the sampling rate, the more accurate the measurement.

Procedure (2/2)

PHYWE

1. Move the magnet into the coil with the north pole first.
2. Move the magnet out of the coil again.
3. Move the magnet into the coil with the south pole first.
4. Move the magnet out of the coil again.
5. Move the magnet in and out of the coil faster.
6. Let the magnet rest in the coil.
7. Rotate the magnet in the coil around the longitudinal axis.
8. Replace the 400 turns coil with a 1600 turns coil and repeat the previous steps.

PHYWE



Report

Evaluation (1/3)

PHYWE

Movement	Observation	Movement	Observation
North pole into coil	<input type="text"/>	Faster movement of the magnet	<input type="text"/>
North pole out of coil	<input type="text"/>	Magnet rests in the coil	<input type="text"/>
South pole into coil	<input type="text"/>	Rotation of the magnet around the longitudinal axis	<input type="text"/>
South pole out of coil	<input type="text"/>	As 1. to 4. for coil with 1600 turns.	<input type="text"/>

Evaluation (2/3)

PHYWE

Drag the words into the correct boxes!

From the results in steps 1 to 6, it is clear that a [] is generated as long as the [] and the [] move relative to each other. However, as step 7 shows, the movement must be such that the [] encompassed by the coil changes in the process. Thus: A voltage is [] in a coil as long as the magnetic field enclosed by the coil changes.

magnetic field

magnet

induced

voltage

coil

☒ Check

Evaluation (3/3)

PHYWE

Drag the words into the correct boxes!

The [] of the induced voltage depends on whether the magnet is moving into or out of the coil and which [] of the magnet is facing the coil. The [] is higher the [] the movement and the higher the number of turns of the induction coil.

pole

induction voltage

faster

direction

☒ Check

Slide

Score / Total

Slide 14: Influence of the movement of the magnet

0/5

Slide 15: Induction voltage

0/4

Total score



Show solutions



Repeat



Export text