


Generation of induced voltage with an electromagnet (DEMO)



Physics

Electricity & Magnetism

Electromagnetism & Induction



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time

20 minutes

This content can also be found online at:



<http://localhost:1337/c/6478ff59d57c800002376ced>


PHYWE

Teacher information



Application

PHYWE



Experimental setup

Transformers are built into many electrical devices. In order to transform current well, the coils must not be too far apart and an iron core is important.

In this experiment, the induced voltage is measured as a function of the distance between the coils and as a function of the iron core.

Other teacher information (1/2)

PHYWE

Prior knowledge



No prior knowledge is required.

Principle



When alternating current flows through a coil, it creates a varying magnetic field, which in turn can induce a current in another coil.

Other teacher information (2/2)

PHYWE

Learning objective



Students should understand how a transformer works.

Tasks



Investigate the functioning of a transformer with and without an iron core.

Theory

PHYWE

The magnetic flux density B of a coil is:

$$B = \mu_r \cdot \mu_0 \cdot I \cdot \frac{n}{l}$$


μ_0 is the magnetic field constant in a vacuum, μ_r is the magnetic permeability, which depends on the material. $\frac{n}{l}$ describes the winding density. With alternating current, the magnetic field also changes, of course. The change in the magnetic field induces a current in the secondary coil.

Lenz's rule:

$$\oint_{\delta A} \vec{E} \cdot d\vec{s} = - \int A \frac{\delta B}{\delta t} \cdot d\vec{A}$$

PHYWE

Student information



Motivation

PHYWE

Transformers are built into many electrical devices. In order to transform current well, the coils must not be too far apart and an iron core is important.

In this experiment, the induced voltage is measured as a function of the distance between the coils and as a function of the iron core.



A transformer station

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Power supply, universal, analog display DC: 18 V, 5 A / AC: 15 V, 5 A	13503-93	1
2	PHYWE Demo Multimeter ADM 3: current, voltage, resistance, temperature	13840-00	1
3	Iron core, I-shaped, laminated, L=300mm	06504-01	1
4	Coil, 600 turns	06514-01	2
5	Two-way switch, single pole	06005-00	1
6	Connecting cord, 32 A, 750 mm, red	07362-01	1
7	Connecting cord, 32 A, 750 mm, blue	07362-04	1
8	Connecting cord, 32 A, 750 mm, black	07362-05	3

Set-up

PHYWE

- Set up the experiment according to Figure 1.
- In the primary circuit, connect the power supply unit, the changeover switch and the coil (as primary coil) in series.
- The secondary circuit consists only of the secondary coil and the demonstration multimeter.




Fig. 1

Procedure (1/2)

PHYWE

- Set the demonstration multimeter to AC voltage and select a suitable measuring range (e.g. 0...15 V~).
- Switch on the power supply and increase the primary voltage one step at a time (e.g. 6 V, 10 V and 15 V) with the switch closed. Observe the secondary voltage in each case.
- Select the primary voltage from $U_P = 15\text{V}$ and increase the distance between the coils step by step according to figure 2: e.g. in steps of 5 cm to approx. 20 cm. Observe the voltage display on the demonstration multimeter at each step.




Fig. 2

Procedure (2/2)

PHYWE




Fig. 3





Fig. 4

- Now push the coils together again in contact. Open the switches and push the iron core into the two coils until it is approximately in the middle (see Fig. 3). Close the switches and observe the secondary voltage for different primary voltages.
- Finally, increase again at a primary voltage of $U_P = 15\text{V}$ the distance between the two coils step by step as before (see Fig. 4) and observe the secondary voltage.

PHYWE

Report



Task (1/3)

PHYWE

The secondary voltage without iron core is...

... independent of the primary voltage.

... anti-proportional to the primary voltage.

... proportional to the primary voltage.

Task (2/3)

PHYWE

Drag the words into the correct boxes!

As shown in the experiment, it is possible to transport without
electrical in the form of wires. The magnetic
 drops rapidly outside the coil.

contact

electrical energy

flux density

✓ Check

Task (3/3)

PHYWE

Drag the words into the correct boxes!

However, the used has a much higher magnetic
and accordingly a much higher , which amplifies the magnetic field and
thus leads to a significantly increased induced in the secondary coil,
which is still clearly detectable even at a greater .

permeability

iron core

voltage

susceptibility

distance

 Check

Slide

Score / Total

Slide 13: Dependence between secondary and primary voltage

0/4

Slide 14: Functionality of the transformer

0/3

Slide 15: Properties of the iron core

0/5

Total score

  0/12 Show solutions Repeat

10/10