Current transformation (DEMO)





http://localhost:1337/c/6478fa50e6e35e00020403a7





Teacher information

Application

PHYWE



Experimental setup

Transformers are needed to transform input currents higher, lower or equal. Current transformers are usually built into the power supply unit of many electrical devices. In this experiment, the dependence of the number of turns and the current strength is investigated.

If a transformer is heavily loaded on the secondary side, the current strength in the secondary circuit depends on the primary current strength and the number of turns of the coils.

To investigate the regularity, different transformers are examined, first keeping the number of turns and current constant on the primary side and then on the secondary side.



www.phywe.de

Other teacher information (1/2)

PHYWE

Prior knowledge




No prior knowledge is required.

Principle

When alternating current flows through a coil, it creates a varying magnetic field, which in turn can induce a current in another coil.

In contrast to the voltages, the currents in a transformer are inversely proportional to the windings. Furthermore, if the primary current flows in a counterclockwise direction through the coil, the current in the secondary coil flows in exactly the opposite direction. By applying the flow theorem, the following applies:

 $I_1 \cdot n_1 = I_2 \cdot n_2$





3/10



Student information

Motivation

Transformers are built into many electrical devices. In order to transform current well, the coils must not be too far apart and an iron core is important.

This experiment investigates the relationship between current strength and number of turns.

Transformer station

PHYWE

PHYWE

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Variable transformer with digital display DC: 020 V, 12 A / AC: 025 V, 12 A	13542-93	1
2	PHYWE Demo Multimeter ADM 3: current, voltage, resistance, temperature	13840-00	2
3	Iron core, U-shaped, laminated	06501-00	1
4	Iron core, I-shaped, laminated	06500-00	1
5	Pins for iron cores, U-shaped	06502-00	1
6	Clamping device for iron cores	06506-00	1
7	Coil, 300 turns	06513-01	2
8	Coil, 600 turns	06514-01	1
9	Coil, 1200 turns	06515-01	1
10	Connecting cord, 32 A, 750 mm, black	07362-05	3
11	Connecting cord, 32 A, 750 mm, blue	07362-04	2



P1434205

Set-up

Set up the experiment according to Fig. 1.

- $\circ~$ Place two coils with 300 turns on the U-core.
- Insert the iron pins into the legs of the U-core and put on the yoke.
- Press the transformer firmly together using the clamping device.
- Connect the primary coil to the AC output of the regulating transformer via a measuring device. Select the measuring range 1 A~.
- $\circ~$ Connect the secondary coil to the second meter. Select the measuring range 1 A~.

Procedure (1/2)

Experiment 1

- $\circ\,$ Adjust the voltage at the regulating transformer so that the primary current is $I_1=1A.$
- $\circ\,$ Read the secondary current strength I_2 and enter it in Table 1 in the report.
- Change the transformer, insert the secondary coil with 600 turns and repeat the measurement. Choose a suitable measuring range.
- $\circ~$ Select the secondary coil with 1200 turns and repeat measurement.

6/10



Fig. 1

PHYWE

PHYWE

PHYWE

Procedure (2/2)

PHYWE


Experiment 2


- Build the transformer with primary coil and secondary coil of 300 turns each.
- $\circ\,$ Adjust the voltage at the regulating transformer so that the secondary current strength is $I_2=500mA$.
- $\circ\,$ Read the primary current I_1 and enter it in Table 2.
- Change the transformer, insert the primary coil with 600 turns and repeat the measurement. Choose a suitable measuring range.
- $\circ\,$ Select the primary coil with 1200 turns and repeat the measurement.



Report

www.phywe.de





PHYWE

Drag the words in	to the correct b	oxes!			
causes an		ber of turns n_2 o	n_1 and a constant print of the secondary coil a	-	ent I ₁ constant increase reduction
ask (4/5)					ЭНУУ
	to the correct b	oxes!			ЭНУУ
ASK (4/5) Drag the words inf Table 1 shows that	to the correct b		halved when the num	ıber of	PHYv primary current

Check

 I_2 .



secondary current

Γask (5/5)		ЭНУЖ
Drag the words into the co	orrect boxes!	
con The following relationship exi	n ₁ is doubled, the only has to be about half as large in order to keep the astant. ists between the current strengths and the number of ily loaded on the secondary side: in the primary and secondary circuits are inversely number of turns.	current strengths primary current number of turns secondary current
Slide		Score/Total
lide 14: Transformer summary		0/3
lide 15: Constant primary currer	nt	0/5
lide 16: Constant secondary curi	rent	0/4
	Total score	0/12

Show solutions

😂 Repeat

Export text



10/10