

The cell wall of the onion



Biology

Microscopy / Cell Biology

Basics of Microscopy & Work Technology

Biology

Microscopy / Cell Biology

Plants & Fungi

Biology

Microscopy / Cell Biology

Cell structure

Nature & technology

From the very small & the very big

Nature & technology

Plants & animals



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time


30 minutes

This content can also be found online at:



<http://localhost:1337/c/5f508d0537fe20003f10204>

PHYWE



Teacher information

Application

PHYWE



Cells of an onion (400x magnification)

With the naked eye we can see the shape of plants and animals and the diversity of these organisms very well. But what they are made up of remains closed to us without any aids. Hundreds of years ago, scientists already tried to explore the structure of organisms with lenses that allowed magnification. The term cell for a closed space was coined by the English scientist Robert Hooke. The year 1838 is regarded as the beginning of cytology. The botanist Matthias Schleiden recognized that plants consist exclusively of cells.

Other teacher information (1/3)

PHYWE

Prior knowledge



Robert Hooke used a simple microscope to find enclosed spaces in the cork of the cork oak and called them cells (lat. cellula = small chamber). The Dutchman Antonie van Leeuwenhoek was able to achieve a higher magnification with his microscope and discovered bacteria. In 1838 and 1839 Matthias Jacob Schleiden and Theodor Schwann determined the composition of all plants and animals from cells.

Scientific principle



Even a simple transmitted light microscope can be used to examine the structure of a cell, because it allows up to 400 times magnification of an object under examination (= micro-preparation).

Other teacher information (2/3)

Learning objective



The students learn to make a simple micro-preparation and view it under a microscope at different magnifications. A detailed drawing of the microscope image will be made and the cell components will be labeled.

Tasks



1. Production of a micropreparation with subsequent rapid staining
2. Creation of a drawing and labelling of the cell components

Other teacher information (2/3)

PHYWE

Learning objective



The students learn to make a simple micro-preparation and view it under a microscope at different magnifications. A detailed drawing of the microscope image will be made and the cell components will be labeled.

Tasks




1. Production of a micropreparation with subsequent rapid staining
2. Creation of a drawing and labelling of the cell components


Other teacher information (3/3)

Colouring


For this preparation, we recommend quick staining with methyl green and Lugol's solution, whereby the effects and advantages of both staining options can be compared.



Without staining (100x)



methyl green (100x)



Lugol solution (100x)

Safety instructions

PHYWE



- Methyl green is corrosive.
- Put on protective goggles!

H and P phrases for methyl green

H314: Causes severe skin burns and severe eye damage

H411: Toxic to aquatic organisms, with long-term effects.


P273: Avoid release into the environment.

P280: Wear protective gloves / protective clothing / eye protection / face protection.

P301 + P330 + P331: If swallowed: Rinse mouth. Do not induce vomiting.


PHYWE

Student Information



Motivation

PHYWE



Cells of an onion
(100x magnification)

In this experiment, you will learn in a simple way and with methods you already know how to create a micro-preparation of plant cells and to look at it under the microscope. With this you can identify the shape of a plant cell, the position of the cell nucleus and other cell components.

Equipment

Position	Material	Item No.	Quantity
1	PHYWE Binocular student microscope, 1000x, mechanical stage	MIC-129A	1
2	Microscopic slides, 50 pcs	64691-00	1
3	Cover glasses 18x18 mm, 50 pcs	64685-00	1
4	Beaker, 100 ml, plastic (PP)	36011-01	1
5	Dropping pipette with bulb, 10pcs	47131-01	1
6	Tweezers, straight, pointed, 120mm	64607-00	1
7	Scalpel holder	64615-00	1
8	Scalpel blades, rounded tip, 10 off	64615-02	1
9	Knife, stainless	33476-00	1
10	Chemicals set for TESS advanced Microscopy	13290-10	1


Equipment

PHYWE

Position	Material	Item No.	Quantity
1	PHYWE Binocular student microscope, 1000x, mechanical stage	MIC-129A	1
2	Microscopic slides, 50 pcs	64691-00	1
3	Cover glasses 18x18 mm, 50 pcs	64685-00	1
4	Beaker, 100 ml, plastic (PP)	36011-01	1
5	Dropping pipette with bulb, 10pcs	47131-01	1
6	Tweezers, straight, pointed, 120mm	64607-00	1
7	Scalpel holder	64615-00	1
8	Scalpel blades, rounded tip, 10 off	64615-02	1
9	Knife, stainless	33476-00	1
10	Chemicals set for TESS advanced Microscopy	13290-10	1

Procedure (1/4)

PHYWE




Make the preparation (1)

One onion is quartered and a single scale is removed.

Procedure (2/4)

Preparation of the preparation (2)

On the inside of a scale lies the thin, silvery onion skin. Make a 2 x 2 mm incision with the scalpel and place the preparation directly into a drop of water.



Procedure (3/4)

PHYWE




Microscopy

Microscopy with the lowest magnification (40x magnification). When you see a nice cell, you move this interesting spot right in the middle of the field of view. Now set the medium magnification (100x magnification) by turning the revolving nosepiece.

Procedure (4/4)

PHYWE




Staining and microscopy again

Pipette a drop of dye solution (methyl green and/or Lugol solution; a micropreparation must be prepared for each dye) next to the preparation and draw it under the cover glass using blotting paper. Microscope again!

PHYWE

Report



Task 1 + 2

PHYWE

Which cell components can you recognize in your micropreparation?

- ☐ Nuclei
- ☐ Mitochondria
- ☐ Cell walls
- ☐ Chloroplasts

☒ Check

How are the cells in your micropreparation arranged?

- ☐ Like a tiled floor
- ☐ Like the cells of a tumor
- ☐ Like a disorganized cluster of cells

☒ Check

Task 3

PHYWE

Which cell components are not found in plant cells?


- ☐ Endoplasmic reticulum
- ☐ Mitochondria
- ☐ Fat cells
- ☐ Golgi apparatus

☒ Check

Task 4

PHYWE

Magnification of a cell



✓ Check

Task 5

PHYWE

Draw a drawing of the stained specimen on a sheet of drawing paper at 100x magnification, i.e. several cells and their arrangement relative to each other. Label the cell wall and nucleus.

