


Plasma currents in the water plantain



Biology

Microscopy / Cell Biology

Plants & Fungi

Biology

Microscopy / Cell Biology

Cell structure

Biology

Plant Physiology / Botany

Physiology of plants



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

30 minutes

This content can also be found online at:



<http://localhost:1337/c/612896714839a1000394d73c>


PHYWE

Teacher information



Application

PHYWE



Water plant (Elodea spec.) (100x)

The cell is a reaction space in which a lot is produced and broken down, i.e. in which a lively metabolism takes place. The constant movement of the plasma in which these substances are dissolved ensures the rapid transport of substances.

Other teacher information (1/5)

PHYWE

Prior knowledge



Students should be familiar with the structure of plant cells. They should also be able to make a fresh preparation and be familiar with the use of a microscope.

Scientific Principle



To familiarize students with the effect of salt on plant cells.

Other teacher information (2/5)

PHYWE

Learning objective



The students should be able to recognise and describe the movement of the cytoplasm as rotating. E.g. the plasma flows around the vacuole in a circle, rotates clockwise or similar.

Tasks




Have students observe the plasma movements in an aquatic plant.

Other teacher information (3/5)

Notes on material procurement

The water plant (Elodea spec.) can either be found in standing water or you can buy it in the aquarium trade as a plant for the aquarium. It can be purchased there at a reasonable price.



The water plant (Elodea spec.)

Other teacher information (4/5)

Information on plasma currents

The cytoplasm of the cell is structured by a protein scaffold, also known as the cytoskeleton. It gives support to the animal cell in particular, which is not supported by the cell wall, and fixes the cell nucleus and some organelles. Transport processes, which can be observed as plasma movements, take place along the fibres of the cytoskeleton. In this process, the protein fibres act like rails on which the organelles are transported with the help of motor proteins.

Other teacher information (5/5)

Notes on implementation

Making the preparation Single-celled algae and somewhat larger sessile organisms also use this aquatic plant as a habitat. The preparations then do not look quite so transparent under the microscope. For this reason the terminal leaves are recommended. They are still fresh and not algal.

Microscopy The plasma currents can be observed very soon in healthy, fresh and non-algal plants. If not too many drawings have been made by the learning group concerned, the water pest cell can also be drawn. It is easy to see that a large area of the cell is apparently empty (vacuole) and the chloroplasts are only present in a narrow area (cytoplasm). However, the dynamics of the process can neither be reproduced in the drawing nor in the photo, but can be described in words.


Safety instructions

PHYWE



- Working with microscopes for too long can lead to physical discomfort (fatigue, headache, nausea), especially when students are untrained.
- To avoid accidents, scalpels should be checked for completeness at the beginning and end of the experiment.
- Microscopes are sensitive. During transport and handling, care should be taken to ensure that everything is done carefully and without rushing.
- The general instructions for safe experimentation in science lessons to be applied to this experiment.


PHYWE



Student Information

Motivation

PHYWE




Water plant (Elodea spec.) (100x)

The cell is a reaction space in which a lot is produced and broken down, i.e. in which a lively metabolism takes place. The constant movement of the plasma in which these substances are dissolved ensures the rapid transport of substances.

Tasks

PHYWE



Elodea spec. (400x)

1. Making the preparation
2. Microscopy

Equipment


Position	Material	Item No.	Quantity
1	PHYWE Binocular student microscope, 1000x, mechanical stage	MIC-129A	1
2	Microscopic slides, 50 pcs	64691-00	1
3	Cover glasses 18x18 mm, 50 pcs	64685-00	1
4	Beaker, 100 ml, plastic (PP)	36011-01	1
5	Dropping pipette with bulb, 10pcs	47131-01	1
6	Tweezers, straight, pointed, 120mm	64607-00	1

Procedure (1/2)

PHYWE

Making the preparation

- Prepare a slide with a drop of water.
- Pluck a leaf from the upper part of the plant. The older leaves are overgrown with algae on the outside and are therefore not well suitable.
- The leaf is placed directly into the water drop and covered with the cover glass.



Plucking a leaf from the upper part of the plant

Procedure (2/2)


Microscopy

The cell plasma (cytoplasm) is very pale without staining and therefore hardly observable. The largest space in the cell is occupied by the vacuole anyway. The narrow area of the plasma contains many green chloroplasts that seem to passively float along in the plasma. Thus, the movement of the plasma can be observed indirectly.

- Observe with ascending magnification.
- If possible, observe leaf areas in which there is only one cell layer. This is mostly in the edge area and on the middle rib of the leaflet.
- You have to give yourself time for the examination. The current is often only visible after a few minutes, when the metabolic processes are activated through the lighting.

PHYWE

Report



Task 1

PHYWE

Which statement is correct?

The largest space in the cell is occupied by the vacuole.

The cytoplasm is motionless.

The largest space in the cell is occupied by the chloroplasts.

The cell is divided into many rooms of equal size.

Task 2

PHYWE

Which of the following statements are correct?

- ☐ The cell is a reaction space in which a lot is produced and broken down. It takes place as metabolism.
- ☐ The movement of the plasma is best observed through the green chloroplasts.
- ☐ The cell is actually a dead space in which no biological processes take place.
- ☐ The cytoplasm has a deep red color.

✓ Check

Task 3

PHYWE

Drag the words to the right places.

The cytoplasm of the cell is structured by a protein scaffold, also known as the . Along the fibers of the cytoskeleton, transport processes take place that can be observed as . In this process, the protein fibers act like rails on which the are transported with the help of motor proteins.

cytoskeleton

plasma movements


organelles

✓ Check

Task 4

PHYWE

Draw the movement of chloroplasts in the cell.



Slide

Score/Total

Slide 16: The structure of the cell	0/1
Slide 17: The cytoplasm	0/2
Slide 18: The plasma movements	0/3

Total  0/6 Solutions Repeat