


# Skeletal muscle (Item No.: P1443001)

## Curricular Relevance



### Difficulty



Easy

### Preparation Time



10 Minutes

### Execution Time



30 Minutes

### Recommended Group Size



1 Student

### Additional Requirements:

- Physiological saline solution
- Beef
- Heart

### Experiment Variations:

### Keywords:

## Task and equipment

## Information for teachers

### Information

Your skeletal muscles make it possible that you can move, hold your body in an upright position when you are standing or sitting, or your face can smile. The muscles of all mammals possess a similar basic structure.

### Information on obtaining material

The study material can be easily obtained from a butcher's shop. The experiments described here were tested with beef and bovine heart. However, skeletal muscles from any other mammal species may also be experimented with, if it should be better available.

### Information on the structure of muscles

We distinguish striated musculature (skeletal musculature) and smooth musculature which is found in the intestines. Muscles are attached to bones by means of tendons. Muscles are composed of muscular fiber bundles which are surrounded by connective tissue and which, again, are composed of numerous muscle fibers. Muscle fibers have developed by the fusion of muscle cells and therefore contain many nuclei. Each muscle fiber consists of several hundred thin myofibrils. The striation recognizable under the microscope is induced by the arrangement of the proteins myosin and actin that are lined up like a string of beads. The striation effect is produced by the myofibrils lying next to each other in one and the same plane. Actin and myosin overlap somewhat and can slide over each other yet more, an action which induces the contraction of the muscle.



### Hazards

- Carmine acetic acid is highly corrosive!
- Put on protective glasses!

### Hazard and Precautionary statements

Carminic acid:

H314: Causes severe skin burns and eye damage.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P260: Do not breathe vapour.

P301 + P330 + P331: IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P302 + P352: IF ON SKIN: Wash with plenty of soap and water.


P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P309 + P310: IF exposed or if you feel unwell: Immediately call a POISON CENTER or doctor/physician.

## Information on how to perform the experiment

### ad 1 and 2: Preparation of muscle tissue

Muscle tissue preparation should proceed in physiological saline because otherwise the surrounding membrane would detach itself. You have to observe that the students produce very thin and transparent slides and, finally, that so much sodium-chloride solution is added that microscopy will be possible free of air bubbles. The striated structure should be explained to the students or, depending on their learning potential, they should elaborate on it by themselves.




Skeletal muscle, unstained, 400x

### ad 3 and 4: Examining the heart muscle

The tissue of the heart muscle is primarily darker than skeletal muscle. The striation is similarly well recognizable as is the case with skeletal muscle. Staining with carminic acid can be easily managed by the students and it enhances the structures.



Heart muscle, unstained, 400x




Heart muscle, stained with carminic acid, 400x

# Skeletal muscle (Item No.: P1443001)

## Task and equipment

### Task

Skeletal muscle is also referred to as striated (striped) muscle. You will explore how this name may be explained.



### Equipment

Position No.	Material	Order No.	Quantity
1	Euromex BioBlue BB.4250 microscope	EUR-BB-4250	1
2	Microscopic slides, 50 pcs	64691-00	1
3	Cover glasses 18x18 mm, 50 pcs.	64685-00	1
4	Scissors, straight, pointed, l 110mm	64623-00	1
5	Dissecting needle, pointed	64620-00	1
6	Dissecting needle, lancet-shaped	64621-00	1
7	Tweezers, straight, pointed, 120mm	64607-00	1
8	Chemicals set for TESS advanced Microscopy	13290-10	1

## Set-up and procedure

### Hazards

- Carmine acetic acid is highly corrosive!
- Put on protective glasses!




### Information

Your skeletal muscles make it possible that you can move, hold your body in an upright position when you are standing or sitting, or your face can smile. The muscles of all mammals possess a similar basic structure.

### Methods and observations

#### 1. Preparation of a muscle

a) Pipette two drops of physiological saline onto a slide (see P1440901)




b) Cut a small piece of meat alongside the muscle fibers using the dissecting scissors. Tear the meat in the salt solution until you obtain very thin fibers.



#### 2. Microscopical examination

If necessary, add another drop of physiological saline before you start with microscopy.



View the specimen in three steps up to the highest power.  
Describe the observations you make in the report.

## **Additional tasks**

### 3. Heart musculature

Examine a piece of heart muscle in the same way.

### 4. Staining with carmine acetic acid

Proceed as described under Task 1, then apply some carmine acetic acid to the mounted slide and view the specimen under the microscope after 3 minutes.

Describe your observations ensuing Tasks 2 and 3 in the report.

## Report: Skeletal muscle

### Result - Observations 1

View the specimen in three steps up to the highest power.

Describe the observations you make! Explain why these muscles are called "striped".

---

---

---

---

### Result - Observations 2

Note your observations when examining a piece of heart muscle in the same way.

---

---

---

---

## Result - Observations 3

Describe your observations ensuing staining.

.....

.....

.....

.....