

Joule-Thomson effect



P2320600

Physics

Thermodynamics

Temperature & Heat



Difficulty level

hard



Group size

-



Preparation time

10 minutes




Execution time

30 minutes

This content can also be found online at:


<http://localhost:1337/c/6475b4cadd89390002768e25>

PHYWE



Informations

Introduction

PHYWE




Fig. 1: Experimental set-up (without PC)

A stream of gas is fed to a throttling point, where the gas (CO_2 or N_2) undergoes adiabatic expansion. The differences in temperature established between the two sides of the throttle point are measured at various pressures and the Joule-Thomson coefficients of the gases in question are calculated. The principle of the Joule-Thomson effect is also applied, for example, in snow cannons for freezing water.

Other information (1/2)

PHYWE

Prior knowledge



The students must be familiar with units such as pressure, temperature, mass and volume and be able to perform calculations using them. They must also be familiar with general good laboratory practice and general laboratory safety regulations.

Principle



The Joule-Thomson effect has various applications, particularly in gas processing and cooling:

- Natural gas purification by condensin impurities after cooling; Cooling of natural gas directly at the extraction site; Production of highly pure gases such as oxygen and nitrogen through air seperation; Air cooling in air conditioners and refrigerators; Production of liquid nitrogen in cryogenics

The Joule-Thomson effect is most effective in gases with high inversion temperatures, as the cool upon expansion at room temperature.

Other information (2/2)

PHYWE

Learning Goal



In this experiment, students become familiar with the different behaviors of gases and deepen their knowledge of physical equations. By going through the different parts with the same volume of air while changing the external influences, they learn about the relationships between pressure, temperature and volume. It is a simple introduction to thermodynamics.

Tasks



The experiment aims to impart students with a deeper understanding of the thermodynamic principles and technical applications of the Joule-Thomson effect through practical experience and measurements.

Determination of the Joule-Thomson coefficient of CO₂ and N₂

Safety instructions 1/2

PHYWE



Fig. 2: safety instructions

Carbon dioxide

H280: Contains gas under pressure; may explode if heated.

P403: Store in a well ventilated place.


Nitrogen

H280: Contains gas under pressure; may explode if heated.

P403: Store in a well ventilated place.

PHYWE

Setup and Procedure



Safety instructions 2/2

PHYWE

The set-up of the experiment is shown in Fig. 1.

If necessary, screw the reducing valves onto the steel cylinders and check the tightness of the main valves. Secure the steel cylinders in their location. Attach the PVC-tubing between the reducing valve and the Joule-Thomson apparatus with hose tube clips. On each side of the glass cylinder, introduce a temperature probe up to a few millimetres from the frit and attach with the union nut. Connect the temperature probe on the pressure side to inlet 1 and the temperature probe on the unpressurised side to inlet 2 of the temperature measurement apparatus.

Safety instructions 2/2

PHYWE

Important

The experimenting room and the experimental apparatus must be in a thermal equilibrium at the start of the measurement. The experimental apparatus should be kept out of direct sunlight and other sources of heating or cooling. Set the temperature measurement apparatus at **temperature difference** measurement. The temperature meter should be switched on at least 30 min before performing the experiment to avoid thermal drift.

Safety instructions 2/2

PHYWE

Important

Read the operating instructions for further explanations of the temperature meter. Open the valves in the following order: steel cylinder valve, operating valve, reducing valve, so that an initial pressure of 100 kPa is established. Reduce the pressure to zero in stages, in each case reading off the temperature difference one minute after the particular pressure has been established. Perform the measurement for both gases, and determine the atmospheric pressure and ambient temperature.

Ausrüstung

Position	Material	Item No.	Quantity
1	Joule-Thomson apparatus	04361-00	1
2	Universal power supply, 600mA 3/4.5/5/6/7.5/9/12V, incl. 9 adaptors	11078-99	1
3	Reducing valve for CO ₂ / He	33481-00	1
4	Reducing valve f.nitrogen	33483-00	1
5	Wrench for steel cylinders	40322-00	1
6	Gas-cylinder Trolley for 2 Cyl.	41790-20	1
7	Steel cylinder,CO ₂ , 10l, full	41761-00	1
8	Steel cylinder,nitrogen,10l, full	41763-00	1
9	Hose clip f.12-20 diameter tube	40995-00	2
10	PVC tubing webbed, inner dia. = 8 mm, l = 1 m	47528-01	2
11	Cobra SMARTsense High Precision Temperature - Sensor for measuring temperature, -50 ... +150 °C (Bluetooth + USB)	12950-00	1
12	measureLAB, multi-user license	14580-61	1

Zusätzliche Ausrüstung

PHYWE

Position	Material	Menge
1	PC with Windows XP® oder higher	1

Theory and evaluation

PHYWE




Fig. 2: Throttling and the Joule-Thomson effect

The state for real gases is given by the van der Waals equation

$$\left(p + \frac{a}{V^2}\right) * (V - b) = RT$$

where p is the pressure, V the molar volume and T the temperature of the gas.

R is the universal gas constant, a and b are the

characteristic Van der Waals coefficients of the gas. The additional

pressure by intermolecular forces of attraction is described by a , b

represents the volume of molecules

Theory and evaluation

PHYWE

In real gases the intrinsic energy U is composed of a thermocinetic and a potential component.

The total change of the intrinsic energy U of a real gas therefore depends not only on the temperature and the molar heat C_V of a gas but also on the volume. The potential energy $-a/V$ in the gas is given by the work against the intermolecular forces

$$U = C_v * T - \frac{a}{V}$$

Theory of evaluation

PHYWE

The total didfferential of the intrinsic energy is given by

$$dU = \left(\frac{\partial U}{\partial T} \right)_V \cdot dT + \left(\frac{\partial U}{\partial V} \right)_T \cdot dV$$

with

$$\left(\frac{\partial U}{\partial T} \right) = C_V \text{ for } \Delta V = 0$$

and

$$\left(\frac{\partial U}{\partial V} \right) = +\frac{a}{V^2} \text{ for } \Delta T = 0 .$$

Theory of evaluation

PHYWE

The expansion of the gas at the throttle-point is adiabatic ($\Delta Q=0$). If external heat losses and friction during the flow of the gas are excluded, the total energy H of the process is constant:

$$H = U + p \cdot V = U_1 + p_1 \cdot V_1 = U_2 + p_2 \cdot V_2$$

see Fig. 2

For further calculations, pressure p can be substituted by using van der Waals equation:

$$H = C_V \cdot T - \frac{a}{V} + \left(\frac{RT}{V-b} - \frac{a}{V^2} \right) \cdot V$$

Theory of evaluation

PHYWE

After calculating the differential of (1) and substituting the total differentials by finite differences again, one obtains a relation between ΔT and ΔV :

$$\Delta H = C_V \cdot \Delta T + \frac{2a}{V^2} \cdot \Delta V + \frac{RV}{V-b} \cdot \Delta T - \frac{RTV}{(V-b)^2} \Delta V + \frac{RT}{V-b} \cdot \Delta V$$

$$\Delta H = \left(C_V + \frac{RV}{V-b} \right) \cdot \Delta T + \left(\frac{2a}{V^2} - \frac{RTb}{(V-b)^2} \right) \cdot \Delta V$$

$\Delta H = 0$ because H is constant.

Theory of evaluation

PHYWE

In the next step, b is considered to be small compared to V . The molar heat at constant volume C_V can be replaced with the molar heat at constant pressure C_p via the ideal gas equation and under the assumption of even gas particle masses m :

$$C_{p,m} = C_{V,m} + R$$

Putting (3) into equation (2) results in

$$\frac{\Delta T}{\Delta V} = \frac{RTb - 2a}{V^2 \cdot C_{p,m}}$$

Theory of evaluation

PHYWE

At the throttle point a pressure gradient $p_1 - p_2$ and a temperature gradient $T_1 - T_2$ are established. This effect is named the Joule-Thomson effect and is described by the coefficient

$$\mu = \frac{T_1 - T_2}{p_1 - p_2} = \frac{\Delta T}{\Delta p}$$

To derive a formula of this coefficient with Van der Waals forces involved, the difference ΔV in (4) has to be expressed in terms of Δp . This can be approximately done via the ideal gas equation again under the assumption of even gas particle masses m :

$$V = \frac{nRT}{p}$$

Theory of evaluation

PHYWE

After differentiating this equation to get an expression for dV , the total differentials are substituted by the finite differences, which results in the following:

$$\Delta V = \frac{RT}{p^2} \Delta p$$

Solving equation (4) for ΔT gives

$$\Delta T = \frac{\Delta V}{V^2 \cdot C_{p,m}} \cdot (RTb - 2a)$$

After replacing ΔV in (7) with the term in (6)

$$\Delta T = \frac{-RT\Delta p}{p^2} \cdot \frac{p^2}{R^2 \cdot T^2} \cdot \frac{RTb - 2a}{C_{p,m}}$$

Theory of evaluation

PHYWE

Putting this into (4), the equation for the Joule-Thomson coefficient results in

$$\mu = \frac{\Delta T}{\Delta p} = \frac{1}{C_{p,m}} \cdot \left(\frac{2a}{RT} - b \right)$$

Safety instructions 1/2

PHYWE




Fig. 3 shows the temperature differences measured at various pressures. The measurement values in Fig. 3 give the straight line gradients.

$$\mu_{\text{CO}_2} = (1.084 \pm 0.050) \cdot 10^{-5} \frac{\text{K}}{\text{Pa}}$$

and

$$\mu_{\text{N}_2} = (0.253 \pm 0.030) \cdot 10^{-5} \frac{\text{K}}{\text{Pa}}$$

where N_2 is an approximation for air.

Theory of evaluation

PHYWE

The two temperature probes may give different absolute values for the same temperature. But this is no problem, as only the temperature difference is important for the determination of the Joule-Thomson coefficients.

The literature values are

$$\mu_{\text{CO}_2} = 1.16 \cdot 10^{-5} \frac{\text{K}}{\text{Pa}}$$

at 20 °C and 10^{-5} Pa and

$$\mu_{\text{air}} = 0.23 \cdot 10^{-5} \frac{\text{K}}{\text{Pa}}$$

at 20 °C and 10^{-5} Pa.

Theory of evaluation

PHYWE

For CO₂, with

$$a = 3.60 \text{ Pa m}^6 / \text{mol}^2$$

$$b = 42.7 \text{ cm}^3 / \text{mol}$$

$$c_p = 366.1 \text{ J/mol K}$$

equation (8) gives the coefficient

$$\mu_{\text{CO}_2} = 0.795 \cdot 10^{-5} \frac{\text{K}}{\text{Pa}}$$

Theory of evaluation

PHYWE

For air, with

$$a = 40 \text{ Pa m}^6 / \text{mol}^2$$

$$b = 39,1 \text{ cm}^3 / \text{mol}$$

$$c_p = 9 \text{ J/mol K}$$

equation (8) gives the coefficient

$$\mu_{\text{air}} = 0.387 \cdot 10^{-5} \frac{\text{K}}{\text{Pa}}$$

Theory of evaluation

PHYWE

Remarks

The formula for the Joule-Thomson coefficient indicates the condition for a cooling process

$$T_i = \frac{2a}{R \cdot b}$$

where T_i is the inversion temperature. If $\Delta p < 0$ the temperature T has to be lower than T_i to cool the gas with the process:

$\Delta p < 0$ and $T < T_i$ then $\Delta T < 0$.

For air (N_2) or CO_2 cooling is observed at room temperature. H_2 or He have a positive coefficient at room temperature, for H_2 this can be hazardous because self ignition is possible.

Measurement 1/2

PHYWE




Fig.3: measureLAB software interface

- start measureLAB Software .
- load experiment P2320600.

Measurement 2/2

PHYWE




Fig.4: measureLAB software interface with the high precision sensor and the two temperature sensors

- start the sensor and press it for at least 3 seconds.
- the sensor will be displayed on the left side of the measureLAB control panel. Please checkmark it.
- start the measurement using this button 
- end the measurement using this button. 