


Fine structure, one-electron and two-electron spectra

(Item No.: P2510600)

Curricular Relevance



Difficulty



Difficult

Preparation Time



1 Hour

Execution Time



2 Hours

Recommended Group Size



2 Students

Additional Requirements:

Experiment Variations:

Keywords:

Diffraction spectrometer, spin, angular momentum, spin-orbital angular momentum interaction, multiplicity, energy level, excitation energy, selection rules, doublets, parahelium, orthohelium, exchange energy, singlet series, triplet series, forbidden transitions

Overview

Short description

Principle

A level energy calculation with the Schrödinger equation predicts a level energy degeneracy for all electron states inside an atom belonging to the same shell having main quantum number n , and same orbit with orbital angular momentum quantum number l . This ignores the effects of magnetic spin-orbit interactions and relativistic effects which shift the energy of levels in atoms causing the fine structure thus lifting this level degeneracy.

By observation of the wavelength of light from transitions between atomic levels this energy shift can be determined. To determine the wavelength, a spectrometer-goniometer is used. It consists of a dispersive element, which deflects light through an angle depending on the wavelength, and a goniometer to determine the angle of deflection of the light beam. As dispersive elements a grating and a flint prism are used.



Fig. 1: Set up of the grating spectrometer

Materials

Position No.	Material	Order No.	Quantity
1	PHYWE Spectrometer-Goniometer with vernier SG1 PRO	35635-03	1
2	Diffraction grating, 600 lines/mm	08546-00	1
3	Prism, 60 degrees, height = 36.4 mm, flint	08237-00	1
4	Spectral lamp He, pico 9 base	08120-03	1
5	Spectral lamp Na, pico 9 base	08120-07	1
6	Spectral lamp Hg 100, pico 9 base	08120-14	1
7	Spectral lamp Cd, pico 9 base	08120-01	1
8	Spectral lamp Zn, pico 9 base	08120-11	1
9	Magnifier with handle, 4x, d=50 mm	87004-03	1
10	Power supply for spectral lamps	13662-97	1
11	Lamp holder pico 9 for spectral lamps	08119-00	1
12	Tripod base PHYWE	02002-55	1

Tasks

One-electron spectra:

1. Set up the spectrometer as a grating spectrometer and determine the grating constant using literature values for the visible spectral lines from a helium lamp.
2. With the determined grating constant measure the wavelengths of the observable sodium spectral lines belonging to the split 3^2P level, determine the energy difference between the $3^2P_{1/2}$ and $3^2P_{3/2}$ level

Two-electron spectra:

1. Set up the spectrometer as a prism spectrometer and determine the dispersion function of the flint prism using literature values for the visible spectral lines from a helium lamp.
2. Determine the wavelengths of spectral lines of zinc, cadmium, and mercury. Observe the increase of triplet-singlet transitions with increasing atomic number.

General remarks:

Darken the room in order to see fainter lines of the spectral lamps and to perceive lines on the red or blue end of the visible spectrum. It is advisable to have a separate lamp for reading the angle scales that can be easily switched off when observing the spectral lines through the ocular telescope.

Align the grating or the prism on the spectrometer table according to the operating instructions of the spectrometer-goniometer 35635-03. Align the slit vertical and the reticule parallel to the slit image. For the grating be sure to check that the vertical position of the slit image is equal on both sides of the undeflected beam as else deviations from verticality of the slit produce angle errors. Also for the prism the slit image vertical position should be constant throughout the whole spectrum. Be aware that because of the dispersion of the lenses the slit image has to be focused for every spectral line.

To avoid lateral play of the focussing mechanics after exerting lateral force to slit or eyepiece, turning the focussing knob back and forth brings the movable part to its reproducible neutral position.

measured angle left			right			calculated grating const.		
8	18	15.30	337	39	15.35	438.92	601.2	603.1
8	37	15.62	337	20	15.67	447.15	602.0	603.9
9	27	16.45	336	27	16.55	471.31	600.8	604.4
10	14	17.23	335	41	17.32	492.19	601.9	604.8
10	34	17.57	335	21	17.65	501.71	601.6	604.3
10	40	17.67	335	15	17.75	504.91	601.1	603.8
13	41	20.68	332	14	20.77	587.73	601.0	603.3
16	44	23.73	329	15	23.75	668	602.5	602.9
average:								602.7

Table 3: Measurement example for He, angle calculated with zero order position at 353.00°

The NIST wavelengths are used to calculate the grating constant.

With the grating constant of Table 3 the sodium lines evaluate as Table 4:

measured angle left			right			calculated wavelength			
<i>1st order</i>									
10	25	17.42	335	30	17.50	496.7	499.0	497.8	497.9
10	26	17.43	335	29	17.52	497.1	499.4	498.3	498.3
11	2	18.03	334	53	18.12	513.7	516.0	514.8	514.9
11	3	18.05	334	52	18.13	514.1	516.4	515.3	515.3
13	0	20.00	332	57	20.05	567.5	568.9	568.2	568.3
13	1	20.02	332	56	20.07	568.0	569.3	568.7	568.8
13	48	20.80	332	11	20.82	589.2	589.7	589.5	589.0
13	50	20.83	332	10	20.83	590.1	590.1	590.1	589.6
14	45	21.75	331	12	21.80	614.9	616.2	615.5	615.4
14	46	21.77	331	11	21.82	615.3	616.7	616.0	616.1
<i>2nd order</i>									
			316	8	36.87		497.8		497.9
			316	6	36.90		498.1		498.3
			314	40	38.33		514.6		514.9
			314	39	38.35		514.8		515.3
36	20	43.33	309	50	43.17	569.3	567.6	568.5	568.3
36	22	43.37	309	47	43.22	569.7	568.1	568.9	568.8
38	21	45.35	307	52	45.13	590.2	588.0	589.1	589.0
38	25	45.42	307	49	45.18	590.9	588.5	589.7	589.6
			305	12	47.80		614.6		615.4
			305	9	47.85		615.1		616.1

Table 4: Measured angles with the sodium lamp, calculated wavelengths with grating constant from Table 3 and NIST data for comparison

Evaluation

In atoms with more than one electron the $n - l$ degeneracy is already lifted because the electrons see a screened potential of the atomic nucleus because of the other electrons. The resulting deviation of the central potential from the $1/r$ -law lifts the level energies of orbits with higher l .

For single electrons outside completely filled orbits it can be assumed that there is a magnetic interaction energy between the spin magnetic moment and the orbital magnetic moment. For orbits with zero angular momentum (S -orbits) the orbital magnetic moment is zero. For P -orbits the e.g. z -component of spin and orbit magnetic moments can be in the same direction, so there is a positive interaction energy increasing the level energy - the magnetic fields of both moments add and the outer field is larger, so the energy stored in magnetic field integrated over the whole space is larger. Else the spin and orbit magnetic moments can be opposed, so the magnetic field weakened and the level energy lowered. Thus the energy of the ${}^2P_{3/2}$ level should be expected to be higher than of the ${}^2P_{1/2}$ level. For all orbits with an orbital magnetic moment this situation is the same and because radiative transitions are allowed from both upper and lower sub-level, the spectrum appears as doublet system.

Since the electron spin magnetic moment is fixed, the size of the spin orbit interaction depends on the magnetic moment of the orbit. For a magnetic field the "speed" of the charge is essential, that is proportional to the square root of the kinetic energy. For orbits with lower angular momentum thus lower l the average kinetic energy is higher, so the corresponding magnetic moments can be stronger though the angular momentum is smaller. Even more so in atoms with several electrons where the inner electrons screen the electric field of the atomic nucleus - they deform the central potential such that the potential is lower at the outside.

This causes in sodium the splitting of the 2D -levels to be far smaller than of the 2P -levels.

The technique to calculate this is perturbation theory taking into account the screening of the electric central potential by the filled inner electron shells. There are many different possible formalisms here so for the calculation refer to literature, starting e.g. with Wikipedia,

https://en.wikipedia.org/wiki/Spin%E2%80%93orbit_interaction

https://en.wikipedia.org/wiki/Fine_structure

From the measured wavelength λ data of Table 4 the wavelength splitting can be calculated to an energy splitting assuming only the $3\,{}^2P_{1/2}$ and $3\,{}^2P_{3/2}$ levels differ in energy of the levels involved with the transition.

For calculation of energies E the wavelength in air is to be corrected for air refractive index of $n = 1.00028$ as is done by NIST in their data, so with speed of light c and Planck constant h the energy is $E = h \cdot c / (n \cdot \lambda)$ and to be divided by electron charge e for results in eV. Table 5 shows the result. The result of 1.8 meV is to be compared with NIST value of 2.13202 meV for the splitting of Na $3\,{}^2P_{1/2}$ and $3\,{}^2P_{3/2}$ level. The large deviation here is due to the angle precision of $0^\circ 1'$ corresponding to an energy of 0.6 meV at 589 nm in second order and 1.6 meV at 589 nm in first order.

λ [nm]	ΔE [meV]
497.8	
498.3	2.3
514.8	
515.3	2.1
568.2	
568.7	1.7
589.5	
590.1	2.4
615.5	
616.0	1.5
497.8	
498.1	1.9
514.6	
514.8	0.9
568.5	
568.9	1.7
589.1	
589.7	2.1
614.6	
615.1	1.6
Av. [meV]:	1.8

Table 5: Measured energy difference between the Na $3\,{}^2P_{1/2}$ and $3\,{}^2P_{3/2}$ level

Fine structure of two electron spectra - prism spectrometer

Set up and procedure

a) Determination of the dispersion function

Insert the helium (He 10) lamp into the lamp holder, connect this to the power supply and turn on.

Align the gonio-spectrometer according to the operating instructions and record the reticule angle where the slit image appears without prism.

Mount and align the diffraction prism according to the operating instructions.

Measure the minimum angle of deviation for visible He spectral lines reading the reticule angle position on the outer scale and calculating the angle difference with the reticule angle without prism. The table angle data read on the inner scale are far less precise and needn't be used for evaluation because at the angle of least deflection also the change of deflection angle with change of table angle is smallest. But they can help to check if there are errors like unintentional change of prism position or badly fixed table and the like. The table angle change should be half the ocular angle change.

Use a spread sheet to plot the angle of least deflection over the wavelength and to identify the He spectral lines, refer to Table 1. You can assume a smooth dispersion curve to identify the lines. Use this plot as calibration curve for wavelength measurement of the spectral lines of the other lamps.

Calculate the refractive index of the prism for each wavelength with the spread sheet and plot the dispersion curve from the prism in use.

b) Determination of two-electron spectral line wavelengths

Turn off the power supply and let cool down the lamp, then exchange the He lamp with the zinc (Zn 10), cadmium (Cd 10) and mercury (Hg 100) lamp. For each lamp record the angle of least deflection for the visible lines. Determine the wavelength from the angle using the calibration curve.

Compare the results with the theoretical values.



Fig. 4: Set up as prism spectrometer

Data and results

If a prism with inner angle α_{prism} at its corner deflects a light beam of wavelength λ at the angle of least deflection through the angle β_λ , then the refractive index of the prism n_λ for this light is according to Snell's Law

$$n_\lambda = \frac{\sin((\alpha_{prism} + \beta_\lambda)/2)}{\sin(\alpha_{prism}/2)}$$

A measurement example is shown in Table 6. For the helium spectral line and transition data refer to the data of Table 1 and Fig. 2.

deg.	min.	β_λ [°]	λ [nm]	n_λ
286	38	66.30	396.47	1.784
287	37	65.32	402.62	1.777
288	20	64.60	412.08	1.771
288	27	64.48	414.49	1.770
289	35	63.35	438.92	1.761
289	57	62.98	447.15	1.757
290	49	62.12	471.31	1.750
291	28	61.47	492.19	1.745
291	43	61.22	501.71	1.743
291	47	61.15	504.91	1.742
293	25	59.52	587.73	1.728
294	20	58.60	668.00	1.720
294	42	58.23	706.72	1.716
294	51	58.08	728.34	1.715

Table 6: Minimum deviation angle for He spectral lines with undeflected light at angle of 352.93°






Fig. 5: Wavelength over angle of minimum deviation β_λ and inserted spline, data of Table 6

A spline may be used for interpolation of the wavelength over angle data as it can be assumed that the function is a smooth convex one.

Then this diagram can be used to read off the wavelength from the measured angle data for the visible spectral lines of the other spectral lamps.

Tables 7 through 9 and Fig. 7 through 9 show the expected data for zinc, cadmium and mercury.



colour	λ [nm]	upper (1)	lower (2)	A [s^{-1}]	E_1 [eV]	E_2 [eV]
blue	462.98	$3d^{10}4s5d, ^1D_2$	$3d^{10}4s4p, ^1P_1$	1.4×10^6	8.4729	5.7957
blue	468.01	$3d^{10}4s5s, ^3S_1$	$3d^{10}4s4p, ^3P_0$	1.4×10^7	6.6545	4.0061
blue	472.22	$3d^{10}4s5s, ^3S_1$	$3d^{10}4s4p, ^3P_1$	1.9×10^7	6.6545	4.0297
blue	481.05	$3d^{10}4s5s, ^3S_1$	$3d^{10}4s4p, ^3P_2$	1.9×10^7	6.6545	4.0779
green	507.0	$3d^{10}4s9p, ^3P_{0,1,2}$	$3d^{10}4s5s, ^3S_1$	5×10^5	9.0998	6.6545
green	518.20	$3d^{10}4s6s, ^1S_0$	$3d^{10}4s4p, ^1P_1$	1×10^7	8.1876	5.7957
green	530.9	$3d^{10}4s8p, ^3P_{0,1,2}$	$3d^{10}4s5s, ^3S_1$	1×10^6	8.9890	6.6545
red	636.23	$3d^{10}4s4d, ^1D_2$	$3d^{10}4s4p, ^1P_1$	4.7×10^7	7.7439	5.7957

Fig. 7: Zn transitions - NIST data




Fig. 7: Cd Grotrian diagram

colour	λ [nm]	upper (1)	lower (2)	A [s^{-1}]	E_1 [eV]	E_2 [eV]
UV	340.37	$4d^{10}5s5d, ^3D_1$	$4d^{10}5s5p, ^3P_0$	7.7×10^7	7.3753	3.7337
UV	346.62	$4d^{10}5s5d, ^3D_2$	$4d^{10}5s5p, ^3P_1$	1.2×10^8	7.3768	3.8009
UV	361.05	$4d^{10}5s5d, ^3D_3$	$4d^{10}5s5p, ^3P_2$	1.3×10^8	7.3790	3.9460
blue	441.30	$4d^{10}5s6s, ^1S_0$	$4d^{10}5s5p, ^3P_1$	2×10^5	6.6096	3.8009
blue	466.24	$4d^{10}5s6d, ^1D_2$	$4d^{10}5s5p, ^1P_1$	5.5×10^6	8.0757	5.4172
blue	467.81	$4d^{10}5s6s, ^3S_1$	$4d^{10}5s5p, ^3P_0$	1.3×10^7	6.3832	3.7337
blue	479.99	$4d^{10}5s6s, ^3S_1$	$4d^{10}5s5p, ^3P_1$	4.1×10^7	6.3832	3.8009
green	508.58	$4d^{10}5s6s, ^3S_1$	$4d^{10}5s5p, ^3P_2$	5.6×10^7	6.3832	3.9460
green	515.47	$4d^{10}5s7s, ^1S_0$	$4d^{10}5s5p, ^1P_1$	3×10^5	7.8218	5.4172
red	632.52	$4d^{10}5s5d, ^3D_2$	$4d^{10}5s5p, ^1P_1$	6×10^6	7.3768	5.4172
red	633.00	$4d^{10}5s5d, ^3D_1$	$4d^{10}5s5p, ^1P_1$	2×10^6	7.3753	5.4172
red	643.85	$4d^{10}5s5d, ^1D_2$	$4d^{10}5s5p, ^1P_1$	5.9×10^7	7.3423	5.4172

Table 8: Cd transitions - NIST data




Fig. 8: Hg Grotrian diagram

colour	λ [nm]	upper (1)	lower (2)	A [s^{-1}]	E_1 [eV]	E_2 [eV]
purple	365.3	$5d^{10}6s6d, {}^3D_{1,3}$	$5d^{10}6s6p, {}^3P_2$	1.3×10^8	8.8500	5.4606
blue	404.66	$5d^{10}6s7p, {}^3S_1$	$5d^{10}6s6p, {}^3P_0$	2.1×10^7	7.7305	4.6674
blue	407.78	$5d^{10}6s7s, {}^1S_0$	$5d^{10}6s6p, {}^3P_1$	4.0×10^6	7.9261	4.8865
blue	410.81	$5d^{10}6s8s, {}^1S_0$	$5d^{10}6s6p, {}^1P_1$	3.0×10^6	9.7209	6.7037
blue	433.92	$5d^{10}6s7d, {}^3D_2$	$5d^{10}6s6p, {}^1P_1$	2.9×10^6	9.5601	6.7037
blue	434.75	$5d^{10}6s7d, {}^1D_2$	$5d^{10}6s6p, {}^1P_1$	8.4×10^6	9.5547	6.7037
blue	435.83	$5d^{10}6s7p, {}^3S_1$	$5d^{10}6s6p, {}^3P_1$	5.6×10^7	7.7305	4.8865
cyan	491.61	$5d^{10}6s8s, {}^1S_0$	$5d^{10}6s6p, {}^1P_1$	5.8×10^6	9.2250	6.7037
green	546.08	$5d^{10}6s7p, {}^3S_1$	$5d^{10}6s6p, {}^3P_2$	4.9×10^7	7.7305	5.4606
yellow	576.96	$5d^{10}6s6d, {}^3D_2$	$5d^{10}6s6p, {}^1P_1$	2.4×10^7	8.8520	6.7037
yellow	578.97	$5d^{10}6s6d, {}^3D_1$	$5d^{10}6s6p, {}^1P_1$	2.4×10^7	8.8445	6.7037
red	623.44	$5d^{10}6s9p, {}^1P_1$	$5d^{10}6s7s, {}^1S_0$	5.3×10^5	9.9143	7.9261
red	690.75	$5d^{10}6s8p, {}^3P_2$	$5d^{10}6s7p, {}^3S_1$	2.8×10^6	9.5249	7.7305

Table 9: Hg transitions -NIST data

Table 10 and 11 show a comparison of the wavelength determined by reading off the wavelength from the interpolated data measured with the prism assuming the He wavelength known (Fig. 5) with the reference data from NIST.

Zn deg. min.	β_λ [°]	λ , Fig. 5 [nm]	λ , NIST [nm]	Cd deg. min.	β_λ [°]	λ , Fig. 5 [nm]	λ , NIST [nm]
288	25	64.52	414	289	14	63.70	429
289	9	63.78	427	289	42	63.23	439
289	11	63.75	428	290	40	62.27	466
290	32	62.40	462	290	42	62.23	467
290	43	62.22	468	291	7	61.82	480
290	52	62.07	474	291	55	61.02	512
291	9	61.78	481	292	3	60.88	518
291	54	61.03	511	293	42	59.23	612
292	7	60.82	520	293	43	59.22	613
292	24	60.53	534	294	0	58.93	636
294	2	58.90	640	294	0	58.93	636
				294	7	58.82	646

Table 10: Measurement data for zinc (Zn) and cadmium (Cd)

Hg	deg.	min.	β_λ [°]	λ , Fig. 5 [nm]	λ , NIST [nm]
287	48	65.13		406	404.7
287	58	64.97		408	407.8
288	9	64.78		410	410.8
289	21	63.58		432	433.9
289	24	63.53		433	434.7
289	28	63.47		434	435.8
291	27	61.48		492	491.6
292	43	60.22		550	546.1
293	14	59.70		580	577.0
293	15	59.68		581	579.0
293	52	59.07		626	623.4
294	40	58.27		684	690.8

Table 11: Measurement data for mercury (Hg)

Evaluation

While for helium eight lines were observable and identifiable with the grating, with the prism it was 14 lines. With the grating for sodium the averaged absolute value of deviation of first order measured wavelength from literature values was 0.16 nm.

Compared to that this value is for the prism for Zn, Cd and Hg with 2 nm more than ten times larger.

For zinc with atomic number 30 no transitions between singulett and tripplett states are observable.

For cadmium with atomic number 48 transitions at 441 nm and 633 nm are just observable, which are transitions between singulett and tripplett states.

For mercury with atomic number 80 the yellow doublett at 578 nm is one of the most prominent lines and other transitions between singulett and tripplett states are well observable, too.

This means that for light elements with low atomic number different selection rules apply for transitions than for heavy elements.

For light elements the selection rule for spin quantum number S is $\Delta S = 0$. LS or Russel-Sounders coupling applies: The two spins s_i couple to form a total spin S and the two orbital moments l_i couple to form a total orbital angular moment L , $L = \sum_i l_i$, $S = \sum_i s_i$, $J = L + S$.

For heavy atoms for the outermost electrons the interaction between the individual spin s_i and orbit l_i is stronger than between different spins so first each $s_i + l_i = j_i$ couple and $J = \sum_i j_i = \sum_i (l_i + s_i)$. This is referred to as jj -coupling. Here the selection rule $\Delta S = 0$ is lifted.

Refer to:

https://en.wikipedia.org/wiki/Angular_momentum_coupling

https://en.wikipedia.org/wiki/Spin%E2%80%93orbit_interaction

Appendix

Atomic level data (source: NIST)

He levels			Na levels			Zn levels			Zn levels		
config.	term	E [eV]	config.	term	E [eV]	config.	term	E [eV]	config.	term	E [eV]
<i>singulett</i>			<i>doublett</i>			<i>singulett</i>			<i>triplett</i>		
1s ²	¹ S ₀	0.0000	2p ⁶ 3s	² S _{1/2}	0.0000	3d ¹⁰ 4s ²	¹ S ₀	0.0000	3d ¹⁰ 4s4p	³ P ₀	4.0061
1s2s	¹ S ₀	20.6158	2p ⁶ 3p	² P _{1/2}	2.1023	3d ¹⁰ 4s4p	¹ P ₁	5.7957	3d ¹⁰ 4s4p	³ P ₁	4.0297
1s2p	¹ P ₁	21.2180	2p ⁶ 3p	² P _{3/2}	2.1044	3d ¹⁰ 4s5s	¹ S ₀	6.9170	3d ¹⁰ 4s4p	³ P ₂	4.0779
1s3s	¹ S ₀	22.9203	2p ⁶ 4s	² S _{1/2}	3.1914	3d ¹⁰ 4s5s	¹ D ₂	7.7439	3d ¹⁰ 4s5s	³ S ₁	6.6545
1s3d	¹ D ₂	23.0741	2p ⁶ 3d	² D _{5/2,3/2}	3.6170	3d ¹⁰ 4s4d	¹ P ₁	7.7999	3d ¹⁰ 4s5p	³ P ₀	7.5938
1s3p	¹ P ₁	23.0870	2p ⁶ 4p	² P _{1/2}	3.7526	3d ¹⁰ 4s5p	¹ S ₀	8.1876	3d ¹⁰ 4s5p	³ P ₁	7.5971
1s4s	¹ S ₀	23.6736	2p ⁶ 4p	² P _{3/2}	3.7533	3d ¹⁰ 4s5d	¹ D ₂	8.4729	3d ¹⁰ 4s5p	³ P ₂	7.6041
1s4p	¹ P ₁	23.7421	2p ⁶ 5s	² S _{1/2}	4.1164	3d ¹⁰ 4s6p	¹ P ₁	8.5062	3d ¹⁰ 4s4d	³ D ₁	7.7823
1s4d	¹ D ₂	23.7363	2p ⁶ 4d	² D _{5/2,3/2}	4.2835	3d ¹⁰ 4s4f	¹ F ₃	8.5344	3d ¹⁰ 4s4d	³ D ₂	7.7827
1s4f	¹ F ₃	23.7370	2p ⁶ 4f	² F _{7/2,5/2}	4.2882	3d ¹⁰ 4s7s	¹ S ₀	8.6794	3d ¹⁰ 4s4d	³ D ₃	7.7834
1s5s	¹ S ₀	24.0112	2p ⁶ 5p	² P _{1/2}	4.3445	3d ¹⁰ 4s6d	¹ D ₂	8.8091	3d ¹⁰ 4s6s	³ S ₁	8.1126
1s5d	¹ D ₂	24.0428	2p ⁶ 5p	² P _{3/2}	4.3448	3d ¹⁰ 4s7p	¹ P ₁	8.8300	3d ¹⁰ 4s6p	³ P ₀	8.4397
1s5f	¹ F ₃	24.0432	2p ⁶ 6s	² S _{1/2}	4.5096	3d ¹⁰ 4s5f	¹ F ₃	8.8446	3d ¹⁰ 4s6p	³ P ₁	8.4409
1s5p	¹ P ₁	24.0458	2p ⁶ 5d	² D _{5/2,3/2}	4.5920	3d ¹⁰ 4s8s	¹ S ₀	8.9214	3d ¹⁰ 4s6p	³ P ₂	8.4435
1s6s	¹ S ₀	24.1690	2p ⁶ 5f	² F _{7/2,5/2}	4.5946	3d ¹⁰ 4s7d	¹ D ₂	8.9910	3d ¹⁰ 4s5d	³ D _{1,2,3}	8.5030
1s6d	¹ D ₂	24.2092	2p ⁶ 6p	² P _{1/2}	4.6242	3d ¹⁰ 4s8p	¹ P ₁	9.0045	3d ¹⁰ 4s4f	³ F _{2,3,4}	8.5343
<i>triplett</i>			2p ⁶ 6p	² P _{3/2}	4.6243	3d ¹⁰ 4s9s	¹ S ₀	9.0584	3d ¹⁰ 4s7s	³ S ₁	8.6474
1s2s	³ S ₁	19.8196	2p ⁶ 7s	² S _{1/2}	4.7129	3d ¹⁰ 4s8d	¹ D ₂	9.0999	3d ¹⁰ 4s7p	³ P _{0,1,2}	8.8013
1s2p	³ P _{0,1,2}	20.9641	2p ⁶ 6d	² D _{5/2,3/2}	4.7594	3d ¹⁰ 4s9p	¹ P ₁	9.1090	3d ¹⁰ 4s6d	³ D _{1,2,3}	8.8293
1s3s	³ S ₁	22.7185	2p ⁶ 6f	² F _{7/2,5/2}	4.7610				3d ¹⁰ 4s5f	³ F _{2,3,4}	8.8445
1s3p	³ P _{0,1,2}	23.0071	2p ⁶ 7p	² P _{1/2}	4.7784				3d ¹⁰ 4s5g	³ G _{3,4,5}	8.8492
1s3d	³ D _{1,2,3}	23.0737	2p ⁶ 7p	² P _{3/2}	4.7785				3d ¹⁰ 4s8p	³ P _{0,1,2}	9.8990
1s4s	³ S ₁	23.5940	2p ⁶ 8s	² S _{1/2}	4.8315				3d ¹⁰ 4s6f	³ F _{2,3,4}	9.0125
1s4p	³ P _{0,1,2}	23.7079							3d ¹⁰ 4s9p	³ P _{0,1,2}	9.0998
1s4d	³ D _{1,2,3}	23.7361							3d ¹⁰ 4s6g	³ G _{3,4,5}	9.0149
1s4f	³ F _{2,3,4}	23.7370							3d ¹⁰ 4s7g	³ G _{3,4,5}	9.1150
1s5s	³ S ₁	23.9720									
1s5p	³ P _{0,1,2}	24.0282									
1s5f	³ F _{2,3,4}	24.0432									
1s5d	³ D _{1,2,3}	24.0427									

Cd levels		Cd levels	
config.	term E [eV]	config.	term E [eV]
<i>singulett</i>		<i>triplett</i>	
4d ¹⁰ 5s ²	¹ S ₀ 0.0000	4d ¹⁰ 5s5p	³ P ₀ 3.7337
4d ¹⁰ 5s5p	¹ P ₁ 5.4172	4d ¹⁰ 5s5p	³ P ₁ 3.8009
4d ¹⁰ 5s6s	¹ S ₀ 6.6096	4d ¹⁰ 5s5p	³ P ₂ 3.9460
4d ¹⁰ 5s5d	¹ D ₂ 7.3423	4d ¹⁰ 5s6s	³ S ₁ 6.3832
4d ¹⁰ 5s6p	¹ P ₁ 7.4276	4d ¹⁰ 5s6p	³ P ₀ 7.2396
4d ¹⁰ 5s7s	¹ S ₀ 7.8218	4d ¹⁰ 5s6p	³ P ₁ 7.2483
4d ¹⁰ 5s6d	¹ D ₂ 8.0757	4d ¹⁰ 5s6p	³ P ₂ 7.2699
4d ¹⁰ 5s7p	¹ P ₁ 8.1211	4d ¹⁰ 5s5d	³ D ₁ 7.3753
4d ¹⁰ 5s8s	¹ S ₀ 8.2952	4d ¹⁰ 5s5d	³ D ₂ 7.3768
4d ¹⁰ 5s7d	¹ D ₂ 8.4109	4d ¹⁰ 5s5d	³ D ₃ 7.3790
4d ¹⁰ 5s8p	¹ P ₁ 8.4383	4d ¹⁰ 5s7s	³ S ₁ 7.7569
4d ¹⁰ 5s9s	¹ S ₀ 8.5300	4d ¹⁰ 5s7p	³ P ₀ 8.0585
4d ¹⁰ 5s8d	¹ D ₂ 8.5917	4d ¹⁰ 5s7p	³ P ₁ 8.0621
4d ¹⁰ 5s9p	¹ P ₁ 8.6094	4d ¹⁰ 5s7p	³ P ₂ 8.0706
4d ¹⁰ 5s10s	¹ S ₀ 8.6634	4d ¹⁰ 5s6d	³ D ₁ 8.1028
4d ¹⁰ 5s9d	¹ D ₂ 8.6999	4d ¹⁰ 5s6d	³ D ₂ 8.1035
4d ¹⁰ 5s10p	¹ P ₁ 8.7120	4d ¹⁰ 5s6d	³ D ₃ 8.1045
4d ¹⁰ 5s11s	¹ S ₀ 8.7411	4d ¹⁰ 5s4f	³ F _{2,3,4} 8.1320
4d ¹⁰ 5s10d	¹ D ₂ 8.7701	4d ¹⁰ 5s8s	³ S ₁ 8.2675
4d ¹⁰ 5s10f	¹ F ₃ 8.8566	4d ¹⁰ 5s8p	³ P ₀ 8.4098
		4d ¹⁰ 5s8p	³ P ₁ 8.4113
		4d ¹⁰ 5s8p	³ P ₂ 8.4155
		4d ¹⁰ 5s7d	³ D _{1,2,3} 8.4302
		4d ¹⁰ 5s5f	³ F _{2,3,4} 8.4425
		4d ¹⁰ 5s9s	³ S ₁ 8.5155
		4d ¹⁰ 5s9p	³ P _{0,1,2} 8.5959
		4d ¹⁰ 5s8d	³ D _{1,2,3} 8.6049
		4d ¹⁰ 5s6f	³ F _{2,3,4} 8.6115
		4d ¹⁰ 5s10s	³ S ₁ 8.6549
		4d ¹⁰ 5s10p	³ P _{0,1,2} 8.7039
		4d ¹⁰ 5s7f	³ F _{2,3,4} 8.7133
		4d ¹⁰ 5s11s	³ S ₁ 8.7411

Hg levels	
config.	term E [eV]
<i>singulett</i>	
5d ¹⁰ 6s ²	¹ S ₀ 0.0000
5d ¹⁰ 6s6p	¹ P ₁ 6.7037
5d ¹⁰ 6s7s	¹ S ₀ 7.9261
5d ¹⁰ 6s7p	¹ P ₁ 8.8394
5d ¹⁰ 6s6d	¹ D ₂ 8.8442
5d ¹⁰ 6s8s	¹ S ₀ 9.2250
5d ¹⁰ 6s8p	¹ P ₁ 9.5298
5d ¹⁰ 6s7d	¹ D ₂ 9.5547
5d ¹⁰ 6s5f	¹ F ₃ 9.5767
5d ¹⁰ 6s9s	¹ S ₀ 9.7209
<i>triplett</i>	
5d ¹⁰ 6s6p	³ P ₀ 4.6674
5d ¹⁰ 6s6p	³ P ₁ 4.8865
5d ¹⁰ 6s6p	³ P ₂ 5.4606
5d ¹⁰ 6s7p	³ S ₁ 7.7305
5d ¹⁰ 6s7p	³ P ₀ 8.6190
5d ¹⁰ 6s7p	³ P ₁ 8.6370
5d ¹⁰ 6s7p	³ P ₂ 8.8286
5d ¹⁰ 6s6d	³ D ₁ 8.8445
5d ¹⁰ 6s6d	³ D ₂ 8.8520
5d ¹⁰ 6s6d	³ D ₃ 8.8563
5d ¹⁰ 6s8s	³ S ₁ 9.1700
5d ¹⁰ 6s8p	³ P ₀ 9.4782
5d ¹⁰ 6s8p	³ P ₁ 9.4807
5d ¹⁰ 6s8p	³ P ₂ 9.5249
5d ¹⁰ 6s7d	³ D ₁ 9.5573
5d ¹⁰ 6s7d	³ D ₂ 9.5601
5d ¹⁰ 6s7d	³ D ₃ 9.5628
5d ¹⁰ 6s5f	³ F ₂ 9.5762
5d ¹⁰ 6s5f	³ F ₃ 9.5764
5d ¹⁰ 6s5f	³ F ₄ 9.5823