

# Qualitative X-ray fluorescence analysis of powder samples



Physics

Modern Physics

Production &amp; use of X-rays



Difficulty level

hard



Group size

2



Preparation time

45+ minutes



Execution time

45+ minutes

This content can also be found online at:

<http://localhost:1337/c/5f7dd0072c24d60003b46120>

PHYWE

# General information



## Application

PHYWE



Setup

Most applications of X rays are based on their ability to pass through matter. Since this ability is dependent on the density of the matter, imaging of the interior of objects and even people becomes possible. This has wide usage in fields such as medicine or security.

## Other information (1/2)

PHYWE

**Prior****knowledge****Main****principle**

The prior knowledge for this experiment is found in the Theory section.

Various powder samples are subjected to polychromatic X-rays. The energy of the resulting fluorescence radiation is analysed with the aid of a semiconductor detector and a multichannel analyser. The energy of the corresponding characteristic X-ray fluorescence lines is determined. The elements of the samples are identified by comparing the line energies with the corresponding table values.

## Other information (2/2)

PHYWE

**Learning  
objective****Tasks**

The goal of this experiment is to get to investigate the spectra of fluorescence radiation.

1. Calibrate the semiconductor energy detector with the aid of the characteristic radiation of the molybdenum X-ray tube.
2. Record the fluorescence spectra that are produced by the samples.
3. Determine the energy values of the corresponding fluorescence lines and compare them with the literature values in order to identify the powder components.

## Safety Instructions

PHYWE



When handling chemicals, you should wear suitable protective gloves, safety goggles, and suitable clothing. Please refer to the appendix for detailed safety instructions.

## Theory

PHYWE

When X-rays interact with matter, they lose energy. In the range of energy that is available during this experiment, the photoelectric effect plays the most important role. This means that on one of the lower shells of the atom, an electron is ejected due to the absorbed photon energy. The now free space is taken by an electron from one of the higher shells. The energy that is produced during this process can be used for the ejection of another electron from one of the higher shells (Auger effect) or for the generation of a photon (fluorescence radiation). Since the energy of the energy levels that are involved in this process is atom-specific, the type of the emitting atom can be determined based on the energy of the fluorescence radiation.


In order to determine the type of atom, the experimental energy values are then compared to the corresponding table values (e.g. "Handbook of Chemistry and Physics"- CRC-Press, Inc. USA). During the assignment of the fluorescence lines, it must be taken into consideration that the relaxations that follow the primary ionisation process can only take place if they fulfil the quantum-mechanical selection rules  $\Delta j = 0, \pm 1$  and  $\Delta l = \pm 1$  ( $j$  = total angular momentum,  $l$  = orbital angular momentum).

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">XR 4.0 expert unit, 35 kV</a>	09057-99	1
2	<a href="#">XR 4.0 X-ray goniometer</a>	09057-10	1
3	<a href="#">XR4 X-ray Plug-in Cu tube</a>	09057-51	1
4	<a href="#">XR 4.0 X-ray material upgrade set</a>	09165-88	1
5	<a href="#">XR 4.0 X-ray Chemical set for edge absorption</a>	09056-07	1

PHYWE

# Setup and Procedure



## Setup (1/2)

PHYWE

- Screw the adapter ring onto the inlet tube of the energy detector and connect the signal and supply cables to the corresponding ports of the detector with the aid of the right-angle plugs.
- Connect the signal and supply cables to the corresponding ports in the experiment chamber of the X-ray unit. In Figure 1, the port for the signal cable is marked in red and the port for the supply cable is marked in green. Connect the external X RED ports of the x-ray unit (see Fig. 2) to the multi-channel analyser (MCA). Connect the signal cable to the "Input" port and the supply cable to the "X-Ray Energy Det." port of the MCA.




Fig. 1:  
Connectors in  
the  
experiment  
chamber



Fig. 2: Connection of the  
multi-channel analyser

## Setup (2/2)

PHYWE

- Secure the energy detector in the holder of the swivel arm of the goniometer. Lay the two cables with sufficient length so that the goniometer can be swivelled freely over the entire range.
- Connect the multi-channel analyser and computer with the aid of the USB cable.




Fig. 3: Goniometer set-up

## Procedure (1/4)

PHYWE

- Bring the goniometer block and the detector to their respective end positions on the right.
- Insert the tube with the 1mm-aperture into the exit tube of the x-ray tube.
- With the X-ray unit switched on and the door locked, bring the detector to the 0° position. Then, shift the detector by some tenths degree out of the zero position in order to reduce the total rate.
- Operating data of the x-ray tube: Select an anode voltage  $U_A = 25$  kV and an anode current  $I_A = 0.02$  mA and confirm these values by pressing the "Enter" button.
- Switch on the x-radiation




Fig. 4: calibration of the multi-channel analyser

## Procedure (2/4)

PHYWE

- In the MEASURE program, select "Multi channel analyser" under "Gauge". Then, select "Settings and calibration". After the "Calibrate" button has been clicked, a spectrum can be measured. The counting rate should be  $< 300$  c/s. Energy calibration settings: - 2-point calibration, - Unit = keV, Gain = 2 – Set the offset so that low-energy noise signals will be suppressed (usually a few per cent are sufficient), See Fig 4.
- Measuring time: 5 minutes. Use the timer of the X-ray unit.
- Make the two coloured calibration lines congruent with the line centres of the two characteristic X-ray lines. The corresponding energy values (see e.g. P2544705)  $E(L_3M_5/L_3M_4) = 8,41\text{keV}$  and  $E(L_2N_4) = 9,69\text{keV}$  are entered into the corresponding fields, depending on the colour. (Note: Since a separation of the lines  $L_3M_4$  and  $L_3M_5$  Lines is not possible, the mean value of both lines is entered as the energy of the line).
- Name and save the calibration.

## Procedure (3/4)

PHYWE

### Spectrum recording






- Insert the tube with the 2-mm-aperture.
- Bring the goniometer block and the detector to their respective end positions on the left. Bring the detector to the  $90^\circ$  position in the 1:2 coupling mode.
- Insert the metal sample with the universal crystal holder (sample at  $45^\circ$ ).
- Operating data of the x-ray tube: Adjust an anode voltage  $U_A = 35\text{ kV}$  and an anode current so that the counting rate is  $\leq 300$  c/s.
- Measuring time: 3 minutes (use the timer of the X-ray unit).



## Procedure (4/4)


PHYWE

### Evaluation of the measurement curves

- In order to determine the line energy, switch from the bar display to the curve display. To do so, click “Display options” and then “Interpolation and straight lines”.
- Extend the relevant line section with the aid of the zoom function 
- Then, select the curve section with  Open the window “Function fitting”  Then, select “Scaled normal distribution” and confirm.
- Find the line centroid of the normal distribution with “Peak analysis”  or determine it with the function “Survey” 

PHYWE

## Evaluation



## Evaluation (1/4)

PHYWE

### Evaluation of the $K_\alpha$ and $K_\beta$ fluorescence lines of the elements with an atomic number $30 < Z < 38$

Figure 5 shows an overview of all of the fluorescence spectra of the elements with an atomic number  $30 < Z < 38$ . Within this group of elements, only the high-intensity characteristic  $K_\alpha$  and  $K_\beta$  radiation can be clearly identified. The energy of the characteristic fluorescence radiation of the sample elements oxygen and sulphur of the oxide or sulphate compounds or of potassium are below the sensitivity limit of the energy detector and partly only have a very low intensity, which is why they cannot be identified in this experiment.




Fig. 5: Total representation of the  $K_\alpha$  and  $K_\beta$  fluorescence lines of the elements with an atomic number  $30 < Z < 38$ .

## Evaluation (2/4)

PHYWE

Figure 5 clearly shows that at the higher atomic numbers of the elements, the energy of the characteristic lines as well as the energetic spacing between the  $K_{\alpha_1}$  and  $K_{\beta_2}$  lines increases.

Table 1 shows the evaluation of the spectra of Fig. 5.

Atomic number Z	Element	$E_{\text{exp.}} / \text{keV}$	$E_{\text{lit.}} / \text{keV}$	Transition
30	Zn	8.60 / 9.52	8.64 / 9.57	$K_\alpha / K_\beta$
32	Ge	9.86 / 10.95	9.89 / 10.98	$K_\alpha / K_\beta$
34	Se	11.19 / 12.46	11.22 / 12.49	$K_\alpha / K_\beta$
35	Br	11.88 / 13.24	11.92 / 13.29	$K_\alpha / K_\beta$
37	Rb	13.34 / 14.93	13.39 / 14.96	$K_\alpha / K_\beta$
38	Sr	14.12 / 15.82	14.16 / 15.83	$K_\alpha / K_\beta$

Table 1 shows the evaluation of the spectra of Fig. 5.

## Evaluation (3/4)

PHYWE

### Evaluation of the L fluorescence lines of elements with a higher atomic number

Figure 6 shows the spectra of the powder samples with a higher atomic number. The ionisation energy of the K level of these elements is above the available excitation energy of the primary radiation of 35 keV max. This is why only the L fluorescence lines can be excited in this case.

Table 2 shows the evaluation of the spectra of Fig.6.




Fig. 6: L fluorescence lines of tungsten, lead and bismuth

## Evaluation (4/4)

PHYWE

The  $L_{\beta_1}$  lines, which are slightly widened, also include parts of the  $L_{\beta_{2,3,4}}$  lines that are of a lower intensity. The  $L_{\gamma_1}$  line clearly becomes asymmetric towards the high-energy edge. This is due to the lower-intensity  $L_{\gamma_{2,3}}$  lines with a slightly higher energy. The spectra of lead and bismuth include the additional low-intensity line (line 3). The fluorescence line at the end of the spectra ( $E = 17.36$  keV) is caused by the primary molybdenum  $K_{\alpha}$  radiation that is scattered at the samples.

Element	Line	$E_{exp.}$ / keV	$E_{lit.}$ / keV	Transition
W, Z = 74	1	7.40	7.39	$L_{\beta_1}M_1 - L_I$
	2	8.36	8.35 / 8.42	$L_{\beta_2}M_{4,5} - L_{\alpha 1,2}$
	3	9.67	9.69	$L_{\beta_3}M_4 - L_{\beta_1}$
	4	11.26	11.28	$L_{\beta_4}N_4 - L_{\gamma_1}$
Pb, Z = 82	1	9.15	9.18	$L_{\beta_1}M_1 - L_I$
	2	10.49	10.45 / 10.55	$L_{\beta_2}M_{4,5} - L_{\alpha 1,2}$
	3	11.30	11.35	$L_{\beta_3}M_4 - L_{\beta_1}$
	4	12.55	12.61	$L_{\beta_4}N_4 - L_{\gamma_1}$
Bi, Z = 83	5	14.73	14.76	$L_{\beta_5}N_4 - L_{\gamma_1}$
	1	9.36	9.42	$L_{\beta_1}M_1 - L_I$
	2	10.77	10.73 / 10.83	$L_{\beta_2}M_{4,5} - L_{\alpha 1,2}$
	3	11.67	11.71	$L_{\beta_3}M_4 - L_{\beta_1}$
	4	12.95	13.02	$L_{\beta_4}N_4 - L_{\gamma_1}$
	5	15.20	15.25	$L_{\beta_5}N_4 - L_{\gamma_1}$

Table 2: Elements of the powder samples with L fluorescence radiation


PHYWE



# Appendix

## Security Information (1/7)

PHYWE

Hazard symbol, signal word	Hazard statements	Precautionary statements
Rubidium chloride (RbCl)		
-	--	--
Germanium(IV) oxide (GeO <sub>2</sub> )		
	H302: Harmful if swallowed H332: Harmful if inhaled	--

## Security Information (2/7)

PHYWE

Hazard symbol, signal word Hazard statements

Precautionary statements

Silver (I) chloride (AgCl)



H400: Very toxic to aquatic life

P273: Avoid release to the environment.

Zinc

-

--

--

## Security Information (3/7)

PHYWE

Hazard symbol, signal word

Hazard statements

Precautionary statements

Selenium



H301: Toxic if swallowed

H331 Toxic if inhaled

H373: Causes damage to organs through prolonged or repeated exposure

H413: May cause long lasting harmful effects to aquatic life

--.

## Security Information (4/7)

PHYWE

Hazard symbol, signal word	Hazard statements	Precautionary statements
<b>Potassium bromide (KBr)</b> 	H315: Causes skin irritation H319: Causes serious eye irritation H335: May cause respiratory irritation	P261: Avoid breathing dust/fume/gas/mist/vapours/spray. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
<b>Strontium sulphate (SrSO<sub>4</sub>)</b>		
-	--	--

## Security Information (5/7)

PHYWE

Hazard symbol, signal word	Hazard statements	Precautionary statements
<b>Lead (IV) oxide (PbO<sub>2</sub>)</b> 	H272: May intensify fire; oxidiser H302: Harmful if swallowed H332: Harmful if inhaled H360: May damage fertility or the unborn child H373: Causes damage to organs through prolonged or repeated exposure H410: Very toxic to aquatic life with long lasting effects	P201: Obtain special instructions before use. P220: Keep/Store away from clothing/.../combustible materials. P273: Avoid release to the environment. P308 + P313: IF exposed or concerned: Get medical advice/attention.

## Security Information (6/7)

PHYWE

Hazard symbol, signal word	Hazard statements	Precautionary statements
<b>Tungsten(IV) oxide (WO<sub>2</sub>)</b>		
	H335: May cause respiratory irritation	--



H335: May cause respiratory irritation

--

## Security Information (7/7)

PHYWE

Hazard symbol, signal word	Hazard statements	Precautionary statements
<b>Bismuth(III) oxide (Bi<sub>2</sub>O<sub>3</sub>)</b>		
	H315: Causes skin irritation H319: Causes serious eye irritation H335: May cause respiratory irritation	P261: Avoid breathing dust/fume/gas/mist/vapours/spray. P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Bismuth(III) oxide (Bi<sub>2</sub>O<sub>3</sub>)

H315: Causes skin irritation  
 H319: Causes serious eye irritation  
 H335: May cause respiratory irritation

P261: Avoid breathing dust/fume/gas/mist/vapours/spray.  
 P305 + P351 + P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.