

Dissociation constants



Students will be able to experimentally determine the extinction (absorbance) of an aqueous solution of thymol blue (thymolsulphonephthalein) in dilute HCl, NaOH and a buffer of known pH value as a function of wavelength between 400 and 700 nm at constant concentration and constant temperature.

Chemistry

General Chemistry

Substances mixtures & separation

Chemistry

General Chemistry

Chemical reactions

Chemical balance



Difficulty level

medium



Group size

2



Preparation time

10 minutes



Execution time

10 minutes


This content can also be found online at:



<http://localhost:1337/c/60c10d3672987c000430de5b>


PHYWE

General information



Application

PHYWE



Experimental setup

In this experiment, students learn how to use and work with the dissociation constant by using Photometric measurements.

Photometric measurements in the visible spectral range can therefore be used to advantage to determine the position of the K_a and pK_a values of the indicator which characterize dissociation equilibrium.

The students will experimentally determine the extinction (absorbance) of an aqueous solution of thymol blue (thymolsulphonephthalein) in dilute HCl and NaOH and calculate the dissociation constant.

Other information (1/2)

PHYWE

Prior knowledge



Students should be familiar with the true and potential electrolytes, Strong and weak acids, Law of mass action, Dissociation constants and pK_a values, Henderson-Hasselbach equation, UV-visible spectrometry, Lambert-Beer's Law and photometry in general.

Scientific principle



Photometric measurements in the visible spectral range can therefore be used to advantage to determine the position of the K_a and pK_a values of the indicator which characterize dissociation equilibrium.

Other information (2/2)

PHYWE

Learning objective



Students will be able to experimentally determine the extinction (absorbance) of an aqueous solution of thymol blue (thymolsulphonephthalein) in dilute HCl, NaOH and a buffer of known pH value as a function of wavelength between 400 and 700 nm at constant concentration and constant temperature.

Tasks



The students will experimentally determine the extinction (absorbance) of an aqueous solution of thymol blue (thymolsulphonephthalein) in dilute HCl, NaOH and a buffer of known pH value as a function of wavelength between 400 and 700 nm at constant concentration and constant temperature. They are to calculate the dissociation constant (indicator constant) K_a from the measurement results.

Safety instructions

PHYWE



- Use protective glasses !
- For this experiment the general instructions for safe experimentation in science lessons apply.
- For H- and P-phrases please consult the safety data sheet of the respective chemical.

Theory

PHYWE

The coloured indicator thymol blue is a weak acid that is partially dissociated in aqueous solution, whereby non-ionized and ionized forms show absorption maximums at different wavelengths in the visible range.


Photometric measurements in the visible spectral range can therefore be used to advantage to determine the position of the K_a and pK_a values of the indicator which characterize dissociation equilibrium.

The HA form that exists alone in acidic solutions absorbs in the blue spectral range ($\lambda_{max} \approx 430nm$) and appears as the complementary colour, yellow. On the other hand, the blue salt form A^- that exists alone in basic solutions absorbs in the yellow spectral range ($\lambda_{max} \approx 595nm$). In buffered solutions a mixed colour is observed because of the existence of HA and A^- , both of which absorb in the visible spectral range.

Equipment

Position	Material	Item No.	Quantity
1	UV-VIS spectrophotometer T70 with monitor 190 - 1100 nm	35658-99	1
2	Cell for spectrophotometer, optical glass	35664-00	2
3	Weighing dishes, square shape, 84 x 84 x 24 mm, 500 pcs.	45019-50	1
4	Volumetric flask, Borosilicate, 50 ml, IGJ12/21	36547-00	3
5	Volumetric flask 1000ml, IGJ24/29	36552-00	3
6	Volumetric pipette, 1 ml	36575-00	1
7	Volumetric pipette, 5 ml	36577-00	1
8	Volumetric pipette, 10 ml	36578-00	1
9	Pipettor	36592-00	1
10	Pipette dish	36589-00	1
11	Graduated cylinder, Borosilicate, 250 ml	36630-00	1
12	Beaker, Borosilicate, tall form, 150 ml	46032-00	3
13	Pasteur pipettes, 250 pcs	36590-00	1
14	Rubber caps, 10 pcs	39275-03	1
15	Microspoon, steel	33393-00	1
16	Funnel, glass, top dia. 50 mm	34457-00	1
17	Wash bottle, plastic, 500 ml	33931-00	1
18	Thermometer -10...+50 °C	38034-00	1
19	Buffer solution, pH 9 1000 ml	30289-70	1
20	Thymol blue indicator 5 g	31896-02	1
21	Hydrochloric acid, 0.1M 1000 ml	48452-70	1
22	Caustic soda sol., 0.1M 1000 ml	48328-70	1
23	Ethyl alcohol, absolute 500 ml	30008-50	1
24	Water, distilled 5 l	31246-81	1

PHYWE



Setup and procedure

Setup and Procedure (1/2)

PHYWE

- The spectrophotometer that is required for this experiment is shown in Fig. right.
- Pipette 10 ml of 0.1 molar sodium hydroxide and 1 ml of 0.1 molar hydrochloric acid into a 1000 ml volumetric flask and fill up to the mark with distilled water.
- Completely dissolve 0.145 g ($3 \cdot 10^{-4} \text{ mol}$) of thymol blue (thymolsulphonephthalein, $\text{C}_{27}\text{H}_{30}\text{O}_5\text{S} \cdot \text{H}_2\text{O}$) in 200 ml of ethanol in a 1000 ml volumetric flask and dilute up to the mark with distilled water.



Spectrophotometer


Setup and Procedure (2/2)

PHYWE

- Pipette 5 ml of this $3 \cdot 10^{-4}$ molar stock solution into each of three 50 ml volumetric flasks: Fill the first flask up to the mark with sodium hydroxide ($c = 1 \cdot 10^{-3} \text{ mol} \cdot \text{l}^{-1}$), the second one with HCl ($c = 1 \cdot 10^{-4} \text{ mol} \cdot \text{l}^{-1}$), and the third flask with buffer solution of pH 9.
- After correcting the zero line of the photometer using a waterfilled cell, record the absorption spectra of the three $3 \cdot 10^{-5}$ molar thymol blue solutions in the visible spectral range between 700 and 400 nm at a slow recording speed.
- Read off the extinction values from the spectra displayed by the monitor in 5 nm steps and plot a graph of them as a function of wavelength.

PHYWE

Evaluation



Evaluation (1/13)


PHYWE

Results

The absorption spectra recorded at $T = 299\text{ K}$ are shown in Fig. right. Absorption spectra of thymol blue ($c_0 = 2 \cdot 10^{-5} \text{ mol} \cdot \text{l}^{-1}$) in $1 \cdot 10^{-4}$ molar HCl (\square) $1 \cdot 10^{-3}$ molar NaOH (\circ) and a buffer solution of $\text{pH} = 9.00$ (\triangle) at $T = 299\text{ K}$.

Note

The graphical evaluation of the measured values can be very easily carried out by means of 'Measure' software. A downloadfile of this software is available as freeware for use in evaluating and graphically representing measured values under URL "www.phywe.com". Fig. right was created by this software.



Absorption spectra

Evaluation (2/13)

PHYWE

Evaluation (1/8)

The coloured indicator thymol blue, which is used in analytical practice, is present in aqueous solutions in the partly dissociated weak acid form:



where $A^- = [C_{27}H_{29}O_5S]^-$

The position of the dissociation equilibrium is quantitatively characterized by the acid or indicator constant K_a or the pK_a value derived from it:

$$K_a = \frac{a_{A^-} \cdot a_{H^+}}{a_{HA}} \approx \frac{c_{A^-} \cdot c_{H^+}}{c_{HA}} \quad (1)$$

where a_i Activity of ion i.

Evaluation (3/13)

PHYWE

Evaluation (2/8)

At ideal dilution, the activity a_i is identical to the concentration c_i .

$$pK_a = -\lg \frac{K_a}{\text{mol} \cdot \text{l}^{-1}} \quad (2)$$

From (1), considering formulation (2) and the analogue definition of the pH value and taking logarithms, we obtain the Henderson-Hasselbach equation (3). This describes, at a given acid strength, the connection between the pH value and the composition (c_{HA} / c_{A^-}) of the buffer system, and so the share of the two forms in the total concentration c_0 of the weak acid.

$$pK_a = pH + \lg \frac{c_{HA}}{c_{A^-}} \quad (3)$$

$$c_0 = c_{HA} + c_{A^-} \quad (4)$$

Evaluation (4/13)

PHYWE

Evaluation (3/8)

We so have, in the acid form ($1 \cdot 10^{-4}$ molar HCl, $pH = 4 < pK_a$), practically only the non-ionized acidic form HA, so that (4) becomes (4.1).

$$c_0 = c_{HA} \quad (4.1)$$

In contrast to this, in basic milieu ($1 \cdot 10^{-3}$ molar NaOH, $pH = 11 > pK_a$), equilibrium is shifted almost completely in the direction of the ionized salt form A^- , and we have

$$c_0 = c_{A^-} \quad (4.2)$$

In buffered solutions of $pH \approx pK_a$ the ionized and non-ionized forms are present in practically the same concentration. These equilibrium concentrations, and so the constants K_a and pK_a , for thymol blue, can be advantageously measured via photometric measurements as, because of their different atomic structures, the acid and salt forms give different absorption spectra which cross each other at an isosbestic point.

Evaluation (5/13)

PHYWE

Evaluation (4/8)

The HA form that exists alone in acidic solutions absorbs in the blue spectral range ($\lambda_{max} \approx 430nm$) and appears as the complementary colour, yellow. On the other hand, the blue salt form A^- that exists alone in basic solutions absorbs in the yellow spectral range ($\lambda_{max} \approx 595nm$). In buffered solutions a mixed colour is observed because of the existence of HA and A^- , both of which absorb in the visible spectral range. In analytical practice, the intensity of absorption is usually quantified by the extinction that is defined in equation (5):

$$E_{\lambda} = \lg \frac{I_0}{I} \quad (5)$$

where I_0, I Intensity of the radiation used before and after passage through an absorbing medium.

Evaluation (6/13)

PHYWE

Evaluation (5/8)

Its dependence on the concentration c_i of an absorbing substance i and the layer thickness d at constant wavelength is given by the Lambert-Beer law:

$$E_{\lambda} = \epsilon_i c_i d \quad (6)$$

where

E_{λ} Extinction at wavelength λ

ϵ_i Decadic molar extinction coefficient of the substance i at wavelength λ

c_i Concentration of the substance i

d Layer thickness in the cell

which, with the simultaneous presence of two absorbing substances (here HA and A^-), takes on the form:

$$E_{\lambda} = \epsilon_{HA} c_{HA} d + \epsilon_{A^-} c_{A^-} d \quad (6.1)$$

In acidic solution ($c_{A^-} = 0, c_{HA} = c_0$) or basic solution ($c_{HA} = 0, c_{A^-} = c_0$), equation (6.1) can be simplified to (6.1.1) and (6.1.2).

$$E_{\lambda} = \epsilon_{HA} \cdot c_0 \cdot d = E_{\lambda, HA} \quad (6.1.1)$$

$$E_{\lambda} = \epsilon_{A^-} \cdot c_0 \cdot d = E_{\lambda, A^-} \quad (6.1.2)$$

Evaluation (7/13)

PHYWE

Evaluation (6/8)

Insertion of $c_{HA} = c_0 - c_{A^-}$ (4) in (6.1) gives:

$$c_{A^-} = \frac{E_\lambda - E_{HA} \cdot c_0 \cdot d}{d(A_{A^-} - E_{HA})} \quad (7)$$

Equation (7) is one possibility for calculating the concentration of the ionized form A^- . In this equation E_λ , represents the extinction of the buffered solution that should be measured at a wavelength ($\epsilon_{HA} \neq \epsilon_{A^-}$) that is sufficiently well away from the isosbestic point. The extinction coefficients can be calculated from equations (6.1.1) and (6.1.2) from the extinctions of the acid and basic thymol blue solutions of $c_0 = 3 \cdot 10^{-5}$ molar at the same wavelength.

With c_{A^-} known, then c_{HA} is accessible via equation (4). Insertion of the values for c_{A^-} and c_{HA} the concentration determined by the buffer of $c = 1 \cdot 10^{-9} \text{ mol} \cdot \text{l}^{-1}$ in the law of mass action (1) gives the indicator constant K_a .

Evaluation (8/13)

PHYWE

Evaluation (7/8)

For the compound examined here, an alternative evaluation procedure can be used. As at wavelengths above 625 nm practically only the conjugated base A^- absorbs ($\epsilon_{HA} = 0$), equation (6.1) changes to:

$$E_\lambda = \epsilon_{A^-} \cdot c_{A^-} \cdot d \quad (6.1.3)$$

The quotient from the relationships (6.1.3) and (6.1.2) is, according to definition, equal to the degree of dissociation α ,

$$\alpha = \frac{c_{A^-}}{c_0} = \frac{E_\lambda}{A_{\lambda, A^-}} \quad (8)$$

Evaluation (9/13)

PHYWE

Evaluation (8/8)

This can be calculated from the extinctions E_λ and E_{λ,A^-} measured in buffered solution or NaOH at constant wavelength. The indicator constant K_a is then accessible from the relationship derived from the law of mass action (1) and the relationships (4) and (8):

$$K_a = \frac{c_{H^+} \cdot \alpha}{1 - \alpha} \quad (9)$$

The negative decadic logarithm of its numerical value is, according to definition (2), equal to the pK_a value.

Evaluation (10/13)

PHYWE

Data and Results

The absorption spectra recorded at $T = 299 \text{ K}$ are shown in the graphic. From the extinctions measured at $\lambda = 460 \text{ nm}$ for the acid solution ($E_{\lambda,HA} = 0.319$), the basic solution ($E_{\lambda,A^-} = 0.073$) and the buffered solution ($E_\lambda = 0.195$), the acid constant can be obtained from the relationships (6.1.1), (6.1.2), (7), (4) and (1) given in the text above, and we obtain $K_a = 1.02 \cdot 10^{-9} \text{ mol} \cdot \text{l}^{-1}$, from which, using (2), it follows that $pK_a = 8.99$ (evaluation procedure 1). The extinctions $E_{\lambda,A^-} = 0.600$ and $E_\lambda = 0.300$ can be taken from the spectra at $\lambda = 460 \text{ nm}$. From these and acc. to (8), a degree of dissociation of $\alpha = 0.50$ is given, which can be inserted in (9) to give $K_a = 1.00 \cdot 10^{-9}$ ($pK_a = 9.00$) (evaluation procedure 2).

The average value of several calculations at various wavelengths is $pK_a = 9.01$. This lies in the region of values given in the literature, that fluctuate around $pK_a = 9.00$ ($T = 293 \text{ K}$).

Evaluation (11/13)

PHYWE

Drag the words into the correct boxes!

The coloured [] thymol blue is a weak acid that is partially dissociated in aqueous solution, whereby non-ionized and ionized forms show absorption maximums at different [] in the visible range. Photometric measurements in the visible spectral range can therefore be used to advantage to determine the position of the K_a and pK_a values of the indicator which characterize [].

dissociation equilibrium

wavelengths

indicator

☒ Check

Evaluation (12/13)

PHYWE

Drag the words into the correct boxes!

In buffered solutions of $\text{pH} \approx 9$ pK_a the [] and non-ionized forms are present in practically the same []. These equilibrium concentrations, and so the constants K_a and pK_a , for [], can be advantageously measured via photometric measurements as, because of their different atomic structures, the acid and salt forms give different [] spectra which cross each other at an isosbestic point.

concentration

absorption

thymol blue

ionized

☒ Check

Evaluation (13/13)

PHYWE

What is given by the Lambert-Beer law?

- ☐ The Lambert-Beer law gives you nothing.
- ☐ Non of the answers is correct.
- ☐ The Lambert-Beer law gives the pH-Value of the buffered solutions.
- ☐ The dependence on the concentration c_i of an absorbing substance i and the layer thickness d at constant wavelength is given by the Lambert-Beer law.

✓ Überprüfen

Slide	Score/Total
Slide 22: Summary of the experiment	0/3
Slide 23: Equilibrium concentrations	0/4
Slide 24: Lambert-Beer law	1/1

Total Score  1/8

👁 Show solutions

🔄 Retry