

# Photosynthesis and cellular respiration in plants with Cobra SMARTsense



Biology

Plant Physiology / Botany

Photosynthesis



Difficulty level

medium



Group size

-



Preparation time

10 minutes



Execution time


30 minutes

This content can also be found online at:



<http://localhost:1337/c/5f45a80a9a658b00033e035f>

PHYWE



## General Information

### Application

PHYWE



Experimental setup with reaction chamber, plant and sensor

Plants are known to produce carbon dioxide ( $\text{CO}_2$ ) and water ( $\text{H}_2\text{O}$ ) to oxygen ( $\text{O}_2$ ) and glucose ( $\text{C}_{12}\text{O}_6\text{H}_{12}$ ) to be implemented. This is called photosynthesis, which means "light composition". As this translation already suggests, this type of metabolism only works with the help of light.

So that plants do not "suffocate" in the dark, it is important that they continue to convert substances into others to produce energy. This is why the reverse reaction occurs without light; plants breathe cells just like animals.

This experiment serves to study photosynthesis and cellular respiration in plants.

## Other information (1/2)

PHYWE

### Prior knowledge



Plants have two types of gas exchange: they carry out photosynthesis when exposed to light and cell respiration when it is dark.

### Scientific principle



This experiment shows that plants produce oxygen in light but consume it in darkness.

## Other information (2/2)

PHYWE

### Learning objective



The aim of this experiment is to investigate the metabolism of plants with regard to the oxygen concentration in a hermetically sealed vessel.

### Tasks



The pupils and students should carry out the following partial experiments in this experiment:

- Determination of the O<sub>2</sub>-production during photosynthesis
- Determination of the O<sub>2</sub>-consumption during cellular respiration

Alternatively, the experiment is very well suited as a demo experiment.

## Safety instructions

PHYWE

The general instructions for safe experimentation in science lessons to be applied to this experiment.

## Theory

PHYWE

Plants, like all living things, are made up of cells. These cells need a constant influx of energy to maintain their metabolism and not die.

During the day, the energy of the sun is used for this purpose, to absorb CO<sub>2</sub> and water to O<sub>2</sub> and dextrose. The corresponding reactions take place in the chloroplasts, where the light of the sun is captured by chlorophyll, the green pigment of the plants, and other colour pigments such as carotenoids, and its energy is used to produce 6 x CO<sub>2</sub> and 6 x H<sub>2</sub>O to 6 x O<sub>2</sub> and 1 x C<sub>6</sub>H<sub>12</sub>O<sub>6</sub> (dextrose). So the energy is actually used to produce glucose, which is then used to produce energy.

In addition, plants breathe cells, which also occurs in animals. To a certain extent, cell respiration always takes place, only that during the day more O<sub>2</sub> as CO<sub>2</sub> is produced. In this process, glucose is "burned" with oxygen to form carbon dioxide and water. All this happens in the mitochondria of plants.

## Equipment

Position	Material	Item No.	Quantity
1	<a href="#">Cobra SMARTsense Oxygen - Sensor for measuring the oxygen content 0 ... 20 mg/l (Bluetooth + USB)</a>	12933-01	1
2	<a href="#">Experiment chamber, 29 cm (11.4"), for Cobra SMARTsense sensors</a>	64837-00	1
3	<a href="#">measureAPP - the free measurement software for all devices and operating systems</a>	14581-61	1
4	<a href="#">Cobra SMARTsense CO2 - Sensor for measuring the carbon dioxide content 0 ... 100000 ppm (Bluetooth + USB)</a>	12932-01	1

PHYWE

## Set-up and procedure



### Set-up (1/2)

PHYWE



Experimental setup with reaction chamber, plant and sensor

First the reaction chamber is built up by inserting the rubber plugs into the corresponding holes. One of these rubber plugs has a small hole through which the probe of the measuring instrument can be inserted. Afterwards the plant is placed in the chamber.

When all this is done, start the measureAPP and connect the sensor to a mobile device or a laptop with Windows 10, close the chamber and the measurement can start.

Note: For significant results, both parts of the experiment should last at least 30 minutes and be of approximately equal length.

## Set-up (2/2)

PHYWE

The Cobra SMARTsense and measureAPP are required for measuring the oxygen concentration in air. Check that "Bluetooth" is enabled on your device (tablet, smartphone, PC with Windows 10) (the app can be downloaded for free from the App Store - QR codes below). Now open measureAPP on your device.



measureAPP for

Android operating systems



measureAPP for

iOS operating systems



measureAPP for

Tablets and PCs with Windows 10

## Procedure (1/2)

PHYWE



Darkened reaction chamber



It is important that the reaction vessel is sealed airtight. If the rubber stopper and lid are not tight enough, you can help with adhesive tape or similar.

The measurement should be set to continuous and it is recommended to set the sensor to zero.

For the second part of the measurement, the chamber should be darkened as shown in the picture on the left. This serves to measure the oxygen consumption as it takes place at night.


## Procedure (2/2)

When the measurements are complete, the results should look something like this. Left is the measurement in light, right is the dark measurement. To calculate theoretically the total oxygen production, one would have to add the oxygen consumed in the dark, because it is also consumed in the light.



## Evaluation (1/3)

PHYWE



Olive leaf under water

Why do plants also need to breathe?

Because otherwise not enough energy would be available at night, as photosynthesis cannot be carried out.

So that not too much O<sub>2</sub> is produced, which is toxic for plants.

Because plants are derived from animals, they can switch when needed.

To avoid collapsing under the load of too much stored glucose.



## Evaluation (2/3)

PHYWE

Plants always breathe

☐ True☐ Wrong☒ Check

Plants always carry out photosynthesis

☐ True☐ Wrong☒ Check

## Evaluation (3/3)

PHYWE

In which part of the cell does photosynthesis take place?

In chlorophyll

In the mitochondria

In the cell nucleus

In chloroplast

In which part of the cell does cellular respiration take place?

In the plasma membrane of the cell


In the mitochondria

In the cell nucleus

In the vacuole

Slide	Score / Total
Slide 13: Metabolism	0/1
Slide 14: Multiple tasks	0/2
Slide 15: Multiple tasks	0/2

Total points



Show solutions



Repeat