


Propagation of sound in solid bodies



Physics

Acoustics

Sound generation & propagation



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time


10 minutes

This content can also be found online at:



<http://localhost:1337/c/616415b2374f4f00038cbacd>

PHYWE



Teacher information

Application

PHYWE



Experiment setup

Transmission (or conduction) of sound through a solid means: one end of the solid is set into vibration (e.g. by the application of a tuning fork) and this vibration propagates through the solid to the other end of the solid. In this experiment, the transmission through the solid is heard as a vibration of the other end directly or through the air.

The experiments will also show that at greater distances (approx. 20 cm) the sound of the tuning fork can no longer be heard through the air, or can only be heard very softly. However, via conduction through solids, the sound is transmitted so well at comparable distances that it can be heard clearly.

Other teacher information (1/3)

PHYWE

Previous knowledge



Students should be familiar with the basic concepts of waves and oscillations.

Scientific principle



The propagation of sound in solids is investigated using three examples:

- Ruler (comparison with railway tracks)
- Conduction over the head (bone conduction)
- Silk thread (comparison with "cord telephone")

Other teacher information (2/3)

Learning objective



Students learn how sound waves propagate through a solid.

Tasks



Sound can be generated, for example, by the vibrations of a tuning fork and normally reaches the ear through the air. Here, the pupils use various examples to investigate whether sound propagates not only in air but also in solid bodies.

1. Is sound transmitted through a ruler made of plastic?
2. Can the human head transmit sound to the ear?
3. Is a thin thread sufficient for sound transmission?

Other teacher information (2/3)

PHYWE

Learning objective



Students learn how sound waves propagate through a solid.

Tasks



Sound can be generated, for example, by the vibrations of a tuning fork and normally reaches the ear through the air. Here, the pupils use various examples to investigate whether sound propagates not only in air but also in solid bodies.

1. Is sound transmitted through a ruler made of plastic?
2. Can the human head transmit sound to the ear?
3. Is a thin thread sufficient for sound transmission?

Other teacher information (3/3)

PHYWE

Notes on structure and implementation

The sound is transmitted through the vibrating lower end (the ball) of the tuning fork to the solid body (e.g. the ruler, the head, the thread or a table).

The easiest way to achieve this in the sound transmission in experiment 3 is to have the tuning fork hanging down by the thread.

Safety instructions


PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE

Student Information



Motivation

PHYWE



A corded telephone

If you connect two cans that are open on one side with a cord on their closed side, you can talk into one can and be understood at the other can. This so-called cord telephone enables communication over the length of the cord, even over distances where normal conversation would no longer be possible.

But how is it possible that you speak into a closed container and the words are still transmitted more understandably than if you speak directly to the other person over the air?

Equipment

| Position | Material | Item No. | Quantity |
|----------|--|----------|----------|
| 1 | Ruler, plastic, 200 mm | 09937-01 | 1 |
| 2 | Tuning fork 440 Hz | 03424-00 | 1 |
| 3 | Impact hammer, rubber | 03429-00 | 1 |
| 4 | Silk thread, l = 200 m | 02412-00 | 1 |
| 5 | Funnel, diameter = 75 mm, plastic (PP) | 46895-00 | 1 |

Equipment

PHYWE

| Position | Material | Item No. | Quantity |
|----------|--|----------|----------|
| 1 | Ruler, plastic, 200 mm | 09937-01 | 1 |
| 2 | Tuning fork 440 Hz | 03424-00 | 1 |
| 3 | Impact hammer, rubber | 03429-00 | 1 |
| 4 | Silk thread, l = 200 m | 02412-00 | 1 |
| 5 | Funnel, diameter = 75 mm, plastic (PP) | 46895-00 | 1 |

Set-up

PHYWE

1. Get all the materials you need to experiment at your workstation.
2. Now read carefully the implementation.

Procedure (1/5)

PHYWE



Figure 1

Experiment with ruler

1. Strike the tuning fork using the striking hammer (Fig. 1).
2. First hold the tuning fork in the air and then move it away from your ear until you can no longer hear it. Estimate the distance you have reached and write it down.

Procedure (2/5)

PHYWE



Figure 2

3. Strike the tuning fork again.
4. Place it on one end of the ruler and then bring the other end to your ear (Fig. 2).
5. Compare your observations from both experiments.

Implementation (3/5)

PHYWE



Figure 3

Experiment: Put tuning fork on the head

1. Strike the tuning fork using the striking hammer.
2. Place the tuning fork quickly one after the other at different places on your head (Fig. 3): at the top, at the back, directly at the ear,
3. Describe how the tone of the tuning fork changes when you put it on your head.
4. Note where on your head you need to hold the tuning fork to hear the sound loudest.

Procedure (4/5)

PHYWE



Figure 4

Experiment with silk thread

1. First connect the funnel and the tuning fork with a silk thread (approx. 1 m) as follows: Pass the thread through the funnel from below and fasten it to the eyelet on the upper edge with a knot. (Fig. 4). Attach the tuning fork to the other end of the thread with another knot.
2. Strike the tuning fork with the hammer and hold it as far away from you as possible without stretching the thread. Note whether you can still hear the tuning fork.

Procedure (5/5)

PHYWE




Fig. 5

3. Then strike the tuning fork again with the striking hammer.
4. Let the tuning fork hang loosely from the thread and hold the funnel to your ear (Fig. 5). Note whether you can hear the tuning fork.

Report

PHYWE

Report



Task 1

PHYWE

Drag the words into the correct gaps

How well a sound wave is transmitted through a medium depends on its . ρ In general, solids are denser than , which in turn are denser than gases.

If a sound wave can propagate through a denser medium, it is intelligible to the human over greater distances. A good example of this is whale song, which can be heard in even over hundreds of miles.

ear

liquids

density

water

☒ Check

Task 2

PHYWE

What describes the physical quantity density ρ ?

- ☐ The physical density ρ describes the mass m per volume V , so $\rho = \frac{m}{V}$.
- ☐ The physical density ρ describes the product between mass m and volume V , so $\rho = m \cdot V$.
- ☐ The physical density ρ describes the number of particles n per volume V , so $\rho = \frac{n}{V}$.
- ☐ The physical density ρ describes the force F per area A , so $\rho = \frac{F}{A}$.

☒ Check