Rusting- "slow combustion"



In this student experiment, the example of "rusting iron wool" clearly shows that oxidation also takes place over longer periods of time and in an "unspectacular" way. The experiment shows that iron corrodes and rusts when it is exposed to humid air for a long time.

Chemistry	Inorganic chemistry	Acids, ba	ses, salts
Difficulty level easy	QQ Group size	Preparation time	Execution time
This content can also be found online at:	回 然 梁父		

http://localhost:1337/c/633b1a8dbdf6080003a941d0





Teacher information

Application

PHYWE



Rusted iron chain

Rust is formed from iron or steel by oxidation with oxygen in the presence of water.

Rust is a hydrous oxide of iron, a chemical compound that belongs to the oxides and also contains water and hydroxide ions. High temperatures are not required.

Rust is porous and does not protect against further decomposition, unlike the oxide layer of many metallic materials such as chrome, aluminium or zinc.

The weathering of ferrous materials in air and water to rust causes damage worldwide every year.







Safety instructions

PHYWE



- Use protective goggles!
- The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE



Student information



Motivation

PHYWE



Wire brushes can clean rusted iron parts

Every iron or steel component is affected by rust in the long term. It is not possible to remove the rust, as it could be counterproductive in the long term. To restore affected iron parts to a good condition, you can remove the rust with a wire brush or sandpaper.

In industry, sandblasters are used to remove paint residues and rust. Usually a protective layer prevents direct contact with oxygen and water. This can consist of paint, grease, oil or another metal. For bicycle chains, lubricants are used to protect against corrosion and wear.

Tasks

PHYWE

 Investigate the rusting of iron wool. 	What processes take place during rusting?
 Note down your observations and answer the questions in the report. 	Ferrous metals rust. False True



Equipment

Position	Material	Item No.	Quantity
1	Beaker, Borosilicate, tall form, 250 ml	46027-00	1
2	Test tube, 180x18 mm,100pcs	37658-10	1
3	Test tube brush w. wool tip,d20mm	38762-00	1
4	Protecting glasses, clear glass	39316-00	1
5	Glass rod, boro 3.3, I=200mm, d=5mm	40485-03	1
6	Iron wool 200 g	31999-20	1
7	Butane burner with cartridge, 220 g	32180-00	1
8	Wood splints, package of 100	39126-10	1



Procedure (1/2)

PHYWE

Loosely insert a suitable iron wool ball into the test tube and push it to the bottom with the glass rod without squeezing it.

Moisten the iron wool slightly and fill the beaker halfway with water.



Procedure (2/2)

PHYWE

Put the test tube in with the opening facing down.

Place the beaker with contents in a storage place (cupboard) where it can remain for a week.

Place some iron wool next to the beaker for comparison.

In the next week, remove the test tube, closing the opening with your thumb.

Light a wood chip and hold it in the test tube.















Task 2	PHYWE
	In winter, the rusting of vehicles is accelerated by Road salt and water Sugar Nitrogen

Task 3

PHYWE

How can iron parts be protected against rusting?		
is the most widely produced and used metal. Whe	en and ox	ygen Iron
attack together, it rusts. The surface must be thoroughly	and the first	cleaned
must be applied immediately afterwards, otherw	ise the iron will	water
again straight away. The paints contain pigment and fillers. The pigments the		
material and protect it from mechanical influences and	new rust formatio	on. prevent
		paint
		colour



Slide	Score/Total
Slide 8: Properties of metals	0/2
Slide 14: Grate	0/4
Slide 15: Corrosion protection	0/7
	Total 0/13
Solutions	t Export text

