

Water, an oxide



Chemistry

Inorganic chemistry

Water



Difficulty level

easy



Group size

2



Preparation time

10 minutes



Execution time

10 minutes

This content can also be found online at:

<http://localhost:1337/c/634122751fa0ea00031e40ce>

PHYWE

Teacher information



Application

PHYWE



The experimental setup

Water plays an essential role in everyday life. Water is also used in many ways in chemistry. The question of whether pure water is an element, a pure substance or a compound can be answered experimentally by the fact that when substances containing hydrogen (hydrocarbons) are burnt, water is produced in addition to carbon dioxide.

The experiment then makes it clear that water is a compound, namely an oxide.

Other teacher information (1/3)

PHYWE

Prior knowledge



The students should have knowledge of the periodic table. Furthermore, the students should already know the terms reduction and reducing agent, as well as the difference between elements and compounds, which can be clearly understood and deepened through this experiment.

Students should know the basics of working with chemicals and be able to work with a butane gas or Bunsen burner.

Principle



In the experiment, the students observe the combustion of substances containing hydrogen and conclude from their observations which processes produce water.

Other teacher information (2/3)

PHYWE

Learning objective



1. The students learn that water is not an element, but a compound.
2. Since water is formed by the combustion of hydrogen-containing substances, water is an oxide.

Tasks



1. Carry out burns of a) butane or natural gas and b) methylated spirits. Write down the observations.
2. Infer from the observations whether water is a compound or an element and how it is formed.

Other teacher information (3/3)

PHYWE

Notes on the set-up and procedure

- Instead of ice water, you can also use chilled water (from the refrigerator). When using water at room temperature as a coolant, the condensation effect is very low.
- The use of a water jet pump is not absolutely necessary. Combustion gas can also be sucked in with a pipetting ball or peplus ball or negative pressure can be dispensed with altogether. Without negative pressure, however, only the smallest amounts of condensed water are produced, which are hardly visible to the naked eye, but would require detection reactions, e.g. with copper sulphate.

Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

Methylated spirits are highly flammable. Extinguish all open flames!


Put on protective goggles!

Make rubber-glass joints slippery with glycerine. Do not use force!

For H and P phrases, please refer to the safety data sheet of the respective chemical.

PHYWE

Student information



Motivation

PHYWE



Water: an element?

Whether in a lake, river, the sea, rain, snow or morning dew, water is omnipresent to us. Even our bodies consist to a large extent of water. People were already aware of the importance of water in ancient times, and since then it has been one of the classical elements, along with air, fire and earth.

But is water really an element or is it a compound?

Tasks

PHYWE



The experimental setup

1. Carry out burning of

- Butane or natural gas
- Methylated spirits.

Note down the observations.

2. Deduce from the observations whether water is a compound or an element and how it is formed.

Equipment

| Position | Material | Item No. | Quantity |
|----------|---|--------------|----------|
| 1 | Support base, variable | 02001-00 | 1 |
| 2 | Support rod, stainless steel, l=370 mm, d=10 mm | 02059-00 | 1 |
| 3 | Boss head | 02043-00 | 2 |
| 4 | Porcelain dish, 75ml, d = 80 mm | 32516-00 | 1 |
| 5 | Crucible tongs, 200 mm, stainless steel | 33600-00 | 1 |
| 6 | Funnel, glass, top dia. 80 mm | 34459-00 | 1 |
| 7 | Beaker, Borosilicate, tall form, 250 ml | 46027-00 | 1 |
| 8 | Test tube, 180x20 mm, side arm, PN19 | MAU-17080301 | 1 |
| 9 | Graduated cylinder, 10 ml, plastic | 36636-00 | 1 |
| 10 | Glass tubes, right-angled | MAU-10030701 | 1 |
| 11 | Test tube rack f. 6 tubes, wood | 37685-10 | 1 |
| 12 | Universal clamp | 37715-01 | 2 |
| 13 | Test tube brush w. wool tip, d20mm | 38762-00 | 1 |
| 14 | Rubber stopper, d = 22/17 mm, 1 hole | 39255-01 | 1 |
| 15 | Rubber tubing, i.d. 6 mm | 39282-00 | 1 |
| 16 | Protecting glasses, clear glass | 39316-00 | 1 |
| 17 | Glycerol, 250 ml | 30084-25 | 1 |
| 18 | Denaturated alcohol (spirit for burning), 1000 ml | 31150-70 | 1 |
| 19 | Wood splints, package of 100 | 39126-10 | 1 |
| 20 | Butane burner with cartridge, 220 g | 32180-00 | 1 |
| 21 | Water jet pump, plastic | 02728-00 | 1 |
| 22 | Rubber tubing, vacuum, i.d. 6mm | 39286-00 | 1 |

Set-up (1/2)

PHYWE

1. Assemble the support with the support base, the support rod and the two sleeves according to Fig. 1 - Fig. 3.

2. Then attach the universal clamps at right angles to each other (Fig. 4).



Figure 1




Figure 2



Figure 3



Figure 4

Set-up (2/2)

PHYWE

3. Put a drop of glycerine on the end of the angle tube and turn it carefully and deeply into the hole of the stopper (Fig. 5). Do not use force!

4. Close the test tube with the stopper (Fig. 6) and clamp it and the trichter each in a universal clamp.

Connect the glass funnel and the angled tube with a piece of tubing to form a closed apparatus. (Fig. 7)

5. Immerse one third of the test tube in a beaker filled with ice water (Fig. 8).




Figure 5




Figure 6



Figure 7



Figure 8

Procedure (1/2)

PHYWE



Figure 9



Figure 10

1. Connect the water jet pump via the vacuum hose to the attachment nozzle of the test tube (Fig. 9) and set a mild negative pressure.
2. Hold the pilot light of the burner under the funnel and let the combustion products flow through the apparatus for about 5 minutes (Fig. 10).
3. Pour out the ice water and dry the beaker. Extinguish the burner flame.

Procedure (2/2)

PHYWE



Figure 11




Figure 12

4. Fill the evaporating dish with approx. 3 ml methylated spirits, remove the filling bottle and the measuring cylinder.
5. Light the wood chip, extinguish all other flames and light the methylated spirit with the wood chip (Fig. 11).
6. Hold the bottom of the dry and still cool beaker over the methylated spirit flame with the crucible tongs (Fig. 12).

PHYWE

Report



Task 1

PHYWE

Write down your observations.

Bunsen burner flame:

Methylated spirits flame:

Task 2

PHYWE

According to the results of the experiment, to which substance class does water belong?

☐ Oxides☐ Noble gases☐ Halogens☐ Alkanes

Which of these substances is found among the products of these oxidation reactions?

☐ Water☐ Ammonia☐ Helium☐ Carbon monoxide

Task 3

PHYWE

Water is a chemical element.

☐ True☐ False☒ Check

The chemical formula of water is H_2O .

☐ True☐ False☒ Check

| Slide | Score / Total |
|--------------------------|---------------|
| Slide 17: Multiple tasks | 0/2 |
| Slide 18: Multiple tasks | 0/2 |

Total  0/4



Solutions



Repeat



Export text