


The cracking of liquid paraffin



Chemistry

Industrial Chemistry

Petrochemistry



Difficulty level

easy



Group size

2



Preparation time

10 minutes



Execution time


10 minutes

This content can also be found online at:



<http://localhost:1337/c/638b6a558ff8b60003be7e2d>

PHYWE



Teacher information

Application

PHYWE



The experimental setup

Crude oil is a mixture of more than 500 different organic compounds.

Among them are many long-chain hydrocarbons, which are not as commercially in demand as short-chain hydrocarbons.

This is why the industry often uses oil cracking, the process by which long-chain alkanes are broken down into short-chain alkanes.

Other teacher information (1/5)

PHYWE

Prior knowledge



Students should already have basic knowledge about saturated and unsaturated hydrocarbons and the composition of petroleum.

Further, students should know the basics of working with chemicals and be able to work with a butane gas or Bunsen burner.

Principle



Paraffin oil is split by heat with the help of a catalyst into short-chain, mainly gaseous hydrocarbons. The proportion of unsaturated compounds increases greatly, which can subsequently be detected by the Baeyer reagent.

Other teacher information (2/5)

PHYWE

Learning objective



The students should learn that higher-boiling alkanes can be decomposed into lower-boiling hydrocarbons by catalytic cracking. The products of this process are alkenes and alkanes.

Tasks



Produce low-boiling hydrocarbons from paraffin oil.

Other teacher information (3/5)

PHYWE


Note on preparation

Prepare soda alkaline potassium permanganate solution. To do this, add a few small potassium permanganate crystals to 50 ml of water, add a spoonful of sodium carbonate.

Unsaturated hydrocarbons can also be detected by bromine water. Due to its carcinogenic properties, the Baeyer reagent is recommended for student experiments.

Note on the theory

In this experiment, the long-chain hydrocarbons of paraffin oil are split to short-chain alkenes and unsaturated hydrocarbons. For example:



Other teacher information (4/5)

PHYWE

Notes on the student experiments

Make sure that the bead catalyst is heated sufficiently. The paraffin oil may only be added drop by drop, otherwise the flask may crack due to the temperature differences that occur. If you use catalysts other than the bead catalyst, for example steel wool or activated carbon, these alternatives must be heated much more.

The burner must not be held to the outlet of the second test tube with preparation nozzle, as flammable gases may escape from there!

Other teacher information (5/5)

PHYWE

Methodological remarks

The test, which is quite demanding in terms of experiment, draws on the "ethene representation" test for the detection methods, which may have to be preceded by this test.


Bunsen burners can still be used instead of butane burners for this experiment.

Disposal

Put the contents of all containers into the collection container for combustible organic substances.

Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.


For H- and P-phrases please consult the safety data sheet of the respective chemical.

Dangers

- Paraffin oil is flammable. Extinguish all open flames during filling!
- Potassium permanganate solution is corrosive. Avoid contact with the skin! Wash off splashes immediately with plenty of water! Use protective gloves and wear safety goggles!
- Make rubber-glass compounds slippery with glycerine!

PHYWE

Student information



Motivation

PHYWE



A petrol pump

The industry's demand for short-chain alkanes is significantly greater than that for long-chain alkanes. This is mainly because short-chain alkanes are needed for a number of essential fuels such as petrol and diesel, while long-chain hydrocarbons are used almost exclusively as fuel oil.

For this reason, long-chain alkanes are often broken down into short-chain alkanes. This so-called "cracking" is examined in more detail in this experiment. Here, the process is supported by a catalyst (Al/Si oxides).

Tasks

PHYWE



The experimental setup

How can higher-boiling petroleum fractions be "refined"?

Produce low-boiling hydrocarbons from paraffin oil.

Equipment

Position	Material	Item No.	Quantity
1	Support base, variable	02001-00	1
2	Support rod, stainless steel, l=370 mm, d=10 mm	02059-00	3
3	Boss head	02043-00	3
4	Spoon, special steel	33398-00	1
5	Round bottom flask w. side arm, 100 ml, PN19	34885-00	1
6	Beaker, Borosilicate, low form, 250 ml	46054-00	1
7	Glass tubes, right-angled	MAU-10030701	1
8	Dropping funnel with drip nozzle, 50ml	36912-00	1
9	Test tube, 180x18 mm, 100pcs	37658-10	1
10	Test tube, 180x20 mm, side arm, PN19	MAU-17080301	2
11	Test tube brush w. wool tip, d20mm	38762-00	1
12	Universal clamp	37715-01	3
13	Rubber stopper, d = 22/17 mm, 1 hole	39255-01	3
14	Rubber tubing, i.d. 6 mm	39282-00	1
15	Protecting glasses, clear glass	39316-00	1
16	Rubber gloves, size M (8), one pair	39323-00	1
17		CHE-881317634	1
18	Glycerol, 250 ml	30084-25	1
19	Potassium permanganate, chem. pur., 250 g	30108-25	1
20	Sodium carbonate, anhyd. 1000 g	30154-70	1
21	Liquid paraffin 1000 ml	30180-70	1
22	Butane burner with cartridge, 220 g	32180-00	1
23	Water, distilled 5 l	31246-81	1

Set-up (1/6)

PHYWE

1. Set up the support according to Fig. 1 to 3.



Figure 1



Figure 2



Figure 3

Set-up (2/6)

PHYWE

2. Clamp the round piston on the left stand rod (Fig. 4).
3. Screw the two angled glass tubes into the pierced rubber stoppers (Fig. 5).



Figure 4




Figure 5

Set-up (3/6)

PHYWE

4. Carefully turn the outlet of the separating funnel as far as it will go into the pierced rubber plug SB 19 (**Attention:** Always lubricate the rubber stopper with glycerine) so that the spout protrudes approx. 5 cm through the rubber stopper (Fig. 6).



Figure 6

Set-up (4/6)

PHYWE

5. Fill the round bottom flask with the bead catalyst until the bottom is well covered (Fig. 7).

6. Close the test tube with the attachment nozzle with the rubber stopper of the angled glass tube (Fig. 8).

7. Connect the neck of the round-bottomed flask to the short leg of one of the angled glass tubes using a piece of tubing and insert them into the test tube.




Figure 7




Figure 8

Set-up (5/6)

PHYWE



Figure 9

8. Close the round-bottomed flask with the rubber stopper on the separating funnel and fill the separating funnel with approx. 5 ml paraffin oil with the tap closed (Fig. 9 and 10).



Figure 10

Set-up (6/6)

PHYWE

9. Connect the neck of the first test tube to the angled glass tube of the second test tube with neck using a piece of tubing.

10. Place the second test tube with the preparation spout in the beaker and fill it with ice (Fig. 11).

11. Attach the test tube with the attachment nozzle to the second stand rod (Fig. 12).



Figure 11



Figure 12

Procedure (1/2)

PHYWE



Figure 13

1. Heat the bead catalyst.
2. Slowly open the tap of the separating funnel so that only a few drops of paraffin oil flow out.
3. Then let the paraffin oil flow in drop by drop over several minutes (Fig. 13).

Procedure (2/2)

PHYWE




Figure 14

4. Finish the experiment and let the apparatus cool down.
5. Detach the second test tube with the preparation nozzle from the apparatus, remove the rubber stopper and add a few drops of soda-alkaline potassium permanganate solution to the condensate (Fig. 14).

PHYWE

Report



Task 1

PHYWE

Write down your observations.

Task 2

PHYWE

Fuels

short-chain alkanes

long-chain alkanes

Petrol

Fuel oil (heavy)

Diesel

Marine fuel

✓ Check

Task 3

PHYWE

Drag the words into the correct boxes!

Long-chain can be broken down into more commercially valuable, short-chain alkanes and by catalytic .

This process usually takes place at very high . Catalytic cracking also uses a that breaks down the long-chain hydrocarbons into short-chain alkanes and alkenes over many reactions.

hydrocarbons

catalyst

alkenes

temperatures

cracking

✓ Check

Task 4

PHYWE

Drag the words into the correct boxes!

In catalytic cracking, are used to bring about a series of complicated reactions.

The of the intermediate steps are usually very reactive and either decompose spontaneously or react further with other intermediates. such as aluminium oxides and silicates are mainly used.

products

catalysts

Aluminium compounds

 Check

Slide

Score/Total

Slide 23: Fuels

0/4

Slide 24: Catalytic cracking

0/5

Slide 25: Cracking catalysts

0/3

Total

  0/12 Solutions Repeat Export text