

From seed to plant



Biology	Plant Physiology / Botany	Physiology of	plants
Biology	Plant Physiology / Botany	Reproduction	in plants
Biology	Plant Physiology / Botany	Germination, growth, development	
Difficulty level	QQ Group size	Preparation time	Execution time
easy	1	10 minutes	10 minutes

This content can also be found online at:



http://localhost:1337/c/5f0d6cd71c41060003916b75





PHYWE



Teacher information

Application PHYWE



It is such a self-evident process in nature that you don't perceive it anymore: seeds, sown in orchards, fields and forests, sprout and grow to become fully developed plants.



Other teacher information (1/2)

PHYWE

Prior knowledge



Scientific principle



While there are kinds of seed, which only germinate in bright light, beans need darkness to sprout. That is why it is vital to encase the beaker thoroughly with black carton. Subsequently the beaker should be located in a bright spot, e. g. in near proximity of a window, to accelerate the plant's growth.

This experiment underlies the common principles of plant growth.

To achieve valid results, the experiment must last at least 10 days.

Other teacher information (2/2)



Learning objective



Tasks



The pupils will experience the development from seed to plant in this experiment.

The students are to plant the beans and observe their further development.



Safety instructions

PHYWE



The common rules of safe experimentation in scientific education to be applied to this experiment.

PHYWE




Student information





Motivation PHYWE



You see it so frequently, you don't realize it anymore with full awareness: small, inconspicuous seeds transform into big trees, pretty flowers, or delicious vegetables. But what happens exactly when a seed becomes a plant? This experiment will show.

Tasks PHYWE

In this experiment you are going plant beans and observe their development. What happens before the plant appears on the surface?




Seed after germination





PHYWE



In this experiment you are going plant beans and observe their development. What happens before the plant appears on the surface?





Equipment

Position	Material	Item No.	Quantity
1	Rubber bands, 50 pieces	03920-00	1
2	Beaker, 250 ml, plastic (PP)	36082-00	1
3	Beaker, Borosilicate, tall form, 600 ml	46029-00	1
4	Petri dish. d 100 mm	64705-00	1





Procedure (1/3)

PHYWE



At first place 10 beans in a petri dish.



Cover them with water afterwards in order that the beans swell.

Procedure (2/3)

PHYWE



Now fill a 600 ml beaker with soil and plant the swollen beans in it, leave about 2 cm free from the top. You need to spread the beans evenly at the edge of the beaker so you can see them grow.



Then encase the beaker with black carton as high as the soil reaches. Fix it with a rubber band.





Procedure (3/3)



Place the beaker in a bright space, e. g. near a window, water it regularly (but not too much) and observe the beaker for at least 10 days.

In order to see what happens beneath the surface, you have to remove the carton once in a while. Encase it again carefully afterwards.


PHYWE




Report













Task 3

Drag the words into the correct boxes

In the beginning the need to swell so that the boll can burst open. Afterwards the beans must be planted in soil in order that and stem can develop.

When the plant appears on the surface it is called , it then begins to grow its first leaves.

beans seedling root

Check



Slide	Score / Total
Slide 15: Colour of the stem	0/1
Slide 16: Multiple tasks	0/2
Slide 17: Bean seed seedlings	0/3
	Total Score 0/6

Show solutions



