


Influence of surface on the absorption of solar energy



Physics

Energy

Renewable energies: Sun



Difficulty level

easy



Group size

1



Preparation time

10 minutes



Execution time

10 minutes


This content can also be found online at:



<http://localhost:1337/c/62e8112c7db41d0003cf9584>

PHYWE

Teacher information



Application

PHYWE



Experimental setup

In this experiment, the white tent in the desert is simulated by a white absorber plate that is illuminated by a halogen lamp for 5 minutes.

The temperature difference before and after irradiation is compared with that of a black absorber plate.

The result is that the white panel absorbs less solar energy and thus heats up less.

This result is to be transferred by the students to practical examples, such as the white tent of the Bedouin and the black skin of the polar bear.

Other teacher information (1/3)

PHYWE

Prior knowledge



Students should be proficient in the use of a power supply unit.

Principle



In this experiment, the absorption behaviour of surfaces is observed and the observed differences are traced back to their physical basis.

Other teacher information (2/3)

PHYWE

Learning objective



The students learn how the absorption behaviour of solar energy depends on the surface.

Tasks



This experiment investigates how black and white absorber plates behave when irradiated with a halogen lamp.

Other teacher information (3/3)

PHYWE

Notes on set-up and procedure

The measured initial temperatures of the absorber plates differ depending on the ambient temperature.

Different results in the temperature increase can be caused by reading inaccuracies on the thermometer or by different orientations of the absorber plate.

Safety instructions

PHYWE



The general instructions for safe experimentation in science lessons apply to this experiment.

PHYWE



Student information

Motivation

PHYWE



A sleeping polar bear

A well-known everyday wisdom says that you should not wear black clothes in the hot summer.

Based on the same physical conditions, polar bears have black skin under their fur and Bedouins prefer to live in white tents.

These circumstances are explained by the absorption behaviour of different surfaces.

In this experiment, a special focus is placed on the colour of a surface and its effect on absorptivity.

Tasks

PHYWE



The experimental setup

This experiment investigates how black and white absorber plates behave when irradiated with a halogen lamp.

Equipment

| Position | Material | Item No. | Quantity |
|----------|--|----------|----------|
| 1 | Halogen lamp with reflector, 12V / 20W | 05780-00 | 1 |
| 2 | Mount for halogen lamp with reflector | 05781-00 | 1 |
| 3 | Solar collector for student experiments | 05760-00 | 1 |
| 4 | Support base, variable | 02001-00 | 1 |
| 5 | Slide mount for optical bench | 09822-00 | 1 |
| 6 | Measuring tape, l = 2 m | 09936-00 | 1 |
| 7 | Support rod, l = 600 mm, d = 10 mm, split in 2 rods with screw threads | 02035-00 | 2 |
| 8 | Digital stopwatch, 24 h, 1/100 s and 1 s | 24025-00 | 1 |
| 9 | Lab thermometer, -10..+110 °C | 38056-00 | 1 |
| 10 | PHYWE Power supply, 230 V, DC: 0...12 V, 2 A / AC: 6 V, 12 V, 5 A | 13506-93 | 1 |

Set-up (1/2)

PHYWE

1. First screw the two-part stand rods together to form two long rods (Fig. 1).

Assemble the stand bench from the variable stand foot and the two rods (Fig. 2 and Fig. 3).



Figure 1



Figure 2



Figure 3

Set-up (2/2)

PHYWE

3. Clamp the lamp in the left part of the stand base and connect it to the power supply unit (12 V~) (Fig. 4).

The power supply unit is switched off.

4. Put the black absorber plate on the stem (Fig. 5). Attach the stem with plate in the tab and place it on the stand bench (Fig. 6).

5. Insert the thermometer into the measuring socket of the plate, move the tab until the distance between the lamp and the plate is 12 cm and align the plate parallel to the tab (Fig. 7).



Figure 4



Figure 5



Figure 6



Figure 7

Procedure

PHYWE




Figure 8




Figure 9

1. Observe the temperature before starting the measurement and wait until it stops changing.
2. Measure the initial angular temperature ϑ_{1black} of the black plate. Enter this under Result in the report.
3. Switch on the lamp (the power supply unit) and start the stopwatch at the same time (Fig. 8).
4. After 5 minutes, switch off the power supply and measure the final temperature of the plate. ϑ_{2black} and note them down.
5. Repeat the experiment with the white plate (Fig. 9).
6. Note the initial temperature ϑ_{1white} and the final temperature ϑ_{2white} .

PHYWE

Report



Task 1

PHYWE

Drag the words into the correct gaps

Light waves come in different [] depending on their energy level.

The higher the [], the smaller the wavelength. When a [] of light is reflected into our eye, our brain interprets the wavelengths as different

[].

Whereby [] is the combination of all colours and [] is the absence of any colour.

colours

energy

black

white

wavelengths

beam

 Check

Task 2

PHYWE

Decide whether or which bracketed word does not belong there.

White materials (reflect / absorb) all wavelengths and the energy they carry is accordingly (not) absorbed.

This causes white materials to heat up (slower / faster).

Black materials (reflect / absorb) all wavelengths and appear black precisely because no light falls from them into our eye. Due to the (reduced / additional) energy absorption, they heat up significantly (more / weaker).

 Check

Task 3

PHYWE

Find all true statements about light

- ☐ During absorption, the energy of the light is converted into kinetic energy.
- ☐ The more wavelengths a material reflects, the faster it heats up due to light irradiation.
- ☐ When light is reflected, the angle of incidence is exactly the same as the angle of reflection.
- ☐ During absorption, the energy of the light is absorbed by the electrons of the material.

 Check

Slide

Score / Total

Slide 15: What is a colour?

0/6

Slide 16: Reflection and absorption

0/6

Slide 17: Light

0/3

Total   0/15 Solutions Repeat