



# leXsolar-EMobility Ready-to-go




Instructions manual

## Layout diagram leXsolar-EMobility Ready-to-go

Item-No.1803

### Bestückungsplan leXsolar-EMobility Ready-to-go

Art.-Nr.1803



- 1** 1800-01 Resistor module (triple) Pro  
1800-01 Widerstandsmodul 3-fach Pro
- 2** 1800-03 Resistor plug element 1 Ohm  
1800-03 Widerstands-Steckelement 1 Ohm
- 3** 2x1800-05 Resistor plug element 10 Ohm  
2x1800-05 Widerstands-Steckelement 10 Ohm
- 4** 1800-04 Resistor plug element 100 Ohm  
1800-04 Widerstands-Steckelement 100 Ohm
- 5** 1100-62 Potentiometer module  
1100-62 Potentiometermodul

Version number  
Versionsnummer

L3-03-167\_02.02.2017

- 6** 1800-15 Distilled water  
1800-15 Destilliertes Wasser
- 7** 1800-13 Lead (Pb)-battery module Pro  
1800-13 Blei-Akkumodul Pro
- 8** L2-04-102 NiZn-battery AAA  
L2-04-102 NiZn-Akku AAA
- 9** 1801-06 LiFePo-battery AAA  
1801-06 LiFePo-Akku AAA
- 10** L2-06-067 Reversible Fuel cell Pro  
L2-06-067 Reversible Brennstoffzelle Pro
- 11** Universal-power supply with **17**  
Stromversorgungsgerät mit **17**
- 12** 2xL2-06-012/013 Test leads black/red  
2xL2-06-012/013 Messleitung schw./rot
- 13** 1800-08 Battery module holder 1xAAA Pro  
with L2-04-021 NiMH battery AAA  
1800-08 Akkualterungsmodul 1xAAA Pro  
mit L2-04-021 NiMH-Akku AAA

- 14** 1118-11 Capacitor modul Pro  
1118-11 Kondensatormodul Pro
- 15** 1118-09 NiMH Battery module 3xAAA Pro  
1118-09 Akkumodul NiMH 3xAAA Pro
- 16** 1800-07 Lithium-polymer-battery module  
1800-07 Lithium-Polymer-Akkumodul Pro
- 17** 9100-13 ChargerModule  
9100-13 ChargerModul
- 18** 9100-03 AV-Module  
9100-03 AV-Modul
- 19** 1801-07 leXsolar Base unit EMobility  
1801-07 leXsolar Grundeinheit EMobility
- 20** 1801-02 Electric model car  
1801-02 Elektro-Modellfahrzeug

# leXsolar-EMobility Ready-to-go

## Instructions manual

### Content

Components.....	5
1.1 Relationship between current, resistance and voltage .....	13
1.2 Series connection of ohmic resistances .....	14
1.3 Parallel connection of ohmic resistances.....	16
2.1 Nominal voltage and capacity of voltage sources.....	17
2.2 Internal resistance of voltage sources .....	20
2.3 Series connection of voltage sources .....	23
2.5 The capacitance of a battery module.....	25
2.5 The energy density of battery modules.....	28
2.6 The $R_i$ efficiency of a battery module.....	32
2.7 The total efficiency of a battery module .....	35
3.1 The charging process of a capacitor.....	39
3.2 The discharge process of a capacitor .....	43
4.1 I-V characteristics of the single NiMH battery module .....	46
4.2 I-V characteristics of the NiZn battery module .....	49
4.3 I-V characteristics of the LiFePo battery module .....	52
4.4 I-V characteristics of the lead battery module.....	55
4.5 I-V characteristics of the lithium-polymer battery module .....	58
4.6 I-V characteristics of the triple NiMH battery module.....	61
5.1 The charging process of the NiMH battery .....	64
5.2 The charging process of the NiZn battery.....	67
5.3 The charging process of the LiFePo battery .....	69
5.4 The charging process of the lead battery.....	72

5.5 The charging process of the lithium-polymer battery .....	75
5.6 The discharging process of a battery module .....	78
6.1 Hydrogen production in the reversible hydrogen fuel cell .....	80
6.2 Characteristic curve of the electrolyzer.....	82
6.3 Hydrogen consumption of a fuel cell.....	84
6.4 Characteristic curve of the fuel cell.....	86
6.5 The efficiency of the hydrogen fuel cell.....	89
7.1 Operation of the electric car with several battery modules .....	91
7.2 Operation of the electric car with the reversible fuel cell.....	94




## I General Information

### Components



The following part contains information about the components of the experimental system. There is sketched a photograph and a small pictogram how the modules are displayed in the experimental setup. Furthermore you get information about the handling and the specifications of the components.

#### Base unit EMobility 1801-07



The base unit is a breadboard where up to 3 components can be plugged in a series and parallel connection. The current flows along the wires on the bottom side. To connect the components on the base unit with other components, there are 4 terminals at the lower end.

#### ChargerModule 9100-13



The ChargerModule is a universal battery charger for all batteries, the capacitor and the reversible fuel cell included in leXsolar-EMobility Ready-to-go. With the additional fixed-voltage outputs constant voltage of 3V or 6V can be applied. To operate the ChargerModule first the power adapter must be plugged in and connected to the input jack on the top right of the module. The charging program is selected by the "+" and "-" button and is displayed by the LEDs. The Power button is used to switch on the ChargerModule. During the charging process, the Power Enable LED flashes once per second and all keys are locked. Pressing the Power Enable button for 0.5s cancels the selected program. When the charging process is complete, there occurs an acoustic signal (3 loud "medium high" beeps, a total of about 2 seconds) and the Power Enable LED is continuously lit.

The ChargerModule provides a constant voltage (cv-mode) or constant current (cc-mode) depending on the charge program. For most battery modules a combined cc/cv-mode is applied. The top LEDs (CC/CV) indicate the applied charging mode.



For open-circuit (for example no battery module is connected to the charger) five high beeps occur and the charging program is terminated immediately. If the voltage of the connected battery module is higher than the maximum charging voltage (for example, if an incorrect battery is connected) or below the specified end-of-discharge voltage the charging program is also terminated. Independent of the connected module the charger switches off after 1 hour to prevent accidental overloading of the battery module. The following charging programs can be selected:

**NiMH (AAA):**

- Only cc-mode (charge current  $I = 250 \text{ mA}$ ) without cv-process
- Upper voltage limit: 1.6V
- Lower voltage limit: 1V

**NiZn (AAA):**

- Starts with cc-mode ( $I=250\text{mA}$ ) up to a switching voltage  $V=1.8\text{V}$
- After reaching the threshold voltage switch to cv-mode, switch-off at a current of 100mA
- Upper voltage limit: 2V
- Lower voltage limit: 1.3V

**LiFePo (AAA):**

- Starts with cc-mode ( $I=200\text{mA}$ ) up to a switching voltage  $V=3.6\text{V}$
- After reaching the threshold voltage switch to cv-mode, switch-off at a current of 100mA
- Upper voltage limit: 3.7V
- Lower voltage limit: 2.8V

**NiMH (triple):**

- Only cc-mode (charge current  $I = 250 \text{ mA}$ ) without cv-process
- Upper voltage limit: 4.8V
- Lower voltage limit: 3V

**Pb:**

- Starts with cc-mode ( $I=500\text{mA}$ ) up to a switching voltage  $V=2.35\text{V}$
- After reaching the threshold voltage switch to cv-mode, switch-off at a current of 200mA
- Upper voltage limit: 2.45V
- Lower voltage limit: 1.8V

**LiPo:**

- Starts with cc-mode ( $I=500\text{mA}$ ) up to a switching voltage  $V=4.1\text{V}$
- After reaching the threshold voltage switch to cv-mode, switch-off at a current of 200mA
- Upper voltage limit: 4.3V
- Lower voltage limit: 3V

**Electrolyzer:**

- Only cv-mode ( $V=2.1\text{V}$ )

**Supercap:**

- Only cv-mode ( $V=2.1\text{V}$ ), switch-off at a current of 50mA
- Upper current limit: 2A
- Switch-off after 10min, independent of current


**3V:**

- Constant voltage of 3V

**6V:**

- Constant voltage of 6V

## AV-Module 9100-03



The AV-Module is a combined voltage and current meter. It holds 3 buttons, whose features are described in the display respectively. By pushing a random button the module will switch on. In the disabled state the display shows the leXsolar emblem. When the display does not show anything or the word „Bat“ is shown, it is necessary to change the batteries in the back (2 x AA batteries 1.2 to 1.5V; Take care of the polarity marked on the bottom of the battery case! Do not touch the button while inserting the batteries).

With the top right button the measuring mode can be switched between voltage mode, current mode or combined voltage-current mode. Both measurement mode and required cable connection will be indicated by the circuit symbols on the display. Take care that in voltage mode no current is applied to the right jack. In the combined mode the voltage can be measured with the right jack as well as with the left one. The influence of the internal resistance of the current measurement is compensated internally. The measured values are signed. When the positive pole is connected to a red jack and the negative pole is connected to the black jack, the value of the voltage will be positive. When current is applied from the left to the right, the current value will be positive, as well. The other way around, the algebraic sign changes.

After 30 min without pushing a button or after 10 min of measuring a constant value, the module will switch off automatically. It can measure voltages up to 12 V and currents up to 2 A. In case of exceeding one of the values, the module interrupts the current flow and shows “overcurrent” or “overvoltage”. This error message can be confirmed by touching a button. The module will resumes measuring, when the values attain acceptable values.

### Specifications:

#### Voltage metering:


- range: 0...12 V
- accuracy: 1 mV
- automatic shutoff in case of overvoltage >12 V

#### Current metering:

- range: 0...2 A
- accuracy: 0,1 mA (0...199 mA) and 1mA (200 mA...1 A)
- automatic shutoff in case of overcurrent >2 A
- internal resistance <0,5 Ohm (0...200 mA); <0,2 Ohm (200 mA...2 A)





### Resistor plug module, triple (1800-01 ) with resistor plug elements



With the resistor plug module and the belonging resistor plug elements parallel connection and series connection of resistors are possible. For parallel connection use one resistor module (triple) with three slots. For series connection use two triple resistor modules. The following resistor plug elements are included:


1 x $R=1\Omega$	1800-03
2 x $R=10\Omega$	1800-05
1 x $R=100\Omega$	1800-04

### Potentiometer module 1100hm Pro 1100-62



The potentiometer module holds a 0-10- $\Omega$ -potentiometer and a 0-100- $\Omega$ -potentiometer. Both are serially connected, so that the potentiometer can attain resistances between 0  $\Omega$  to 110  $\Omega$ . The measuring error amounts to 0.5  $\Omega$  for the small resistor and 3  $\Omega$  at other one. The maximum current amounts to 1A.

### Capacitor module 1118-11




The capacitor module consists of 2 series-connected capacitors. Charging voltages for the capacitor should not exceed 5 V. It is possible to short-circuit the capacitor to discharge, because there are fuses to avoid damages.

#### Specifications:

Capacitance: 5 F  
Maximum voltage: 5,4 V



### NiMH-battery module, single L2-04-021 with mount 1800-08




#### Specifications:

$V=1,0V \dots 1,35V$

End-point voltage: 1V

Max. charging voltage: 1,6V

### NiMH-battery module, triple 1118-09




#### Specifications:

$V=3V \dots 4,05V$

End-point voltage: 3V

Max. charging voltage: 4,8V

### LiFePo-battery module 1801-06




#### Specifications:

$V=3,2V \dots 3,4V$

End-point voltage: 2,8V

Max. charging voltage: 3,6V

### NiZn-battery module 1801-06




#### Specifications:

$V=1,3V \dots 1,8V$

End-point voltage: 1,3V

Max. charging voltage: 1,9V

### Lead-battery module 1800-13




#### Specifications:

$V=1,9V \dots 2,15V$

End-point voltage: 1,9V

Max. charging voltage: 2,35V

### Lithium-Polymer-battery module 1800-07




#### Specifications:

$V=3V \dots 4,2V$


End-point voltage: 3V

Max. charging voltage: 4,2V

### Expansion *(not included in EMobility Ready-to-go)*: Battery adapter cable 1800-09



All battery modules are equipped with an additional connection for the four-point measurement. The adapter cable is connected with the black connector to this port:




To measure the voltage, the red and the black cable are connected with the measurement device. For measuring the resistance the white cable instead of the red is used.

### Reversible fuel cell (L2-06-067)




The reversible fuel cell consists of an electrolyzer and a fuel cell. To fill the reversible fuel cell you should proceed in the following way:

1. Fill the rev. fuel cell with distilled water as shown in the alongside figure.
2. Fill both storage cylinders up to the top of the tubules, which are inside the cylinders.
3. Knock the rev. fuel cell slightly on the table.
4. Continue filling in water until it flows through the tubules.
5. Close the storage cylinders with the plugs and turn over the rev. fuel cell (the plugs must be on the bottom).



**!! Advice:** To charge the reversible fuel cell the applied voltage should not exceed 1.5 V. Otherwise the resulting current could exceed 1 A, which would damage the fuel cell.

### Electric model car with battery adapter 1801-02



The electric model car can be used with the reversible fuel cell or the battery modules. The fuel cell can be plugged directly onto the car. The battery modules can be plugged with the adapter onto the car.

The car will move when both cables are connected with the voltage source. There will be a short circuit when the wires are held during the short circuit.




## 1.1 Relationship between current, resistance and voltage

### Task

Examine the relationship between voltage, current and resistance in a simple electrical circuit.

### Setup



### Equipment required

- base unit
- 1 ChargerModule
- 1 Potentiometer module
- 1 AV-Module
- 1 cable

### Procedure

1. Set up the experiment according to the circuit diagram. The Charger module is plugged into the base unit rotated by  $90^\circ$  (see sketch). Use the ChargerModule with constant voltage mode at 6V. For handling instructions see page 5.
2. Measure current  $I$  and voltage  $V$  for different resistances  $R$  at the potentiometer (for values see table). Use the AV-Module in voltage-current-mode. For handling instructions see page 7.
3. Note your measured values in the table.

### Data

$R (\Omega)$	100	80	60	40	20
$V (V)$					
$I (mA)$					
$V/I (\Omega)$					

### Evaluation

1. Calculate for each measuring point the ratio  $V/I$  and note your values in the table.
2. Which law reflects your findings? Illustrate these principles using data from the table.

---



---



---




## 1.2 Series connection of ohmic resistances

## Task

Examine the series connection of ohmic resistances.

## Setup



## Required devices

- base unit
- 1 ChargerModule
- 1 Potentiometer module
- 1 resistor module, triple
- 1 resistor plug element ( $R=100\Omega$ )
- 1 AV-Module
- 4 cables

*Additionally needed:*

- 1 voltage measurement device

## Execution

1. Set up the experiment according to the circuit diagram. Use the ChargerModule with constant voltage mode at 6V. For handling instructions see page 5.
2. Adjust the resistance  $R$  of the potentiometer to a value of  $R_{Pot}=100\Omega$  and use the resistor plug element of  $R_S=100\Omega$  at the triple resistor module.
3. Measure each voltage  $V$  and current  $I$  over both resistances ( $V_{tot}$ ) and the single voltage ( $V_1, V_2$ ).

**Note:** The AV-Module is plugged into the base unit rotated by  $90^\circ$  (see sketch). It is used in current-mode. If no further measurement device is available, you can use the AV-Module in voltage mode to measure the voltage. For this purpose the slot of the AV module should be electrically bridged using a cable. For handling instructions see page 7.

4. Repeat your measurement for further resistance values at the potentiometer (see table).
5. Note your measured data in the table.

## Measurements

$R_S (\Omega)$	100	80	60	40	20
$V_1 (V)$					
$V_2 (V)$					
$V_{tot} (V)$					
$I (mA)$					
$R_{ges} = V_{tot} / I (\Omega)$					



## 1.2 Series connection of ohmic resistances

## Evaluation

1. Calculate each the ratio  $R_{\text{tot}} = V_{\text{tot}}/I$  and note your values in the table above.
2. Calculate each the sum of the single voltages ( $V_1 + V_2$ ) and compare it the voltage over both resistances ( $V_{\text{tot}}$ ).
3. What is the influence of the resistance on the current  $I$  and the voltages  $V_1 + V_2$ , respectively  $V_{\text{tot}}$ ?
4. What is the connection between the total resistance  $R_{\text{tot}}$  and the single resistances? Formulate a law for the calculation of the total resistance in a series connection of resistances.

2.

	$V_1 + V_2$	$V_{\text{tot}}$
$R_{\text{Pot}} = 100\Omega / R_S = 100\Omega$		
$R_{\text{Pot}} = 80\Omega / R_S = 100\Omega$		
$R_{\text{Pot}} = 60\Omega / R_S = 100\Omega$		
$R_{\text{Pot}} = 40\Omega / R_S = 100\Omega$		

3.

---

---

---

---

---

---

---

---

4.

---

---

---

---

---

---

---

---




## 1.3 Parallel connection of ohmic resistances

## Task

Examine the parallel connection of ohmic resistances.

## Setup



## Required devices

- base unit
- 1 ChargerModule
- 1 resistor module, triple
- 3 resistor plug elements  
(2x  $R=10\Omega$ , 1x  $R=100\Omega$ )
- 1 AV-Module
- 1 cable

## Execution

1. Set up the experiment according to the circuit diagram. The Charger module is plugged into the base unit rotated by  $90^\circ$  (see sketch). Use the ChargerModule with constant voltage mode at 3V. For handling instructions see page 5.
2. Start with 1 x  $10\Omega$  resistance. Measure the voltage and current  $I$ . Use the AV-Module in current-voltage mode. For handling instructions see page 7.
3. Repeat your measurement for the parallel connection of the following resistances and note your measured data in the table:

- $R_1=10\Omega$  /  $R_2=10\Omega$
- $R_1=10\Omega$  /  $R_2=100\Omega$
- $R_1=10\Omega$  /  $R_2=10\Omega$  /  $R_3=100\Omega$

## Measurements

	$R_1=10\Omega$	$R_1=10\Omega$ / $R_2=10\Omega$	$R_1=10\Omega$ / $R_2=100\Omega$	$R_1=10\Omega$ / $R_2=10\Omega$ / $R_3=100\Omega$
V (V)				
I (mA)				
$R_{ges}=V/I$ ( $\Omega$ )				

## Evaluation

1. What is the influence of the resistance on the current  $I$  and the voltage  $V$ ?
2. Formulate a law for the calculation of the total resistance in a parallel connection of resistances.

---



---



---



---




## 2.1 Nominal voltage and capacity of voltage sources

### Task

Determine the open-circuit voltage and the capacity of single cells.

### Setup



### Equipment required

- base unit
- 1 AV-Module
- 1 battery module NiMH, single
- 1 battery module Pb
- 1 battery module LiPo
- 1 battery module NiZn
- 1 battery module LiFePo

### Procedure

1. Set up the experiment according to the circuit diagram. Use the provided battery modules.
2. Measure the respective open-circuit voltages of the voltage sources  $V_0$  and write down your measured values in the table below. Use the AV-Module in voltage mode. For handling instructions see page 7.

### Evaluation

1. Use the open-circuit diagram below to determine the charge state of the voltage sources and note your values (in percent) in the table.
2. Calculate the remaining capacity of each battery using the determined charge states and the indicated maximum capacity. Use the following formula:

$$\frac{\text{Remaining capacity}}{\text{Maximum capacity}} = \frac{\text{Charge State in \%}}{100}$$


3. Calculate the required battery capacity to operate a radio with a power of 20W for a period of 3h at a battery voltage of 12V.
4. A starter battery was loaded for 5h with a capacity of 40Ah. Calculate the discharge current.

### Data

battery module	$V_0$ in V	Charge state in %
NiMH		
NiZn		
LiFePo		
Pb		
LiPo		

## 2.1 Nominal voltage and capacity of voltage sources

## Evaluation



*Diagram 2.1: Determination of the charge state*

battery module	Capacity
NiMH	600mAh
NiZn	550mAh
LifePo	200mAh
Pb	2500mAh
LiPo	980mAh

*Table 2.1: Maximum capacity of the battery modules*

2.

[illegible]

## Evaluation

[illegible]

A full-page grid of small squares, typical of graph paper. The grid consists of 20 columns and 15 rows of squares, creating a uniform pattern across the entire page.


[illegible]

## 2.2 Internal resistance of voltage sources

### Task

Determine the internal resistance of the provided voltage sources.

### Setup



### Equipment required

- base plate
- all battery modules
- 1 AV-Module
- 1 triple resistor module
- 1 resistor plug-element (10Ω)
- cables

### Procedure

1. Set up the experiment according to the circuit diagram.
2. Measure first of all the open-circuit voltage of the voltage sources  $V_0$  without a resistor and enter your data in the table.
3. Now measure the voltage  $V_{load}$  and the current  $I_{load}$  by closing the circuit (plug in the resistor).

### Evaluation

1. Determine the internal resistance  $R_i$  of the voltage sources and enter this value into the table too. The internal resistance of the cells is given by the following formula:

$$R_i = \frac{V_0 - V_{load}}{I_{load}} - R_C$$

**Advice:** The contact resistances  $R_C$  have to be subtracted from the calculated value in order to get the actual internal resistances. The contact resistances are also given in the data table.

2. In what way are the voltage sources different from each other?
3. What percentage of the power consumes the triple NiMH-battery module for its own heating in this example?
4. A starter battery has a voltage  $U_0 = 12V$  and an internal resistance  $R_i = 20m\Omega$ . An external starter of  $60m\Omega$  is then connected.
  - a) Which current is flowing during start?
  - b) Calculate the voltage drop at the clamps during start.

## 2.2 Internal resistance of voltage sources

## Data

	$V_0$ in V	$V_{load}$ in V	$I_{load}$ in mA	$R_i$ in m $\Omega$	$R_c$ in m $\Omega$
NiMH-battery module, single ( $R_{load} = 5 \Omega$ )					50
NiZn-battery module ( $R_{load} = 5 \Omega$ )					50
LiFePo-battery module ( $R_{load} = 5 \Omega$ )					50
Lead-battery module ( $R_{load} = 5 \Omega$ )					5
Lithium-polymer-battery module ( $R_{load} = 10 \Omega$ )					50
NiMH-battery module, triple ( $R_{load} = 10 \Omega$ )					220

## Evaluation

1.

This image shows a full page of blank graph paper. The grid consists of small, equal-sized squares formed by thin gray lines. There are 20 columns and 20 rows of squares, creating a total area of 400 small squares. The grid covers the entire page except for a narrow white margin around the edges.

2.

---

## Evaluation

1.

[illegible]

2.

[illegible]




## 2.3 Series connection of voltage sources

## Task

Investigate the behavior of voltage sources in a single cell setup and series connection respectively!

## Setup

## Equipment required



- base plate
- 1 battery module NiMH, single
- 1 battery module NiMH, triple
- 1 resistor module, triple
- 2 resistor plug elements (2x10Ω)
- 1 AV-Module
- cables

## Procedure

1. Set up the experiment according to the circuit diagram. Use the single NiMH-battery module first. Set up a parallel connection with both resistor plug elements to achieve a load resistance of 5Ω. Do not plug in the triple resistor module yet to avoid a discharge of the battery module.
2. First measure the open-circuit voltage  $V_0$  without the resistor and note your value in the table.
3. Plug in the triple resistor module and measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$ . Use the AV-module in current-voltage-mode.
4. Repeat the experiment with two further NiMH single cells.
5. Now use the examined cells in the adapter for the triple NiMH-module and measure again the above values.

## Data

## Single battery module:

	$V_0$ in V	$V_{\text{Load}}$ in V	$I_{\text{Load}}$ in mA
1st cell ( $R_C = 50 \text{ m}\Omega$ )			
2nd cell ( $R_C = 50 \text{ m}\Omega$ )			
3rd cell ( $R_C = 50 \text{ m}\Omega$ )			

## Triple battery module:

	$V_0$ in V	$V_{\text{Load}}$ in V	$I_{\text{Load}}$ in mA
Triple cell ( $R_C = 220 \text{ m}\Omega$ )			

## Evaluation

- $$R_i = \frac{V_0 - V_{load}}{I_{load}} - R_C$$

2. Why is it better to use a single cell with the same voltage as a triple cell rather than using a comparable battery with several cells connected in series?

A full-page sheet of white graph paper featuring a light gray grid. The grid consists of small, equal-sized squares arranged in a continuous pattern across the entire page. There are no margins, text, or other markings present.

---

---

---

---

---




## 2.5 The capacitance of a battery module

## Task

Determine the capacitance of a battery module.

## Setup



## Equipment required

- base plate
- 1 LiFePo battery module
- 1 battery adapter cable
- 1 resistor plug module, triple
- 2 resistor plug elements (2x10Ω)
- 1 AV-Module
- cables

*Additionally needed (optional):*  
- PC with data analysis software

## Execution

1. Set up the experiment according to the circuit diagram. Set up a parallel connection with both resistor plug elements to achieve a load resistance of 5Ω. Do not plug in the triple resistor module yet to avoid a discharge of the battery module.
2. First measure the open-circuit voltage  $V_0(1)$  without the resistor and note your value in the table.
3. Plug in the triple resistor module and measure 15min the voltage  $V_{Load}$  and current  $I_{Load}$  at intervals of 1min. Use the AV-module in current-voltage-mode.
4. Measure the open-circuit voltage  $V_0(2)$  five minutes after the experiment.

**Advice:** The battery module should have a rest capacity of 50% (corresponds to  $V_0=3.3V$ ). The experiment has to be interrupted as soon as the discharge current drops significantly.

## Measurements

$V_0(1) = \underline{\hspace{2cm}}$

t in min	0	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
$V_{Load}$ in V																
$I_{Load}$ in mA																

$V_0(2) = \underline{\hspace{2cm}}$



## 2.4 The capacitance of a battery module


### Evaluation

1. Enter your values in the given diagram.
2. Determine the charge state and capacity of the battery module before and after the experiment. Use the diagram and table from the experiment “Nominal voltage and capacity of voltage sources”. Estimate from your values the loss of capacity during the experiment.
3. Explain why some parts of the discharge curve can lead to problems with the charge level indication of LiFePo batteries.
4. Transfer your experimentally determined values in a data analysis software. Define with the help of the software a polynomial curve which approximately describes the course of the I-t curve. Then determine the dissipated charge from the integral of the I-t-curve:

$$Q = \int_{t_1}^{t_2} I dt$$

3. Compare the estimated charge Q to the loss of capacity that you determined in task 2.

### Diagram



## Evaluation

2.

This image shows a full page of blank graph paper. The grid consists of small, equal-sized squares formed by thin gray lines. There are 20 columns and 20 rows of squares, creating a total area of 400 small squares. The grid is perfectly aligned and covers the entire page without any margins or additional markings.

3.

---

4.

---



---




## 2.5 The energy density of battery modules

### Task

Determine the energy density of several battery modules.

### Setup



### Equipment required

- base plate
- 1 AV-Module
- All battery modules

### Execution

1. Set up the experiment according to the circuit diagram.
2. Measure each open-circuit voltage of the voltage sources  $V_0$  and note your values in the table. Use the AV-module in voltage mode.

### Evaluation

1. Use the illustration below (diagram 2.6) to determine the charge state of the voltage sources and note the respective percentage values in the table.
2. Calculate the remaining capacity  $Q_R$  with the following formula (maximum capacity is stated in table 2.6):

$$\frac{\text{remaining capacity } Q_R}{\text{maximum capacity } Q_{max}} = \frac{\text{charge state in \%}}{100}$$

3. Calculate the energy content of the various battery modules and enter your values into the table. The energy content can be calculated by the following formula (note units!):

$$E = V \cdot I \cdot t = V_0 \cdot Q_R$$

4. Calculate using the weight of the batteries (given in table 2.6) the mass-based energy  $\omega$  [kJ/kg].
5. Why batteries with relatively low (mass-based) energy density like the lead-acid battery are widely used in various applications, despite their heavy weight?



## 2.5 The energy density of battery modules

## Measurements

battery module	$V_0$ in V	charge state in %	$Q_R$ in mAh	E in kJ
NiMH				
NiZn				
LiFePo				
Pb				
LiPo				

## Evaluation




Diagram 2.6: Determination of charge state

battery module	maximum capacitance $Q_{\max}$	weight m in g
NiMH	600mAh	11.3
NiZn	550mAh	11.2
LiFePo	200mAh	7.8
Pb	2500mAh	177.4
LiPo	980mAh	20.0


Table 2.6: Maximum capacitance and weight of the battery modules




## 2.5 The energy density of battery modules

### Evaluation

2.



3.








## 2.6 The $R_i$ efficiency of a battery module

### Task

Determine the  $R_i$  efficiency of a battery module.

### Setup



### Equipment required

- base plate
- 1 battery module
- 1 potentiometer module
- 1 AV-Module
- cables

### Execution

1. Set up the experiment according to the circuit diagram.
2. Measure first of all the open-circuit voltage  $V_0$  of the voltage sources without the potentiometer module and enter your data in the table.
3. Now plug in the resistor module and measure the voltage  $V_{\text{load}}$  and the current  $I_{\text{load}}$  for different resistances  $R_{\text{pot}}$  at the potentiometer. Use the AV-module in current-voltage-mode. Note your values in the table.

**Advice:** Interrupt the current flow (for example, by removing the cable) after each individual measurement to prevent excessive discharge of the module during the experiment

### Evaluation

1. Determine for each resistance  $R_{\text{pot}}$  at the potentiometer module the  $R_i$  efficiency  $\eta$  of the battery module and enter your values into the table. The  $R_i$  efficiency  $\eta$  can be found using the following formula:

$$\eta = \frac{P_{\text{load}}}{P_0} = \frac{V_{\text{load}} \cdot I_{\text{load}}}{V_0 \cdot I_{\text{load}}}$$

2. Enter your values into the diagram.
3. Describe and explain the behavior of the  $R_i$  efficiency in dependence of the current  $I_{\text{load}}$ .


## 2.6 The $R_i$ efficiency of a battery module

## Measurement

$$V_0 = \underline{\hspace{2cm}}$$

R <sub>Pot</sub> in Ω									
V <sub>Load</sub> in V									
I <sub>Load</sub> in mA									
η in %									

## Diagram



## Evaluation

1.

A blank sheet of graph paper with a grid pattern. The grid consists of small squares formed by thin gray lines. There are 20 columns and 15 rows of squares. A single horizontal line runs across the middle of the page, separating the top 7 rows from the bottom 8 rows. This line is slightly thicker than the other grid lines. The entire grid is enclosed within a rectangular border made of the same thin gray lines.

## Evaluation

This image shows a full page of blank graph paper. The grid consists of small, equal-sized squares formed by thin gray lines. There are no margins, text, or other markings on the page.

3.

---

---

---

---


## 2.7 The total efficiency of a battery module

### Task

Determine the total efficiency of a battery module.

### Setup


#### Part 1: Charging process



### Equipment required

- base plate
- 1 battery module NiMH, single
- 1 ChargerModule
- 1 AV-Module
- 1 potentiometer module
- cables

#### Part 2: Discharging process



### Execution

#### Part 1: Charging process

1. Set up the experiment according to the circuit diagram (Part 1). Use the ChargerModule in NiMH-mode (single). For advices about the handling of the ChargerModule, see page 5. Do not start the charger yet.
2. Measure first of all the open-circuit voltage  $V_0(1)$  of the battery module and enter your data in the table.
3. Switch on the charger and measure for ten minutes in intervals of one minute, the voltage  $V$  and the current  $I$ . Use the AV-Module in current-voltage-mode. Enter your values in the table.
4. Measure 5 minutes after completion of the first part of the experiment again the open circuit voltage  $V_0(2)$ .

**Advice:** The charge state of battery module should be maximal 50% before the beginning of the experiment (corresponds to an open circuit voltage of 1,18V). Optionally, the battery module must be discharged before the experiment with the help of the resistor modules.



## 2.7 The total efficiency of a battery module

## Execution

## Part 2: Discharging process

1. Set up the experiment according to the circuit diagram (Part 2). Do not plug in the potentiometer module yet to avoid the beginning of the measurement without recording the measured data.
2. Measure first of all the open-circuit voltage  $V_0(1)$  of the battery module and enter your data in the table.
3. Plug in the potentiometer module and regulate the value of the discharge current to the value of charge current from part 1 of the experiment.
4. Measure then ten minutes the voltage  $V$  and the current  $I$  in intervals of one minute. Readjust if necessary the resistance of the potentiometer to keep the discharge current constant. Enter your values in the table.
5. Stop after ten minutes the current flow. Measure five minutes the open-circuit voltage  $V_0(2)$  of the battery module immediately after completion of the experiment in intervals of one minute. Enter your values in the table.

## Measurements

## Part 1: Charging process

 $V_0(1) = \underline{\hspace{2cm}}$ 

t in min	0	1	2	3	4	5	6	7	8	9	10
V in V											
I in mA											
W in J											

 $V_0(2) = \underline{\hspace{2cm}}$ 

## Part 2: Discharging process

 $V_0(1) = \underline{\hspace{2cm}}$ 

t in min	0	1	2	3	4	5	6	7	8	9	10
V in V											
I in mA											
W in J											

Open-circuit voltage  $V_0(2)$  after experiment:

t in min	0	1	2	3	4	5
$V_0(2)$ in V						

## 2.7 The total efficiency of a battery module

## Evaluation

1. Calculate each the energy  $W$  that was expended, respectively consumed during the charging/discharging process and enter your values in the table.


$$W = V \cdot I \cdot t$$

2. Enter your values in the given diagram.
3. Determine the electric energy  $W_1$ , that was required during the first part of the experiment for charging the battery ( $t=10\text{min}$ ). Determine further the electric energy  $W_2$ , that was submitted during the second part of the experiment ( $t=10\text{min}$ ). Calculate the total efficiency  $\eta$  of the battery module.

$$\eta = \frac{W_2}{W_1}$$

4. Describe what mainly affects the efficiency of battery modules.
5. Explain the so-called *Peukert-Effect*.

## Diagram



## Evaluation

- 1.

[illegible]

## Evaluation

[illegible]

3.

[illegible]

4.

---

---

---

---

---

---

5.

---

---

---


---

## 3.1 The charging process of a capacitor

### Task

Record the charging curve of a capacitor.

### Setup



### Equipment required

- base plate
- 1 ChargerModule
- 1 AV-Module
- 1 capacitor module
- 1 resistor module, triple
- 2 resistor plug elements ( $2 \times R = 10 \, \Omega$ )
- cables

### Execution

#### Part 1: Supercap-mode

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in Supercap-mode. For handling instructions see page 5.
2. Plug in the resistor module of  $10 \, \Omega$  and measure 90s the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  in intervals of 10s. The AV-module is in current-voltage-mode.
3. Enter all your values in the table.
4. Discharge the capacitor and repeat the experiment for a resistance of  $5 \, \Omega$  (parallel connection of  $2 \times 10 \, \Omega$ ).

#### Part 2: Constant-voltage-mode

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in Constant-voltage-mode with 3V.
2. Record analogously to part 1 the charging curve for different resistances. Measure now 120s voltage  $V_{\text{Load}}$  and Current  $I_{\text{Load}}$  in intervals of 10s and enter your values in the table.

**Advice:** The resistor module should be plugged off before the measurement to avoid starting the experiment without recording the data.

### Evaluation

1. Enter your values in the diagram.
2. Describe the charging behavior of the capacitor.
3. Estimate the period of time after which the capacitor has a charge state of 60% (In constant-voltage-mode, where  $3V \equiv 100\%$ ).
4. Name scopes of application for supercaps.



## 3.1 The charging process of a capacitor

## Measurements: Part 1

 **$R_1 = 10\Omega$ :**

t in s	0	10	20	30	40	50	60	70	80	90
$V_{\text{Load}}$ in V										
$I_{\text{Load}}$ in mA										

 **$R_1 = 5\Omega$ :**

t in s	0	10	20	30	40	50	60	70	80	90
$V_{\text{Load}}$ in V										
$I_{\text{Load}}$ in mA										

## Measurements: Part 2

 **$R_1 = 10\Omega$ :**

t in s	0	10	20	30	40	50	60	70	80	90	100	110	120
$V_{\text{Load}}$ in V													
$I_{\text{Load}}$ in mA													

 **$R_1 = 5\Omega$ :**



t in s	0	10	20	30	40	50	60	70	80	90	100	110	120
$V_{\text{Load}}$ in V													
$I_{\text{Load}}$ in mA													




## 3.1 The charging process of a capacitor

### Diagrams

#### Part 1: Supercap-mode




#### Part 2: Constant-voltage-mode





### 3.1 The charging process of a capacitor

#### Diagrams



#### Evaluation

2.

---

---

---

---

---

3.

$R=10\ \Omega$ :

$R=5\ \Omega$ :

4.

---

---

---

---

---

---

---


---

## 3.2 The discharge process of a capacitor

### Task

Record the discharge curve of the capacitor.

### Setup



### Equipment required

- base plate
- 1 AV-Module
- 1 resistor module, triple
- 2 plug-in resistor elements (2 x  $R=10\Omega$ )
- 1 capacitor module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the resistor module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the resistor module ( $R=10\Omega$ ) and measure 90s the voltage  $V_{Load}$  and current  $I_{Load}$  in intervals of 10s. Use the AV-module in current-voltage-mode.
4. Repeat the experiment for a resistance of  $5\Omega$  (parallel connection of 2 x  $10\Omega$ )

**Advice:** The capacitor should have the same charge state before starting both parts of the experiment. Hence, charge the capacitor after the first part. For handling instructions see experiment "The charging process of a capacitor"

### Evaluation

1. Enter your values in the given charts.
2. What can you conclude from the diagrams on the discharge behavior of the capacitor?
3. Estimate the time after which the capacitor has a charge state of 60%.
4. Calculate the charge state of the capacitor at the beginning and after 90s discharge time for a resistance of  $R=10\Omega$  (capacity  $C=5,0F$ ).
5. The capacitance of a capacitor is given as "n47". Which capacity corresponds to this specification?



## 3.2 The discharge process of a capacitor

## Measurements



 **$R_1 = 10\Omega$ :**

t in s	0	10	20	30	40	50	60	70	80	90
V in V										
I in mA										

 **$R_2 = 5\Omega$ :**

t in s	0	10	20	30	40	50	60	70	80	90
V in V										
I in mA										

## Diagrams



### 3.2 The discharge process of a capacitor

## Evaluation

2.

---

---

---

---

---

---

3.

[illegible]

4.

5.




## 4.1 I-V characteristics of the single NiMH battery module

### Task

Obtain the I-V characteristics of a single NiMH battery module.

### Setup



### Equipment required

- base plate
- 1 NiMH battery module, single
- 1 AV-Module
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the potentiometer module and adjust the resistance to  $50\Omega$ .
4. Measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  and note your values in the table. Use the AV-module in current-voltage-mode.
5. Decrease in several steps the resistance  $R_{\text{Pot}}$  at the potentiometer and measure each the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$ . Note all your values in the table.

**Advice:** Interrupt the current flow (for example by removing the cable) after each single measurement to avoid discharge of the battery module during the experiment.


### Evaluation

1. Enter your values in the given chart.
2. Compare your measured characteristics with the added characteristics and make a statement about the charge state of the cell. Calculate the remaining capacity of the battery module. You find instructions for this in experiment "Nominal voltage and capacity of voltage sources".
3. Name applications for NiMH batteries. Explain those according to their characteristics
4. Give reasons, why NiMH batteries can not be used in safety-related devices such as fire alarms or emergency flashlights.
5. What benefits have NiMH batteries over NiCd batteries?

## Measurements

$$V_0 = \underline{\hspace{2cm}}$$
[illegible]

## Diagram



## Evaluation

2.

[illegible]



## 4.1 I-V characteristics of the single NiMH battery module

### Evaluation

3.

---

---

---

---

---

---

---

---

4.

---

---

---

---

---

---

---

---

---

---

5.

---

---

---

---

---

---

---

---

---

---

---

---




## 4.2 I-V characteristics of the NiZn battery module

### Task

Determine the I-V characteristics of the NiZn battery module.

### Setup



### Equipment required

- base plate
- 1 NiZn battery module
- 1 AV-Module
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the potentiometer module and adjust the resistance to  $60\Omega$ .
4. Measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  and note your values in the table. Use the AV-module in current-voltage-mode.
5. Decrease in several steps the resistance  $R_{\text{Pot}}$  at the potentiometer and measure each the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$ . Note all your values in the table.

**Advice:** Interrupt the current flow (for example by removing the cable) after each single measurement to avoid discharge of the battery module during the experiment.

### Evaluation

1. Enter your values in the given chart.
2. Compare your measured characteristics with the added characteristics and make a statement about the charge state of the cell. Calculate the remaining capacity of the battery module. You find instructions for this in experiment "Nominal voltage and capacity of voltage sources".
3. Why was the NiZn battery practically used first in the 2000 years, although Adison has patented these battery types already in 1901?
4. Name advantages of the NiZn batteries over the NiMH systems, especially in the automotive industry.




## 4.2 I-V characteristics of the NiZn battery module

## Measurements

 $V_0 =$  \_\_\_\_\_


$R_{Pot}$ in $\Omega$									
$U_{Load}$ in V									
$I_{Load}$ in mA									

## Diagram



## Evaluation

2.





## 4.2 I-V characteristics of the NiZn battery module

### Evaluation

3.

---

---

---

---

---

---

---

---

---

---

4.

---

---

---

---

---

---

---

---

---


---

## 4.3 I-V characteristics of the LiFePo battery module

### Task

Determine the I-V characteristics of the LiFePo battery module.

### Setup



### Equipment required

- base plate
- 1 LiFePo battery module
- 1 AV-Modul
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the potentiometer module and adjust the resistance to  $100\Omega$ .
4. Measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  and note your values in the table. Use the AV-module in current-voltage-mode.
5. Decrease in several steps the resistance  $R_{\text{Pot}}$  at the potentiometer and measure each the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$ . Note all your values in the table.

**Advice:** Interrupt the current flow (for example by removing the cable) after each single measurement to avoid discharge of the battery module during the experiment.

### Evaluation

1. Enter your values in the given chart.
2. Compare your measured characteristics with the added characteristics and make a statement about the charge state of the cell. Calculate the remaining capacity of the battery module. You find instructions for this in experiment "Nominal voltage and capacity of voltage sources".
3. Name advantages and disadvantages of LiFePo batteries against other battery types.
4. What are the main fields of application for LiFePo batteries?

## Measurements

[illegible]

The graph shows the relationship between  $V \ln V$  (y-axis, ranging from 2.8 to 3.5) and current  $I$  in mA (x-axis, ranging from 0 to 300). Four data series are plotted: 'Upper limit' (dashed black line at 3.5), 'Charge state 90%' (solid green line), 'Charge state 50%' (solid yellow line), and 'Charge state 10%' (solid red line). All three charge state curves show a decreasing trend as current increases. A 'Lower limit' (dashed black line at 2.8) is also indicated.

I in mA	Upper limit	Charge state 90%	Charge state 50%	Charge state 10%	Lower limit
0	3.5	3.38	3.29	3.14	2.8
50	3.5	3.32	3.24	3.09	2.8
100	3.5	3.29	3.21	3.05	2.8
150	3.5	3.26	3.18	3.02	2.8
200	3.5	3.24	3.15	2.99	2.8
250	3.5	3.22	3.12	2.97	2.8
300	3.5	3.20	3.10	2.95	2.8

[illegible]

### 4.3 I-V characteristics of the LiFePo battery module

## Evaluation


[illegible][illegible]

## 4.4 I-V characteristics of the lead battery module

### Task

Determine the I-V characteristics of the lead battery module.

### Setup



### Equipment required

- base plate
- 1 lead battery module
- 1 AV-Modul
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the potentiometer module and adjust the resistance to  $60\Omega$ .
4. Measure the voltage  $V_{Load}$  and current  $I_{Load}$  and note your values in the table. Use the AV-module in current-voltage-mode.
5. Decrease in several steps the resistance  $R_{Pot}$  at the potentiometer and measure each the voltage  $V_{Load}$  and current  $I_{Load}$ . Note all your values in the table.

**Advice:** Interrupt the current flow (for example by removing the cable) after each single measurement to avoid discharge of the battery module during the experiment.

### Evaluation


1. Enter your values in the given chart.
2. Compare your measured characteristics with the added characteristics and make a statement about the charge state of the cell. Calculate the remaining capacity of the battery module. You find instructions for this in experiment "Nominal voltage and capacity of voltage sources".
3. What is meant when a lead battery suffers from sulfation ?
4. Why are electrodes in the automotive industry often performed highly porous in lead-acid batteries?

#### 4.4 I-V characteristics of the lead battery module

## Measurements

$$V_0 = \underline{\hspace{2cm}}$$
[illegible]

## Diagram



## Evaluation

2.

[illegible]



## 4.4 I-V characteristics of the lead battery module

### Evaluation

3.

---

---

---

---

---

---

---

---

---

---

4.

---

---

---

---

---

---

---

---

---


---

## 4.5 I-V characteristics of the lithium-polymer battery module

### Task

Determine the I-V characteristics of the lithium-polymer battery module.

### Setup



### Equipment required

- base plate
- 1 LiPo battery module
- 1 AV-Modul
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the potentiometer module and adjust the resistance to  $60\Omega$ .
4. Measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  and note your values in the table. Use the AV-module in current-voltage-mode.
5. Decrease in several steps the resistance  $R_{\text{Pot}}$  at the potentiometer and measure each the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$ . Note all your values in the table.

**Advice:** Interrupt the current flow (for example by removing the cable) after each single measurement to avoid discharge of the battery module during the experiment.


### Evaluation

1. Enter your values in the given chart.
2. Compare your measured characteristics with the added characteristics and make a statement about the charge state of the cell. Calculate the remaining capacity of the battery module. You find instructions for this in experiment "Nominal voltage and capacity of voltage sources".
3. Name the main fields of application of LiPo batteries.
4. How can you optimize the life time of a LiPo battery?

## Measurements

R <sub>Pot</sub> in Ω									
V <sub>Load</sub> in V									
I <sub>Load</sub> in mA									

## Diagram



## Evaluation

[illegible]



## Evaluation

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.

[illegible]


## 4.6 I-V characteristics of the triple NiMH battery module

### Task

Determine the I-V characteristics of the triple NiMH battery module.

### Setup

### Equipment required



- base plate
- 1 NiMH battery module, triple
- 1 AV-Modul
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet.
2. Measure the open-circuit voltage  $V_0$  of the capacitor and note your value.
3. Plug in the potentiometer module and adjust the resistance to  $100\Omega$ .
4. Measure the voltage  $V_{Load}$  and current  $I_{Load}$  and note your values in the table. Use the AV-module in current-voltage-mode.
5. Decrease in several steps the resistance  $R_{Pot}$  at the potentiometer and measure each the voltage  $V_{Load}$  and current  $I_{Load}$ . Note all your values in the table.

**Advice:** Interrupt the current flow (for example by removing the cable) after each single measurement to avoid discharge of the battery module during the experiment.

### Evaluation

1. Enter your values in the given chart.
2. Compare your measured characteristics with the added characteristics and make a statement about the charge state of the cell. Calculate the remaining capacity of the battery module. You find instructions for this in experiment "Nominal voltage and capacity of voltage sources".
3. Calculate the total voltage and capacity of a series connection of two batteries with each 12V open circuit voltage and a capacity of 50Ah.




## 4.6 I-V characteristics of the triple NiMH battery module

## Measurements

 $V_0 =$  \_\_\_\_\_


$R_{Pot}$ in $\Omega$										
$V_{Load}$ in V										
$I_{Load}$ in mA										

## Diagram



## Evaluation

2.



## Evaluation

[illegible]


## 5.1 The charging process of the NiMH battery

### Task

Record the charge characteristics of a NiMH battery module.

### Setup

### Equipment required



- base plate
- 1 ChargerModule
- 1 NiMH battery module, single
- 1 AV-Module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in NiMH-mode. For handling instructions of the ChargerModule see page 5. Do not switch on the ChargerModule yet.
2. Measure the open-circuit voltage  $V_0$  of the battery module and note your value.

**Advice:** The battery module should have a charge state of maximum 50% (this corresponds to an open-circuit voltage of 1.18V). If the charge state is over 50% you can discharge the battery module with the resistances or the electric car.

3. Switch on the ChargerModule and measure the voltage  $V_{Load}$  and current  $I_{Load}$  in intervals of 1min and note your values in the table. Use the AV-module in current-voltage-mode.

### Evaluation

1. Enter your values in the diagram.
2. Describe and give reasons for the behavior of voltage and current during the charging process.
3. Explain the so-called *Memory-Effect* and the *Lazy-Effect*. Describe their influence on conventional NiMH batteries.



5.1 The charging process of the NiMH battery

Measurements


$V_o(1) =$  \_\_\_\_\_

t in min													
V in V													
I in mA													

t in min													
V in V													
I in mA													

t in min													
V in V													
I in mA													

Diagram





## 5.1 The charging process of the NiMH battery

### Evaluation

2.

---

---

---

---

---

---

---

---

---

---

3.

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---




## 5.2 The charging process of the NiZn battery

### Task

Record the charge characteristics of a NiZn battery module.

### Setup

### Equipment required



- base plate
- 1 ChargerModule
- 1 NiZn battery module
- 1 AV-Module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in NiZn-mode. For handling instructions of the ChargerModule see page 5. Do not switch on the ChargerModule yet.
2. Measure the open-circuit voltage  $V_0$  of the battery module and note your value.

**Advice:** The battery module should have a charge state of maximum 20% (this corresponds to an open-circuit voltage of 1.4V). If the charge state is over 20% you can discharge the battery module with the resistances or the electric car.

3. Switch on the ChargerModule and measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  in intervals of 10s and note your values in the table. Use the AV-module in current-voltage-mode.

### Evaluation

1. Enter your values in the diagram.
2. Describe and give reasons for the behavior of voltage and current during the charging process.
3. Determine the time after which the transition from cc-mode (constant current) to cv-mode (constant voltage) occurs.
4. Why does the voltage in cv-mode further increase slightly (despite an applied constant voltage)?

### Measurements


$V_0(1) = \underline{\hspace{2cm}}$

t in s													
V in V													
I in mA													



## 5.2 The charging process of the NiZn battery

### Diagram



### Evaluation

2.

---

---

---

---

---

---

---

---

3.

4.

---

---

---

---

---

---

---

---


## 5.3 The charging process of the LiFePo battery

### Task

Record the charge characteristics of a LiFePo battery module.

### Setup

### Equipment required



- base plate
- 1 ChargerModule
- 1 LiFePo battery module
- 1 AV-Module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in LiFePo-mode. For handling instructions of the ChargerModule see page 5. Do not switch on the ChargerModule yet.
2. Measure the open-circuit voltage  $V_0$  of the battery module and note your value.

**Advice:** The battery module should have a charge state of maximum 50% (this corresponds to an open-circuit voltage of 3.3). If the charge state is over 50% you can discharge the battery module with the resistances or the electric car.

3. Switch on the ChargerModule and measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  in intervals of 1min and note your values in the table. Use the AV-module in current-voltage-mode.

### Evaluation

1. Enter your values in the diagram.
2. Describe and give reasons for the behavior of voltage and current during the charging process.
3. Determine the time after which the transition from cc-mode (constant current) to cv-mode (constant voltage) occurs.
4. Describe the influence of the depth of discharge (DOD) on the lifetime of LiFePo batteries.



## 5.3 The charging process of the LiFePo battery

## Measurements


 $V_o(1) = \underline{\hspace{2cm}}$ 

t in min												
V in V												
I in mA												

t in min												
V in V												
I in mA												

t in min												
V in V												
I in mA												

## Diagram





## 5.3 The charging process of the LiFePo battery

### Evaluation

2.

---

---

---

---

---

---

---

---

---

---

3.

4.

---

---

---

---

---

---

---

---

---

---


## 5.4 The charging process of the lead battery

### Task

Record the charge characteristics of a lead battery module.

### Setup

### Equipment required



- base plate
- 1 ChargerModule
- 1 lead battery module
- 1 AV-Module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in lead-mode. For handling instructions of the ChargerModule see page 5. Do not switch on the ChargerModule yet.
2. Measure the open-circuit voltage  $V_0$  of the battery module and note your value.

**Advice:** The battery module should have a charge state of maximum 50% (this corresponds to an open-circuit voltage of 2.03). If the charge state is over 50% you can discharge the battery module with the resistances or the electric car.

3. Switch on the ChargerModule and measure the voltage  $V_{Load}$  and current  $I_{Load}$  in intervals of 1min and note your values in the table. Use the AV-module in current-voltage-mode.

### Evaluation

1. Enter your values in the diagram.
2. Describe and give reasons for the behavior of voltage and current during the charging process.
3. Determine the time after which the transition from cc-mode (constant current) to cv-mode (constant voltage) occurs.



5.4 The charging process of the lead battery

Measurements


$V_o(1) =$  \_\_\_\_\_

t in min													
V in V													
I in mA													

t in min													
V in V													
I in mA													

t in min													
V in V													
I in mA													

Diagram





## 5.4 The charging process of the lead battery

### Evaluation

2.

---

---

---

---

---

---

---

---

---

---

---

---

3.




## 5.5 The charging process of the lithium-polymer battery

### Task

Record the charge characteristics of a lithium-polymer battery module.

### Setup

### Equipment required



- base plate
- 1 ChargerModule
- 1 lithium-polymer battery module
- 1 AV-Module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Use the ChargerModule in lithium-polymer-mode. For handling instructions of the ChargerModule see page 5. Do not switch on the ChargerModule yet.
2. Measure the open-circuit voltage  $V_0$  of the battery module and note your value.

**Advice:** The battery module should have a charge state of maximum 75% (this corresponds to an open-circuit voltage of 3.9). If the charge state is over 75% you can discharge the battery module with the resistances or the electric car.

3. Switch on the ChargerModule and measure the voltage  $V_{\text{Load}}$  and current  $I_{\text{Load}}$  in intervals of 1min and note your values in the table. Use the AV-module in current-voltage-mode.

### Evaluation

1. Enter your values in the diagram.
2. Describe and give reasons for the behavior of voltage and current during the charging process.
3. Determine the time after which the transition from cc-mode (constant current) to cv-mode (constant voltage) occurs.
4. Why does the voltage in cv-mode further increase slightly (despite an applied constant voltage)?
5. Why it is dangerous to apply a purely cc-charge method for Lithium-based batteries?



5.5 The charging process of the lithium-polymer battery

Measurements


$V_o(1) =$ \_\_\_\_\_

t in min														
V in V														
I in mA														

t in min														
V in V														
I in mA														

t in min														
V in V														
I in mA														

Diagram





## 5.5 The charging process of the lithium-polymer battery

### Evaluation

2.

---

---

---

---

---

---

---

---

---

---

3.

4.

---

---

---

---

---

---

---

---

---

---

5.

---

---

---

---

---

---

---

---

---

---




## 5.6 The discharging process of a battery module

### Task

Record the discharging curve of a battery module.

### Setup



### Equipment required

- base plate
- 1 NiMH battery module, single
- 1 AV-Modul
- 1 potentiometer module
- cables

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer module yet to avoid the start of the experiment without recording data.
2. Measure the open-circuit voltage  $V_0(1)$  of the battery module and note your value.
3. Adjust the resistance of the potentiometer so that a discharge current of 250mA is flowing. The AV-module on the base plate is operated in current-voltage-mode.
4. Measure 10min the voltage  $V$  and current  $I$  in intervals of 1min and note your values in the table. Adjust the potentiometer if necessary to keep the discharge current constant. Remove after ten minutes the cable from the base unit.
5. Measure 5min after termination of the experiment again the open-circuit voltage  $V_0(2)$ .

**Advice:** Before the experiment the battery module should have a charge state of minimum 75% (this corresponds to an open-circuit voltage of 1.26V)

### Measurements

$V_0(1) =$  \_\_\_\_\_


$V_0(2) =$  \_\_\_\_\_

t in min	0	1	2	3	4	5	6	7	8	9	10
V in V											

### Evaluation

1. Enter your values in the diagram.
2. Calculate the capacity of the battery module before and after the experiment from the open-circuit voltage. You find instructions in experiment "Nominal voltage and capacity of voltage sources".
3. Name reasons for the deep-discharge of battery modules. Describe possibilities to protect the modules from this process.

## Diagram



## Evaluation

2.

A full-page sheet of white graph paper with a light gray grid. The grid consists of small squares, approximately 10 units wide by 10 units high, covering the entire area of the page. There are no margins or other markings on the paper.

3.

---

---

---

---

---

---




## 6.1 Hydrogen production in the reversible hydrogen fuel cell

### Task

Investigate the hydrogen production in a reversible hydrogen fuel cell.

### Setup



### Equipment required

- base plate
- 1 reversible hydrogen fuel cell
- 1 ChargerModule
- 1 AV-Module
- 1 Potentiometer module
- cables

### Procedure

1. Fill the reversible hydrogen fuel cell with distilled water. You find handling instructions on page 11.
2. Set up the experiment according to the circuit diagram. Pay attention to the polarity of the connections.
3. The Charger module is plugged into the base unit rotated by 90 ° (see sketch). Use the ChargerModule in Electrolyzer-mode. For handling instructions see page 5.
4. Adjust the resistance of the potentiometer to 2Ω and measure voltage V and current I at the fuel cell. The AV-module on the base plate is operated in current-voltage-mode.
5. Repeat the experiment for different resistances at the potentiometer (see table) and measure each the time which is required for the production of 3 ml of hydrogen (H<sub>2</sub>).

**Advice:** Make sure that the circuit is open before the start of each measurement (for ex. by removing a cable) so that the experiment does not start without recording data.

### Evaluation

1. Enter your values in the diagram.
2. How does the oxygen production rate cohere with the current of a fuel cell?
3. Explain the behavior of voltage and current in dependence from the resistance.


### Data

R in Ω	2	4	6	8
V in V				
I in mA				
t (3ml H <sub>2</sub> ) in min				



## 6.1 Hydrogen production in the reversible hydrogen fuel cell

## Diagram



## Evaluation

2.

---

---

---

---

---

3.

---

---

---

---

---

---

---

---

---

---




## 6.2 Characteristic curve of the electrolyzer

### Task

Use the electrolyzer to produce hydrogen and record the corresponding I-V-curve.

### Setup



### Equipment needed

- base plate
- 1 ChargerModule
- 1 AV-Module
- 1 reversible fuel cell
- 1 potentiometer module
- cables

### Execution

1. Fill the reversible hydrogen fuel cell with distilled water. You find handling instructions on page 11.
2. Set up the experiment according to the circuit diagram. Pay attention to the polarity of the connections.
3. The Charger module is plugged into the base unit rotated by 90 ° (see sketch). Use the ChargerModule in constant voltage-mode at 3V. For handling instructions see page 5.
4. Adjust the potentiometer to the maximum resistance of 110Ω and measure the voltage V and current I at the reversible fuel cell. Use the AV-module in current-voltage-mode. Note your values in the table.
5. Decrease the resistance at the potentiometer module in several steps and measure each the voltage V and current I. Note your values in the table.

**Advice:** The current circuit should be open at the beginning (for example by removing a cable) to avoid the start of the experiment without recording data.

### Measurements


V in V									
I in mA									

V in V									
I in mA									

## Evaluation

- Diagram



## Evaluation

2.

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and extend across the width of the page. There are no margins, text, or other markings on the paper.




## 6.3 Hydrogen consumption of a fuel cell

### Task

Determine the hydrogen consumption of the reversible fuel cell.

### Setup



### Equipment needed

- base plate
- 1 AV-Module
- 1 reversible fuel cell
- 1 potentiometer module
- cables

### Execution

1. Before starting the experiment an amount of 8ml hydrogen has to be produced. You find handling instructions in experiment „Hydrogen production in the reversible hydrogen fuel cell”.
2. Set up the experiment according to the circuit diagram. Use the AV-module in current-voltage-mode.
3. Adjust the potentiometer to a resistance of  $4\Omega$  and measure 5min voltage  $V$ , current  $I$  and the hydrogen consumption at the reversible fuel cell in intervals of 1min. Note your values in the table.
4. Refill the fuel cell with 12ml hydrogen and repeat the experiment with a potentiometer resistance of  $2\Omega$ .

### Evaluation

1. Enter your values in the diagram.
2. Describe the behavior of voltage  $V$ , current  $I$  and hydrogen production in the course of the experiment.
3. Describe the influence of the loading resistance on the operation of the fuel cell.

### Measurements

**$R=4\Omega$**

t in min	1	2	3	4	5
V in V					
I in mA					
H <sub>2</sub> in ml					




## 6.3 Hydrogen consumption of a fuel cell

## Measurements

 **$R=2.0\Omega$** 

t in min	1	2	3	4	5
V in V					
I in mA					
H <sub>2</sub> in ml					

## Diagram



## Evaluation

1.

---

---

---

---

3.

---

---

---


---

## 6.4 Characteristic curve of the fuel cell

### Task

Record the I-V-curve of a PEM fuel cell.

### Setup



### Equipment needed

- base plate
- 1 AV-Module
- 1 reversible fuel cell
- 1 potentiometer module
- cables

### Preparation

First, you have to produce appr.10ml of hydrogen. For handling instructions see experiment „Hydrogen production in the reversible hydrogen fuel cell“. Directly after  $H_2$  production the hydrogen fuel cell will behave like a capacitor. For this reason you should decrease its voltage down to approximately 0.9 V before measurement by letting a current of roughly 500 mA flow for 20 seconds.

### Execution

1. Set up the experiment according to the circuit diagram. Do not plug in the potentiometer yet.
2. Measure the open-circuit voltage of the fuel cell and note your value in the table.
3. Plug in the potentiometer and adjust it to the maximum resistance of 110Ω. Measure the voltage V and current I at the reversible fuel cell. Use the AV-module in current-voltage-mode. Note your values in the table.
4. Decrease the resistance at the potentiometer module in several steps and measure each the voltage V and current I. Note your values in the table.

### Evaluation

1. Enter your values in the diagram.
2. Describe the course of the I-V-characteristic.
3. Which area of the curve should be used for the operation of a consumer? Justify your answer.
4. Explain the decrease of voltage at higher currents.




### 6.4 Characteristic curve of the fuel cell

#### Measurements

V in V														
I in mA														

V in V														
I in mA														

#### Diagram



#### Evaluation

2.

---

---

---

---

---

---

---

---

---

---



### 6.4 Characteristic curve of the fuel cell

#### Evaluation

3.

---

---

---

---

---

---

---

---

4.

---

---

---

---

---

---

---

---

---

---




## 6.5 The efficiency of the hydrogen fuel cell

### Task

Determine the efficiency of the reversible hydrogen fuel cell.

### Setup



### Equipment required

- base plate
- 1 reversible hydrogen fuel cell
- 1 potentiometer module
- 1 AV-module
- cables

### Preparation

First, you have to produce appr. 5 ml of hydrogen. For handling instructions see experiment „Hydrogen production in the reversible hydrogen fuel cell“. Directly after  $H_2$  production the hydrogen fuel cell will behave like a capacitor. For this reason you should decrease its voltage down to approximately 0.9 V before measurement by letting a current of roughly 500 mA flow for 20 seconds.

### Procedure

1. Set up the experiment according to the circuit diagram.
2. Adjust the resistance of the potentiometer to  $5\Omega$ . Use the AV-Module in current-voltage-mode.
3. Measure the time it takes for the circuit to use up 2 ml of  $H_2$  and record voltage  $V$  and current  $I$  at this point.

### Data

$V =$  \_\_\_\_\_

$I =$  \_\_\_\_\_

$t =$  \_\_\_\_\_

### Evaluation

1. Calculate the electrical energy which was used up during the experiment. The electrical energy can be found using the following formula:

$$W_2 = V \cdot I \cdot t$$

2. For how long could the current flow with an entire filling of  $H_2$  (12 ml)?




## 6.5 The efficiency of a hydrogen fuel cell

## Evaluation

3. Determine the efficiency of the reversible hydrogen fuel cell. The efficiency of the reversible fuel cell is given by:

$$\eta = \frac{W_2}{W_1}$$

(The lower fuel value of 2ml of  $H_2$  is  $W_1 = 22Ws$ )






### 7.1 Operation of the electric car with several battery modules

#### Task

Observe the driving behavior of the car with different battery modules and conclude the characteristics from it.

#### Setup



#### Equipment required

- Electric car with module plate
- 1 AV-Module
- lead battery module
- NiZn battery module
- NiMh battery module, single
- LiFePo battery module
- LiPo battery module
- capacitor module
- Stop watch

#### Preparation

For the experiment you need enough space (min. 2x2m). Tilt the front axle of the car to the left, so that the car drives a circular path. Mark the starting and the finishing line of the car on the circular path with adhesive tape or something like that. The battery modules should be fully loaded and the capacitor module should be loaded to 5V before starting the experiment.

#### Execution

Execute the instructions for every battery module:

1. Measure the open circuit voltage  $V_{OC}$  of the battery module and record your data in the table.
2. Plug the battery module onto the module plate and first connect only **one** cable.
3. Position the car at the starting line and connect the second cable shortly before putting down the car.
4. Measure the time that the car needs for 4 rounds and repeat the measurement several times without stopping the car. Record your data in the table.
5. Let the car drive for 5 minutes and note your observations.
6. Calculate the difference to the previous round to determine the time for 4 rounds.

Advice: Pay attention to the car. It should not hit something, because the axles could get damaged. Hold the car shortly before starting it, because it could tip otherwise.

#### Evaluation



1. Compare the different battery modules and give reasons for the differences. Which properties of the respective module can you conclude from the differences?
2. Which parameters have influence on the measurement?
3. Which type of battery module would you classify as suitable for the use of an electric car?
4. Why should you load the capacitor to max. 5V to gain reasonable results?

## Data

	4 rounds	8 rounds	12 rounds	16 rounds	20 rounds	Observation after 5min (resp. when the car stops)
--	----------	----------	-----------	-----------	-----------	--

**Lead battery module  $V_0 =$  \_\_\_\_\_**

time in s						
time for 4 rounds						

**NiZn battery module  $V_0 =$  \_\_\_\_\_**

time in s						
time for 4 rounds						

**NiMH battery module  $V_0 =$  \_\_\_\_\_**

time in s						
time for 4 rounds						

**LiFePo battery module  $V_0 =$  \_\_\_\_\_**

time in s						
time for 4 rounds						

**LiPo battery module  $V_0 =$  \_\_\_\_\_**

time in s						
time for 4 rounds						

**Capacitor module  $V_0 =$  \_\_\_\_\_**

time in s						
time for 4 rounds						



## 7.1 Operation of the electric car with several battery modules

### Evaluation

1.

---

---

---

---

---

---

---

---

2.

---

---

---

---

---

---

---

---

3.

---

---

---

---

---

---

---

---

4.

---

---

---

---

---

---

---

---




## 7.2 Operation of the electric car with the reversible fuel cell

### Task

Observe the driving behavior of the car with the fuel cell and conclude the characteristics from it.

### Setup



### Required devices

- 1 Electric car with module plate
- 1 reversible fuel cell with mount
- 1 Stop watch

### Preparation

For the experiment you need enough space (min. 2x2m). Tilt the front axle of the car to the left, so that the car drives a circular path. Mark the starting and the finishing line of the car on the circular path with adhesive tape or something like that. Produce 12ml of  $H_2$  with the fuel cell (see. exp.5.1).

### Execution

1. Measure the open circuit voltage  $V_{OC}$  of the fuel cell after the production of 12ml  $H_2$  and record your data in the table.
2. Plug the fuel cell module onto the car and first connect only **one** cable.
3. Position the car at the starting line and connect the second cable shortly before putting down the car.
4. Measure the time that the car needs for 4 rounds and repeat the measurement several times without stopping the car. Record your data in the table.
5. Let the car drive for 5 minutes and note your observations.
6. Calculate the difference to the previous round to determine the time for 4 rounds.

Advice: Pay attention to the car. It should not hit something, because the axles could get damaged. Hold the car shortly before starting it, because it could tip otherwise.



## 7.2 Operation of the electric car with the reversible fuel cell

## Evaluation

1. Compare the operation of the electric car with the fuel cell to the operation with conventional accumulators like in the prior experiment.
2. Inform yourself about the application of fuel cells in the automotive industry. Which forms of storage of hydrogen are in use?

## Data

	4 rounds	8 rounds	12 rounds	16 rounds	20 rounds	Observation after 5min (resp.when the car stops)
--	----------	----------	-----------	-----------	-----------	---

Fuel cell:  $V_0 =$  \_\_\_\_\_

time in s						
time for 4 rounds						

## Evaluation

1.

---

---

---

---

---

---

---

---

2.

---

---

---

---

---

---

---

---



leXsolar GmbH  
Strehleener Straße 12-14  
01069 Dresden / Germany

Telefon: +49 (0) 351 - 47 96 56 0  
Fax: +49 (0) 351 - 47 96 56 - 111  
E-Mail: [info@lexsolar.de](mailto:info@lexsolar.de)  
Web: [www.lexsolar.de](http://www.lexsolar.de)