PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Telefon +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.de

12921-00 Cobra SMARTsense pH - Sensor zur Messung des pH-Wertes 0 ... 14 (Bluetooth)

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12921-00 Cobra SMARTsense pH - Sensor for measuring the pH value 0 ... 14 (Bluetooth)

Storage solution: 3 M KCl. The probe should be stored permanently in this solution.

Procedure after use: After use, replace the plastic cap with the storage solution on the sensor. The probe remains in the storage solution throughout the entire storage period.

Changing the storage solution (recommendation): Once per school year. Additional regular visual inspection: If the solution is no longer clear, it should be changed earlier.

Calibration: 2-point calibration with any two pH buffer solutions. (e.g. CHE-882395492)

Calibration frequency (recommendation): Once per school year.

12922-00 Cobra SMARTsense Conductivity - Sensor for measuring conductivity in liquids 0 ... 20000 $\mu\text{S}/\text{cm}$ / 0 ... 60°C (Bluetooth)


Storage solution: Distilled water. Fill the transparent plastic cap with this solution. A conductivity value of < 15 $\mu\text{S}/\text{cm}$ is optimal. The probe should be stored in this solution at all times.

Procedure after use: After using the sensor, replace the plastic cap with the storage solution. The sensor remains in the storage solution for the entire storage period.

Changing the storage solution (recommendation): Regularly, about once a month.

Calibration: 1-point calibration with a calibration solution of 1413 $\mu\text{S}/\text{cm}$. (e.g. 47070-02)

Calibration frequency (recommendation): Once per school year.

12933-01 Cobra SMARTsense Oxygen - Oxygen sensor for air and aqueous solutions 0 ... 20 mg/l (Bluetooth + USB)

Storage or electrolyte solution: 3 M KCl.

Filling and use:

The lower end of the electrode can be unscrewed. The entire cavity between the membrane cap and the electrode must be completely filled with 3 M KCl – without air bubbles. The sensor must also be completely filled with the KCl solution during measurement.

Changing the electrolyte/storage solution:

If the sensor is not used for a long period of time, the KCl solution should be removed and the sensor stored in a dry place. This protects the membrane from damage or contamination. Three replacement membranes are included in the scope of delivery.

Note: No recalibration is necessary after refilling; the last calibration is stored in the sensor.

Calibration: 2-point calibration with a zero oxygen solution and ambient air (20.9% O_2). (e.g. 12933-N, 47072-50)

Calibration frequency (recommendation): Once per school year.