

PHYWE Systeme GmbH & Co. KG Robert-Bosch-Breite 10 D–37079 Göttingen

Telefon +49 (0) 551 604-0 Fax +49 (0) 551 604-107 E-mail info@phywe.de Internet www.phywe.de

## Braunsches Rohr mit Fassung

06988-00

## **FAQ**

## The electron beam is not in the centre of the fluorescent screen

Yes, this is to be expected. When the sawtooth is switched off, there is anode
potential at both deflection plates. If you now remove the magnet and put the tube
upright, the point is in the middle of the screen.

## No electron beam is visible on the screen when the connection to the sawtooth generator is disconnected.

 This is to be expected. If the 2 deflector plates are potential-free, the visible beam travels at the height of the deflector plates. This can be seen well in a darkened environment.