

PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Telefon +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.com

Operating instructions



The unit complies
with the correspond-
ing EC guidelines




Fig. 1: 39839-01 Crystal lattice kit: diamond

Function and Applications

This kit includes all needed elements to build a diamant crystall model made of 30 carbon atoms.

Crystal lattice-kit: diamond					
Atoms	Holes (Angle)	Color	Dia. (mm)	Qty.	Drawing
C-Carbon	4(109°28')	Black	23	30	
Links		Color	Length (mm)	Qty.	Drawing
Medium Bond		Grey white	27	40	
ABS plastic box (size: 23.5 × 17 × 3.5 cm)				1	