

PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Telefon +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.com

Operating instructions



The unit complies
with the correspond-
ing EC guidelines



Fig. 1: 39841-01 Crystal lattice kit: fullerene (C60)

Function and Applications

This kit includes all needed elements to build a Buckminster fullerene (C60) model.

Crystal lattice kit: fullerene (C60)					
Atoms	Holes (Angle)	Color	Diam. (mm)	Amount	Drawing
C-Carbon	3(107°)	Black	23	60	
Links		Color	Length (mm)	Amount	Drawing
Medium Bond (single)		Grey white	27	60	
Medium Bond (double)		Purple	27	30	
ABS plastic box (size: 23.5 × 17 × 3.5 cm)				1	

PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Telefon +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.com

Operating instructions



The unit complies
with the correspond-
ing EC guidelines



Fig. 1: 39841-01 Crystal lattice kit: fullerene (C60)

Crystal lattice kit: fullerene (C60)

39841-01