

PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Telefon +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.com

Operating instructions



The unit complies
with the correspond-
ing EC guidelines




Fig. 1: 39848-01 Crystal lattice kit: ice

Function and Applications

This kit includes all needed elements to build an ice crystal model made of 35 water molecules.

Crystal lattice kit: ice					
Atoms	Holes (Angle)	Color	Dia. (mm)	Qty.	Drawing
H-Hydrogen	2(180°)	White	17	70	
O-Oxygen	4(109°28')	Red	23	35	
Links		Color	Length (mm)	Qty.	Drawing
Medium Bond		Purple	27	52	
Short Bond		White	11	70	
ABS plastic box (size: 23.5 × 17 × 3.5 cm)				1	