

Cornelsen Experimenta Molecular Construction Kits 1 and 2

Selection of Important Compounds

With Molecular Construction Kit 1, the following compounds can be built:

Alkanes (C_nH_{2n+2}):

1. Methane – CH₄
2. Butane – C₄H₁₀

Alkenes (Olefins, C_nH_{2n}):

3. Ethene (Ethylene) – C₂H₄

Alkynes (C_nH_{2n-2}):

4. Ethyne (Acetylene) – C₂H₂

Alkanols (Alcohols, C_nH_{2n+2}O):

5. Methanol (Methyl Alcohol) – CH₃OH
6. Ethanol (Ethyl Alcohol) – C₂H₅OH
7. Propanetriol (Glycerin) – CH₂OH–CHOH–CH₂OH

Alkanals (Aldehydes, R–CHO):

8. Methanal (Formaldehyde) – HCHO
9. Ethanal (Acetaldehyde) – CH₃CHO

Amino Acids and Proteins:

10. α-Aminopropanoic Acid (Alanine) – CH₃CHNH₂COOH
11. Formation of a dipeptide from two identical amino acids (glycine)

Polymers:

12. From ethylene, polyethylene is formed.

Halogenoalkanes:

13. Monochloromethane (Methyl Chloride) – CH₃Cl

Amino Acids and Proteins (continued):

14. Aminoethanoic Acid (Glycine) – NH₂CH₂COOH

Carbohydrates:

15. D-Glucose (Grape Sugar) – C₆H₁₂O₆
16. Fructose (Fruit Sugar) – C₆H₁₂O₆
17. Sucrose (Cane Sugar, Beet Sugar) – C₁₂H₂₂O₁₁
18. Starch – (C₆H₁₀O₅)_n
19. Cellulose – (C₆H₁₀O₅)_n

Carboxylic Acids:

20. Methanoic Acid (Formic Acid) – HCOOH
21. Ethanoic Acid (Acetic Acid) – CH₃COOH
22. Propanoic Acid (Propionic Acid) – C₂H₅COOH
23. Butanoic Acid (Butyric Acid) – C₃H₇COOH

General formula for monocarboxylic acids:

R–COOH

(The remainder group R is represented with a universal building block from Molecular Construction Kit 2.)

Selection of Important Compounds – Aromatic and Cyclic

For all the following compounds, both Kits 1 and 2 are required.

Aromatic Compounds:

- Benzene – C_6H_6
- Nitrobenzene – $\text{C}_6\text{H}_5\text{NO}_2$
- Naphthalene – C_{10}H_8

Acids:

- Nitric Acid – HNO_3

Dyes:

- Methyl Orange

Medicinal Compounds:

- 6-Amino-Penicillanic Acid

Materials List:

No.	Component	Quantity
1	Connecting pieces, grey	60
2	Hydrogen (H), monovalent, white	25
3	Oxygen (O), divalent, red	15
4	Chlorine (Cl), monovalent, green	5
5	Nitrogen (N), trivalent, blue	5
6	Carbon (C), tetravalent, black	14
7	Additional connecting pieces	80
8	Universal building blocks, grey	4
9	Benzene rings, black	9
10	Sulfur (S), divalent, yellow	8
11	Oxygen (O), divalent, red	4
12	Nitrogen (N), trivalent, blue	4
13	Carbon (C), tetravalent, black	6
14	Phosphorus (P), pentavalent, violet	4
15	Sulfur (S), hexavalent, yellow	4

Instructions for Use:

- The connecting pieces are pushed firmly into the atom spheres.
- To form double or triple bonds, bend the connectors strongly.
- After use, all parts should be placed back according to the storage diagram to easily check completeness.

Legal Notice:

© 2008 Cornelsen Experimenta, Berlin – All rights reserved.

Any use beyond what is legally permitted requires prior written consent.

This includes scanning or making the content publicly available on networks or school intranets.

Cornelsen assumes no liability for damages resulting from misuse of materials.

Contact:

Cornelsen Experimenta GmbH

Holzhauser Straße 76, Berlin

Service Tel.: 0800 435 90 20

Service Fax: 0800 435 90 22

Email: info@cornelsen-experimenta.de