

PHYWE Systeme GmbH & Co. KG
Robert-Bosch-Breite 10
D-37079 Göttingen

Telefon +49 (0) 551 604-0
Fax +49 (0) 551 604-107
E-mail info@phywe.de

Einraum-Übersicht General map to fill



Nr. No.	Menge Qty.	Deutsch	English	Art. Nr. Art. No.
TESS advanced Applied Sciences Set Akustik 1, AE-1			TESS advanced Applied Sciences set Acoustics 1, AE-1	25289-88
1	1	Lineal, $l = 200$ mm, Edelstahl	Ruler, stainless steel, $l = 200$ mm	09937-01
2	1	Software „measure Acoustics“	Software “measure Acoustics”	14441-61
	0,5	Siliconschlauch, $d_i = 3$ mm	Silicon tubing, $d_i = 3$ mm	39292-00
	1	PVC-Schlauch, $d_i = 7$ mm	PVC tubing, $d_i = 9$ mm	03985-00
3	1	Anschlaghammer, Gummi	Striking hammer	03429-00
4	2	Stimmgabel 440 Hz	Tuning fork 440 Hz	03424-00
5	1	Maßband, $l = 2$ m	Measuring tape, $l = 2$ m	09936-00
6	1	Stimmgabel 880 Hz	Tuning fork 880 Hz	03421-00
7	1	Seidenfaden, Nähseide, auf Röllchen, $l = 200$ m	Silk thread, $l = 200$ m	02412-01
8	2	Trichter, Oben- $d = 75$ mm	Filter funnel, PP, $d = 75$ mm	46895-00
9	1	Rahmentrommel, $d = 20$ cm	Frame drum, $d = 20$ cm	13289-11
	1	Styroporkugel mit Haken	Styrofoam sphere with hook	13289-13
10	1	Laborbecher, PP, niedrige Form, 100 ml	Beaker, 100 ml, low form, plastic	36011-01